Vortex IR

Instruction Manual (v.1)

N2400-0001

Vortex IR

CE

www.starlabgroup.com

S T A R
L A B

ENGLISH..... 3

DEUTSCH 8

FRANCAIS 13

ITALIANO 18

EC DECLARATION OF CONFORMITY

STARLAB certifies that

N2400-0001 Vortex IR

has been manufactured and inspected in accordance

with the following European standards:

CEI EN 61010-2-010 Safety requirements for electrical equipment

for measurement, control and for laboratory use

UL 61010-1 Electrical equipment for laboratory use

CAN/CSA-C22.2 No.61010-1 General requirement, Canadian electrical code

As an ISO 9001 and ISO 13485 certified company,
our working procedures and processes are checked and audited regularly
and the quality of our products is continuously monitored.

STARLAB International GmbH

VORTEX IR

Cat. No. N2400-0001

General Information

Before using the unit, please read this entire Instruction Manual carefully.

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

- 1. This unit must be used for laboratory applications only.** The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.
2. The manufacturer reserves the right to modify the characteristics of its products to constantly improve their quality.
3. Make sure that liquid does not come into contact with the power cable or the electrical parts inside the unit.
4. Use the transformer supplied with the unit.
5. Do not use the unit if it is damaged or malfunctions. If this occurs, contact your local STARLAB office or STARLAB distributor.
6. Personal protective equipment must be compatible with the possible risks posed by the material being processed and the material of the containers being used.

1. Introduction

Mixing by spinning can be carried out quickly and safely by placing the tube containing the sample on the special plastic mixing cup head. The Vortex IR offers the possibility of either automatic or continuous mode: the mixing speed can be adjusted using the knob (1). The instrument automatically starts mixing when the tube crosses the infrared sensor field (3) near the mixing cup head (4). For continuous mode press the key (2), stirring is continuous at the selected speed, no pressure is required on the cup head. The automatic mode indicator light (5) indicates that the instrument is running in automatic mode.

The mixing cup heads have a push fit attachment which allows quick and easy replacement with any of the accessories available for mixing tubes of different sizes. See section 3.1 for instructions on changing the heads.


Fig. 1

2. Assembling and installing

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box should contain:

- Vortex IR complete with cup head
- Plug adapter for power supply
- Switching 100–240 V / 12 V plug
- Instruction manual

2.1 Electrical connection

After unpacking place the unit on the laboratory bench. Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied. Ensure that the socket and the relative cut-off device conform to current safety norms and are within easy reach.

2.2 Start-up

Check that the mixing cup head for holding the tube is correctly inserted in place. Also check that the speed control knob (1) is turned to 0 rpm. Use the ON-OFF switch (6) on the side of the unit to turn the

unit ON. The start-up and the control of the speed are effected with the knob. When switched on, the unit automatically starts up in Sensor mode and the indicator light (5) is lit.

3. Operation controls

Connect the unit to mains and turn it ON using the ON-OFF switch (6). **Each time the unit is switched on it automatically starts in Sensor mode** and the sensor mode indicator light (5) will be on. To select continuous mode, press the Continuous-Sensor button (2) on the front of the unit and the indicator light (5) will go out to indicate the unit is in continuous mode. In Continuous mode, the mixing cup head runs continuously at the speed set.

When Sensor mode is selected, mixing starts automatically (at the speed set) when the tube crosses the sensor field (infrared optical system) as shown in Fig. 2, without the need for any pressure on the mixing cup head. The sensor field has been purposely designed to detect the presence of any kind of tube to be mixed. The Vortex IR allows the mixing of tubes without the need to exert pressure on the cup head, thus reducing effort by the operator. **NOTE:** remember to hold the tube tightly during mixing!


Fig. 2

<p>CONTINUOUS-SENSOR MODE BUTTON (2)</p>	<p>Pushing the Continuous-Sensor button (2) activates the respective modes:</p> <ul style="list-style-type: none"> - Continuous (indicator light off) The instrument runs in continuous mode with several accessories; no pressure is required. - Sensor (indicator light on) Allows automatic operation mode when a tube crosses the IR field; even with different accessories.
<p>SPEED CONTROL KNOB (1)</p>	<p>Use for quick and precise regulation of mixing speed. Range: 0 to 3000 rpm.</p>
<p>ON-OFF SWITCH (6)</p>	<p>Use to power the unit on and off. 0 position is OFF / I position is ON.</p>

Always disconnect the unit from the power supply when not in use.

3.1 Replacing the cup head

To replace the cup head, pull it upwards and off. The other push-fit attachments (see Chapter 6, Ordering Information) can then be inserted.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodic cleaning. In compliance with the product guarantee law, repairs to the unit must be carried out in our factory, unless previously agreed otherwise with local distributors (see Warranty, Chapter 7).

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Specifications

Speed range	0–3000 rpm
Speed setting	Analogue
Speed control	Electronic
Operating temperature range	5 °C to 40 °C (max. humidity 80%)
Operation modes	Two: sensor or continuous
Type of movement	Orbital (4.5 mm Ø)
Maximum load	0.5 kg
Construction materials	Unit: Zinc alloy die cast and technopolymer. Cup head: SEBS Foamed attachments: PE. Feet: MQ/MVQ. Control panel: PET
Power	15 W
Power supply	Input: AC 100–240 V, 50–60 Hz, 0.5 A. Output: DC 12 V, 1.25 A
Safety	Electronic protection degree IP 42-EN 60529. Pollution degree 2-EN61010-1
Weight	2.2 kg
Dimensions (WxDxH)	180 x 220 x 70 mm

Technical specifications subject to change without notice.

6. Ordering information

N2400-0001	Vortex IR
N2400-1008	Replacement Cup for Single Test Tube
N2400-1001	Vortex IR Foam Attachment for 19 x 1.5 ml Microcentrifuge Tubes
N2400-1002	Vortex IR Foam Attachment for 5 x 50 ml Tubes (or any tube up to 16 mm Ø)
N2400-1003	Vortex IR Foam Attachment for 4 x 15 ml Tubes (or any tube up to 29 mm Ø)
N2400-1004	Vortex IR Foam Attachment for Microtitre Plate
N2400-1005	Vortex IR Universal Foam Attachment (50 mm Ø)
N2400-1006	Vortex IR Customisable Soft Foam Attachment (29 mm Ø)

NOTE: Do not exceed a speed of 800 rpm when the microtiter foam attachment (N2400-1004) is being used. High speeds may cause sample spillage. Do not exceed a speed of 1200 rpm when all other foam attachments are being used. Only the cup head (N2400-1008) supplied with the unit should be used up to the maximum speed of 3000 rpm.

Please see our website www.starlabgroup.com for information on our wide range of consumables including centrifuge tubes, microcentrifuge tubes and PCR tubes.

7. Warranty

The Vortex IR is under warranty against defects in materials and workmanship under normal service, for two years from the shipping date to the purchaser. This warranty excludes damages resulting from transport, misuse, carelessness, or neglect. STARLAB's liability under the warranty is limited to the receipt of reasonable proof from the customer that the defect is embraced within the terms of the warranty. All claims made under this warranty must be presented to STARLAB or a STARLAB Distributor within two years following the date of delivery of the product. Our "General terms and conditions" apply.

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STARLAB GmbH
STARLAB (UK), Ltd
STARLAB S.r.l.