This is your **Aspirator** user manual.

USER MANUAL · BEDIENUNGSANLEITUNG · MODE D'EMPLOI

User Manual	3–22
Bedienungsanleitung	23–42
Mode d'emploi	43–61

1	Safety information	4
1.1	Important information	4
1.2	General information	5
1.3	Intended use	5
1.4	Setting up and installing the system	5
1.5	Ambient condition	6
1.6	Operating conditions	6
1.7	Safety during operation	6
1.8	Maintenance and repair	8
2	Technical data	9
2.1	Wetted parts	10
2.2	System parts	11
2.3	On/off switch	11
3	Use and operation	12
3.1	First steps: Installation	12
3.2	Operating the Aspirator	12
3.3	During operation	13
3.4	Filtration	13
3.5	Filter and collection bottle	13
3.6	Storage following use	14
4	Cleaning and decontaminating	15
5	Assembling of components	16
5.1	Replacing the filter	16
5.2	Assembling a second hand controller connection set	16
5.3	Assembling quick coupling bottle - pump unit	17
5.4	Quick coupling sets	17
6	Accessories - spare parts	18
7	Troubleshooting	19
7.1	Replacing the fuse	20
8	Repair - Maintenance - Return	20
8.1	Return form	21

Latest and up-to-date version of the declaration of conformity is available at:
www.starlab.click/doc-aspirator

1 Safety information

1.1 Important information

- > Keep this manual complete and accessible to personnel at all times!
- > Read this manual carefully before installing or operating the equipment. Observe the instructions contained in this manual.
- > Do not modify the equipment without authorization.

This manual is an integral part of the equipment described therein. It describes the safe and proper use of the vacuum pump. Make operating personnel aware of dangers arising from the pump and the pumped substances. Starlab cannot accept any liability for inappropriate use of these pumps or for damage from failure to follow instructions contained in this manual.

This manual is only to be used and distributed in its complete and original form. It is strictly the users' responsibility to check carefully the applicability of the guidance in this manual with respect to this product.

The following signal word panels and safety symbols are used throughout this manual:

This is the **safety alert** symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury and death.

DANGER indicates a hazardous situation which, if not avoided, will result in death or serious injury.

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE is used to address practices not related to personal injury.

Caution! Hot surface!

Disconnect equipment from AC power.

Within the European Community, the disposal of electrical devices is regulated by national regulations based on EU Directive 2012/19/EU pertaining to waste electrical and electronic equipment (WEEE). According to these regulations, any devices supplied after August 13, 2005, in the business-to-business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste. Because disposal regulations may differ from one country to another within the EU, please contact your supplier if necessary.

The document "Safety information for vacuum equipment" is part of this manual! Read the "Safety information for vacuum equipment" and observe the instructions contained therein!

1.2 General information

WARNING

- › Read and comply with this manual before installing or operating the equipment.
- › Before operating the equipment read and comply with the section "Cleaning and decontamination"!

NOTICE

Remove all packing material, remove the product from its packing-box, remove the protective covers from the inlet and outlet ports and keep. Inspect the equipment. If the equipment is damaged, notify the supplier within three days; state the item number of the product together with your order number and the supplier's invoice number. Retain all packing material for inspection. **Do not use the equipment if it is damaged.** Lift and transport the equipment by the handle.

1.3 Intended use

WARNING

- › The system and all system parts must not be used on humans or animals. Prevent any part of the human body from coming into contact with the vacuum.
- › Make sure that the individual components are only connected, combined and operated according to their design and as indicated in the instructions for use. Use only original Starlab accessories.
- › Comply with notes on correct vacuum and electrical connections, see section "Use and operation".

CAUTION

- › The systems are designed for **ambient temperatures** during operation between +10 °C and +40 °C. Check the maximum temperatures if installing the system in a cabinet or a housing and make sure ventilation is adequate. Install an external automatic ventilation system if necessary.

NOTICE

- › Ensure that the equipment and all components are suitable for the intended application. Use the system only for **aspiration of liquids** and **filtration**.

1.4 Setting up and installing the system

DANGER

- › Equipment must be connected only to a **suitable electrical supply** and a suitable ground point. **Failure to connect the motor to ground may result in deadly electrical shock.**

The supply cable may be fitted with a molded European IEC plug or a plug suitable for your local electrical supply. If the plug has been removed or has to be removed, the cable will contain wires color coded as follows: green or green and yellow: earth; blue or white: neutral; brown or black: live. The device has an internal fuse.

WARNING

- › Do not permit any **uncontrolDEL pressurizing** (e.g. make sure that an exhaust tubing cannot become blocked). **Risk of bursting!**
- › Due to the high compression ratio of the pumps, pressure at the outlet port may be generated that is higher than the maximum permissible pressure compatible with the mechanical stability of the system.
- › Always provide a free and pressureless exhaust line.
- › Keep the electrical power cord away from heated surfaces.
- › Keep the electrical power cord away from hot surfaces.

CAUTION

- › Comply with maximum permissible pressures and pressure differences; see section "Technical data". Do not operate the pumping unit with overpressure at the inlet or outlet.
- › Check that line voltage and current conform with the equipment. (See rating plate.)
- › Avoid overpressure of more than 0.2 bar in instances in which inert gas is connected.

NOTICE

Provide a firm level platform for the equipment. Ensure a stable position of the pump without any mechanical contact except of the pump feet. Comply with all applicable **safety regulations**. Check fan regularly for dust/dirt, clean if necessary to avoid reduced ventilation. Avoid high heat supply. If the equipment is brought from cold environment into a room for operation, allow the equipment to warm up. (Watch for water condensation on cold surfaces.) The mains plug is a disconnecting device to separate the pump from the supply voltage. Ensure that the mains plug is easily accessible at all times to allow the separation of the device from the power supply. The diameter of the an outlet tubing should be at the least as large as the diameter of the pump connections. Comply with **all applicable and relevant safety requirements** (regulations and guidelines), **implement the required actions and adopt suitable safety measures.**

1.5 Ambient conditions

NOTICE

To the best of our knowledge the equipment is in compliance with the requirements of the applicable EC-directives and harmonized standards (see "Declaration of conformity") with regard to design, type and model, especially directive DIN EN 61010-1. This directive gives, in detail, conditions, under which the equipment can be operated safely (see also IP degree of protection). Adopt suitable measures in case of differences, e.g. using the equipment outdoors, installation in higher altitudes, conductive pollution or bedewing. Pay attention to the **permissible maximum ambient temperatures** (see "Technical data").

1.6 Operating conditions

DANGER

- > **The devices have no approval for operation in, or for pumping of, potentially explosive atmospheres.**
- > The devices are **not suitable** to pump:
 - **unstable substances** and substances which react explosively under **impact** (mechanical stress) and/or when being exposed to **elevated temperatures** without air,
 - **self inflammable** substances,
 - substances which are **inflammable without air**
 - **explosive substance**
- > The pumps have **no approval** for operation below ground.

WARNING

- > The pumps are **not suitable** for pumping dust.

CAUTION

- > The devices are **not suitable** for pumping substances which may form **deposits** inside the pump. Deposits and condensate in the pump may lead to increased temperatures even to the point of exceeding the maximum permitted temperatures! Increased temperatures may cause ignition of inflammable substances that may have been deposited inside the pump during use.
- > If there is a danger of the formation of **deposits** in the pump chamber (check inlet and outlet of the pump), inspect the pump chamber regularly and clean if necessary.
- > **Take into consideration interactions and chemical reactions of the pumped media.**
- > Ensure that the materials of the wetted parts are compatible with the pumped substances, see section "Technical data".

NOTICE

If pumping **different substances**, it is recommended that the pump be purged with air or inert gas prior to changing the pumped media in order to pump out residues and to avoid reactions of the pumped substances with each other and with the pump materials.

1.7 Safety during operation

DANGER

- > Avoid interactions of media in the collection bottle. Comply with material safety data sheets and notes on safe use of the manufacturer. **Do not mix incompatible disinfectants and/or incompatible reagents / solvents or any unknown substances.**
For example sodium hypochlorite (chlorine bleach)

Incompatible chemicals and agents	Possible results of mixing with sodium hypochlorite (chlorine bleach)
Acids or acidic compounds (e.g. hydrochloric acid, aluminium chloride)	Release of chlorine gas
Ammonia containing compounds (e.g. ammonium hydroxide, quaternary ammonium salts)	Formation of explosive compounds, release of chlorine gas and other hazardous gases
Organic chemicals (e.g. solvents, polymers, amines, oils)	Formation of chlorinated organics, release of chlorine gas and other hazardous gases
Metals (e.g. copper, iron) Hydrogen peroxide	Release of oxygen, overpressure, rupture of a closed system
Reducing agents (e.g. sodium thiosulfate)	Production of heat, boiling
Guanidine salts (e.g. guanidine hydrochloride, guanidine thiocyanate)	Release of toxic gases, e.g. chlorine, chloramine, hydrogen cyanide

 **DANGER**

- > Avoid the formation of dangerous reactions / gases in the Aspirator, especially in the collection bottle. If this is impossible, dispose of dangerous gases at the outlet of the pump appropriately.
- > Adopt suitable measures to prevent the release of dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapors and gases. In such cases, install an appropriate collecting and disposal system and take protective action for pump and environment.
- > If the equipment is used in combination with **hazardous materials** (e.g., in **medical-microbiological laboratories**) check all relevant safety and health requirements and determine the applicability of regulatory limitations prior to use if necessary.
- > Use appropriate safety work materials and methods if necessary, e.g., steam sterilization, sterilization indicators or germicides. For notes concerning the sterilization of wetted parts, see the "Technical Data" section and "Cleaning and decontaminating". It is the user's responsibility to ensure effective sterilization.
- > Never operate a defective or damaged device.
- > The user must take suitable precautions to prevent any formation of explosive mixtures in the expansion chamber or in the housing. In case of a diaphragm crack, for example, mechanically generated sparks, hot surfaces or static electricity may ignite these mixtures. Use inert gas for venting if necessary.
- > Potentially explosive mixtures at the outlet of the pump have to be drained appropriately, evacuated, or diluted with inert gas to non-explosive mixtures.

 **WARNING**

- > Check the collection bottle regularly for cracks. Never use a collection bottle with cracks nor expose it to vacuum.
- > Use the system only in combination with the integrated hydrophobic filter (sterilizable) to protect pump and environment (people).
- > Comply with applicable regulations when disposing of chemicals. Take into consideration that chemicals may be polluted. Take adequate precautions to protect people from the effects of dangerous substances (chemicals, thermal decomposition products of fluoroelastomers), wear appropriate safety-clothing and safety glasses.
- > Use only **Starlab spare parts and accessories**; otherwise safety and performance of the equipment as well as the electromagnetic compatibility of the equipment might be reduced. The CE mark (see rating plate) may become void if original equipment manufacturer parts are not used.
- > In case of leaks in the manifold, aspirated substances might be released into the environment. Comply especially with notes on operation and use and maintenance.
- > Due to the residual **leak rate of the equipment**, there may be an exchange of gas, albeit extremely slight, between the environment and the vacuum system. Adopt suitable measures to prevent contamination of the pumped substances or the environment.

NOTICE

Maximum liquid level in collection bottle: approx. 80 %, depending on the application. (When working with low boiling liquids, or liquids that tend to foam, the maximum liquid level might be reduced).

Ensure that, in case of failure, the system will always return to a safe status. Provide appropriate protective measures (i.e. precautions which allow for the requirements of the respective application) even for the case of failure and **malfunction**. When operating the pump, the loss of pumping or venting must not lead to a critically dangerous situation under any circumstances.

Electronic equipment is never 100 % fail-safe. This may lead to an indefinite status of the equipment. Failure of the pumping unit (e.g., by power failure) or connected components, or change of parameters must not lead to a critically dangerous situation under any circumstances. In case of diaphragm cracks or leaks in the manifold, pumped substances might be released into the environment or into the pump housing. Comply especially with notes on operation and use and maintenance.

- > Pay attention to the safety symbol "hot surfaces" on the equipment. Depending on operation conditions and ambient conditions dangers due to hot surfaces may arise. Hot parts may cause burns if touched. Adopt suitable measures to prevent any danger arising from hot surfaces. Ensure that hot surfaces of the pump do not cause burns. Provide a suitable contact guard if necessary.

In case of overload, the motor is shut down by a **self-hold thermal cutout** in the winding. **Attention:** Manual reset is necessary. Switch off the pump or isolate the equipment from power supply. Identify and eliminate the cause of failure. Wait approx. five minutes before restarting the pump.

1.8 Maintenance and repair

NOTICE

Parts liable to wear have to be replaced regularly. Under typical operating conditions, the projected lifetime of the diaphragms and valves is > 15000 operating hours. Bearings have a typical durability of 40000 h. Motor capacitors have a typical durability in the range of 10000 to 40000 h depending strongly on the operating conditions like ambient temperature, humidity or load.

DANGER

- > Never operate the pump if covers or other parts of the pump are disassembled. Ensure that the pump cannot be operated accidentally.
- > **Isolate equipment from power supply before removing the cover!**
- > Before starting maintenance, unplug the equipment and wait 5 seconds to allow the capacitors to discharge.
- > Attention: The pump might be contaminated with noxious or otherwise dangerous process chemicals that have been pumped during operation. Ensure that the pump is decontaminated before maintenance.

WARNING

- > Take adequate precautions to protect people from the effects of dangerous substances if contamination has occurred. Wear appropriate safety-clothing when you come in contact with contaminated components.
- > **Parts liable to wear have to be replaced regularly.**
- > Never operate a defective or damaged pump.
- > Check every motor capacitor regularly by measuring its capacity and estimating its operation time. Exchange old capacitors early enough to prevent a failure. When motor capacitors fail, they can get hot enough to melt or cause a flame which could be **dangerous for persons and equipment in the vicinity**. Motor capacitors have to be replaced by an electrician.
- > Before starting maintenance, vent the system and allow sufficient cooling of the pump. Drain condensate, if applicable. Avoid the release of pollutants.

NOTICE

Clean polluted surface with a clean, slightly moistened cloth. To moisten the cloth we recommend water or mild soap. Ensure that **maintenance** is carried out by suitably trained technicians. Ensure that the maintenance technician is familiar with the safety procedures which relate to the products processed by the pumping system. In order to comply with law (occupational, health and safety regulations, safety at work law and regulations for environmental protection) vacuum pumps, components and measuring instruments returned to the manufacturer can be repaired only when certain procedures (see section "**Repair - Maintenance - Return**") are followed.

- > Work on the equipment may be performed only by a qualified person
- > Work on electrical equipment in particular may be performed only by qualified electricians.

2 Technical data

Type	Aspirator
Pump	
Maximum pumping speed (ISO 21360) at 50/60 Hz	0.7 / 0.8 m ³ /h
Ultimate vacuum (absolute)	150 mbar
Maximum permissible inlet pressure (absolute)	1.1 bar
Maximum permissible outlet pressure (absolute)	1.1 bar
Rated power	0.04 kW
No-load speed 50/60 Hz	1500 / 1800 rpm
Control range	150 - 850 mbar
Permissible ambient temperature storage / operation	-10 to +60 / +10 to +40 °C
Maximum installation altitude m 2000 above mean sea level	2000 m above mean sea level
Permissible relative atmospheric moisture during operation (no condensation)	30 to 85 %
Maximum permissible range of supply voltage: Attention: Observe specifications of rating plate!	230 V~ ±10% 50/60 Hz
Rated current at: 230 V~ 50/60 Hz	0.8 / 0.6 A
Fuse	2 x 2.5 AT, 0.19 x 0.78 in (5 x 20 mm)
Motor protection	Thermal cutout
Overvoltage category	II
Degree of protection IEC 60529	IP 40
Pollution degree	2
A-weighted emission sound pressure level* (uncertainty K _{pA} : 3 dB(A))	49 dB(A)
System	
Volume of collection bottle	4 l (polypropylene)
Connection at the inlet	Hose nozzle DN 3/5 mm
Connection at the outlet	Hose nozzle DN 8 mm / silencer
Hydrophobic filter	
Type	Midisart® 2000
Manufacturer	Sartorius stedim
Nominal pore diameter	0.2 µm
Dimensions L x W x H approx. Aspirator	408 x 194 x 500 mm
Weight approx. Aspirator	7.7 kg


* Measurement according to EN ISO 2151:2004 and EN ISO 3744:1995 at 230V/50Hz and ultimate vacuum with exhaust tube at outlet. Subject to technical changes!

2.1 Wetted parts

Components	Wetted materials
Pump	
Housing cover insert	PTFE, carbon reinforced
Head cover	ETFE, carbon fibre reinforced
Diaphragm clamping disc	ETFE, carbon fibre reinforced
Diaphragm	PTFE
Valve	PTFE / FFKM
Inlet	ETFE
Outlet	ETFE
Silencer	Silicone rubber
Filter	
Diaphragm	PTFE
Housing	PP
Tubing	Silicone rubber
Collection bottle 4l	
Bottle / screw cap	PP
Hose nozzle at bottle head (connection hand controller)	PP
Hose in bottle	PTFE
Hose nozzle (connection filter)	PPS, glass fibre reinforced
Counter nut	PP
Closing screw	PPS, glass fibre reinforced
Collection bottle 2l (accessory)	
Bottle	Borosilicate glass
Cap insert	PP
Flat seal	EPDM
Hose nozzle at bottle head (connection hand controller)	PP
Hose in bottle	PTFE
Hose nozzle (connection filter)	PPS, glass fibre reinforced
Closing screw	PPS, glass fibre reinforced

2.2 System parts

Position	Designation
1	Pump
2	Mains connection
3	Fuse holder
4	On / Off switch
5	Touch panel
6	Rating plate
7	Outlet
8	Connection tubing
9	Hydrophobic protection filter
10	Connection filter
11	Connection hand controller
12	Closing screw (optional: connection second hand controller)
13	Screw cap / bottle cap with insert
14	Collection bottle
15	Handle

2.3 On/off switch

System switched off


System switched on

Attention: Operate the on/off switch only by hand.

3 Use and operation


3.1 First steps: Installation


- > Unpack the equipment.

- > Read and comply with the document "**Safety information for vacuum equipment**"!

- > Assemble connection tubing with filter.

- > Setting up the Aspirator. Connect to power supply.

Check line voltage and current prior to switching on!


- > Attach aspiration hose of the hand controller to the hose nozzle of the bottle head.

3.2 Operating the Aspirator

The unit has a touch panel so the keys need only be touched gently. Do not press the keys.

The keys "+" and "-" have to be touched > 0.25 seconds to be actuated.

A successful action is confirmed by a blip and the flashing of LEDs.

Position	Designation
1	Key to reduce suction power
2	Display suction power
3	Key to increase suction power

Use the keys to set suction power of the system.

The suction power can be set linearly in a range from 150 mbar (1 LED flashes) to 850 mbar (8 LEDs flash) underpressure (relative to atmospheric pressure).

A **flashing LED** indicates that the actual available suction power differs from the preset suction power.

Attention: If the collection bottle is under vacuum and the vacuum demand is reduced, the existing underpressure inside the bottle remains until the vacuum is reduced (pressure rises) by further aspiration with the hand controller.

NOTICE


Depending on ambient atmospheric pressure (depending on altitude or weather conditions) and on the vapor pressure of the media inside the collection bottle, the setting of a high underpressure (e.g., 850 mbar, 8 LEDs flash) may lead to continuous operation of the pump. In this case, a reduction of the suction power is recommended.

3.3 During operation

DANGER

- › Attention: When removing the hand controller, liquid in the tube may leak! Adopt suitable measures to prevent the release of dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapors and gases. In such cases, use an appropriate collecting and disposal system and take protective action for personnel, equipment and environment.


WARNING

- › **Max. ambient temperature: 40 °C**

CAUTION

- › Use the system only with the **integrated hydrophobic filter** to protect the vacuum supply from aspirated liquids and aerosols, and to protect the environment/user from contamination risk.

Silencer at the outlet

Attention: Dust-laden gases, deposits and condensed solvent vapor can restrict air flow out of the silencer. The resultant back pressure can lead to damage of pump bearings, diaphragms, and valves. Under those conditions, a silencer must not be used. Check the silencer regularly and replace if necessary. In case of harmful / poisonous gases or condensate at the outlet, remove silencer and replace with an exhaust tube.


NOTICE

The hand controller allows the aspiration of excess liquids. See manual "Hand controller".

The aspiration system does not feature an air admittance valve. If the collection bottle is under vacuum and the vacuum demand is reduced, the existing underpressure inside the bottle remains until the vacuum is reduced (pressure rises) by further aspiration with the hand controller.

Attention: Removing the connection tubing from the screw cap leads to immediate venting of the collection bottle. In systems without quick couplings, removal of the tubing at the pump inlet will also vent the system.

3.4 Filtration

Connect the suction flask to the connection of the hand controller at the bottle head.

NOTICE

The pump achieves its **pumping speed, ultimate vacuum** and vapor pumping rate only at operating temperature (after approx. 15 minutes). Prevent internal condensation, transfer of liquids or dust. The diaphragm and valves will be damaged if liquids are transferred through the pump in significant amounts.

In case of excess temperature the motor is shut down by a **thermal cutout** in the winding. Manual reset is necessary. Switch off or unplug the pump. Identify cause of failure and eliminate. Wait approx. five minutes before restarting the pump. After a power failure, the system will automatically start again.

3.5 Storage following use

NOTICE

Short-term:

Empty bottle.

Has the pump been exposed to **condensate**?

Set the hand controller to continuous aspiration and the Aspirator to maximum underpressure and allow the pump to continue to run for a few minutes. Has the pump been exposed to media which may damage the pump materials or form deposits? Check and clean pump head if necessary.

Long-term:

Take measures as described in short-term shutdown section. Close inlet and outlet port (e.g., with transport caps).

Store the pump in dry conditions.

3.6 Filter and collection bottle

NOTICE

The hydrophobic filter retains water and aqueous solutions also in the form of aerosols. When using solvents, or due to evaporation of water, formation of condensate at the outlet of the Aspirator is possible.

WARNING

- › **Important:** Comply with applicable regulations when disposing of solvents.
- › **Attention:** Do not remove connections in case of liquid in the tube. Do not touch contaminated parts when removing the screw cap. Risk of infection! Wear appropriate protective clothing, e.g., gloves.

Filter


Change filter if necessary in case of contamination, weak aspiration or no aspiration.

- › Use only the original spare part (see section "Accessories - Spare parts") and sterilize and dispose of the filter according to appropriate regulations.

Note: After disinfection with alcohol, allow the filter to dry completely because the filter may become clogged by alcohol.

Emptying the collection bottle of the Aspirator


Check liquid level in the collection bottle regularly.

Maximum admissible liquid level in collection bottle: approx. 80 %, depending on the application. (In case of low boiling liquids or in case of suction of liquids which tend to foam, the admissible liquid level might be reduced.)


Switch off the Aspirator.

Vent the collection bottle.

Remove the connecting tube at the pump inlet, isolate connections in case.

Remove screw cap from the collection bottle after venting.

Remove bottle from the support.

Sterilize and dispose of collected liquid according to appropriate regulations.

4 Cleaning and decontaminating

Over time, discoloration and material changes (e.g., resiliency, elasticity / tightness, cracking) due to repeated steam sterilizations / autoclavings and / or chemical disinfections may occur. Check all parts regularly. Replace defective parts.

121°C
2 bar absolute


Autoclaving

The collection bottle with bottle head and screw cap, the quick coupling and the filter are designated for steam sterilization at 121 °C and 2 bar absolute (1 bar overpressure). Time of exposure according to DIN 58946 $t_e = 20$ minutes.

NOTICE

Prior to autoclaving loosen or remove the bottle head from the bottle.

It is the user's responsibility to ensure effective autoclaving.

max.
80%

NOTICE

Autoclaving

Maximum liquid level in collection bottle: 80 %. When working with low boiling liquids or liquids that tend to foam, the maximum liquid level may be less.

121°C
2 bar absolute
max. 20 x

The number of autoclaving cycles can be marked on the plastic disc (Memory Disc) of the filter (**max. 20 autoclaving cycles** according to the manufacturer instruction).

NOTICE

UV disinfection is permitted but may lead to discoloration of the plastic parts.

Attention: UV disinfection acts only at the surface.

NOTICE

Over time, discoloration and material changes (e.g. resiliency, elasticity/tightness, cracking) due to repeated steam sterilizations may occur.

Important notes on use of disinfectants, see also "Safety during operation"

Aspirator with plastic bottle
Sodium hypochlorite (chlorine bleach)

Aggressive disinfectants which release chlorine or oxygen radicals, e.g. sodium hypochlorite (chlorine bleach) or peroxide compounds can corrode the material of the 4l polypropylene collection bottle and other components (e.g. couplings). This may lead to stress cracks and breakage of the 4l polypropylene bottle.

Aspirator with glass bottle, optional
Sodium hypochlorite (chlorine bleach)

- > Therefore, prolonged use of these disinfectants in the 4l polypropylene bottle is strongly discouraged. After brief disinfecting operations in the bottle, rinse bottle thoroughly to avoid leaving residues of disinfectant in the bottle.
- > Alternative: Use the Aspirator with 2l borosilicate glass bottle.

Attention: The quick-coupling accessory sets between the pumping unit and hand controller are not suitable for use with sodium hypochlorite (chlorine bleach).

- > Using corrosive disinfectants may result in damage, malfunction and/or failure of the equipment.
- > Chlorine permeates the hydrophobic filter on top of the collection bottle and may damage the materials of the equipment or the vacuum supply.
- > Leaking of liquid from a damaged collection bottle or suction tube may lead to exposure of personnel and material or to damage/destruction of wetted equipment or laboratory furniture.

Sekusept® Plus

The use of the disinfectant **Sekusept® Plus** (Manufacturer: Ecolab GmbH & Co OHG, Düsseldorf, Germany) in extensive in-house testing did not cause any damage to the collection bottle. Hence the disinfectant Sekusept® Plus - even when used in the collection bottle during suction - is preferred to sodium hypochlorite (chlorine bleach). Comply with use and safety instructions of the manufacturer!

- > Even if using disinfecting solutions **after termination of the suction** it is absolutely necessary to assure the compatibility of disinfectant and the parts to be disinfected.
- > For information about the compatibility with the materials of the pumping unit, ask the manufacturer of the disinfectant.
- > The wetted materials of the system are listed in the section "Technical data".

5 Assembling of components

5.1 Replacing the filter

atmospheric pressure

Vent the collection bottle.

Ensure that there is no liquid in the tube to avoid risk of contamination.

Remove connecting tube from the filter. Remove the filter from the piece of tube at the hose nozzle.

Attach new filter.

Observe flow direction. Position filter with the printed side "FACING DOWN" towards the bottle.

Attach the connecting tube.

5.2 Assembling a second hand controller connection set (with / without coupling) or conversion to quick coupling hand controller - bottle

Empty bottle.

Decontaminate equipment if necessary.

Remove filter with connection tubing from bottle head.


Remove screw cap from the collection bottle.

Unscrew closing screw.

Only Aspirator with polypropylene bottle: Remove counter nut inside the screw cap.

Aspirator with polypropylene bottle



Insert hose connection (1a) with seal ring (2) or coupling (1b) with seal ring (2) in the screw cap.

Secure lead through with counter nut (3). Assemble hose (4) under the screw cap to the lead through.

Aspirator with glass bottle (accessory)

Screw hose connection (1a) with seal ring (2) or coupling (1b) with seal ring (2) into the screw cap.

Assemble hose (3) under the screw cap to the lead through.

Screw cap onto bottle.


Assemble filter.

Assemble hand controller.


5.3 Assembling quick coupling bottle - pump unit

Attention: The extension set "Quick coupling bottle - pump unit" is designated for two different extension versions. Therefore, the set consists of parts which may not be necessary in each individual case. Redundant parts are not credited by Starlab.

Assembling to pump inlet


Remove tubing.
Disassemble hose nozzle.
Hold counter piece with open-ended wrench size 19 mm.

Assemble coupling (2) with seals (1, 3) and hose nozzle (4).
Hold counter piece with open-ended wrench size 19 mm.
Connect tubing.

Assembling to an additional bottle

Assemble coupling part (1) with seal (2) and hose nozzle (3).
Connect tubing.

5.4 Quick coupling sets

Quick coupling set: Hand controller to bottle

Quick coupling made of PVDF, with adapter to connect a hand controller to a collection bottle. When disconnected, the collection bottle closes vacuum tight.

Quick coupling set: Bottle to pump unit

Quick coupling made of PVDF, to connect a collection bottle to a Aspirator. When disconnected, the collection bottle closes vacuum tight.

NOTICE

Attention: The quick-coupling accessory sets between the pumping unit and hand controller are not suitable for use with sodium hypochlorite (chlorine bleach).

6 Accessories - spare parts

4l PP bottle with protection filter and inlet pipe **N2400-9010**

2l glass bottle with protection filter and inlet pipe **N2400-9012**

Note: Quick couplings must be ordered separately!

4l PP bottle with screw cap **N2400-9011**

(w/o filter, connection for hose, blind plug)

Quick coupling hand controller to bottle **N2400-9008**

Attention: The quick-coupling accessory sets between the bottle and the pumping unit as well as between the bottle and the hand controller are not suitable for use with sodium hypochlorite (chlorine bleach).

Quick coupling bottle to pump **N2400-9009**

Attention: The quick coupling accessory sets between the bottle and the pumping unit as well as between the bottle and the hand controller are not suitable for use with sodium hypochlorite (chlorine bleach).

Protection filter with connective tubing **N2400-9007**

Hand controller for Aspirator **N2400-9001**

Tubing for hand controller, 2.5 m **N2400-9005**

Connection Kit for 2nd hand controller **N2400-9006**

(to be assembled at a Starlab Aspirator, does not include hand controller or quick coupling)

For further accessories and spare parts see instructions for use of the hand controller.

Adapter for 2l glass bottle **N2400-9013**

Mobile underframe **N2400-9014**


Fault	Possible cause	Remedy
<ul style="list-style-type: none"> ⊗ Pump fails to start, power switch does not light up. System becomes vented. 	<ul style="list-style-type: none"> > Not plugged in, or failure of electrical supply? > Power switched off? > Internal fuse blown? 	<ul style="list-style-type: none"> ☑ Plug in, or check fuse controlling power outlet. ☑ Switch on power. ☑ Identify cause of failure, replace fuse (under the cover).
<ul style="list-style-type: none"> ⊗ Pump not running, power switch on. 	<ul style="list-style-type: none"> > Motor overloaded? > Red LED glows continuously? > Filter clogged? 	<ul style="list-style-type: none"> ☑ Switch off, allow motor to cool, identify cause of failure and eliminate before switching on again. ☑ Eliminate leaks. ☑ Replace filter.
<ul style="list-style-type: none"> ⊗ Pump not running, red LED keeps coming on. 	<ul style="list-style-type: none"> > Leak at the system? 	<ul style="list-style-type: none"> ☑ Eliminate leaks.
<ul style="list-style-type: none"> ⊗ No aspiration at the hand controller, pipette drips (above causes excluded). 	<ul style="list-style-type: none"> > Adapter or aspiration tube defective? > Filter clogged? > Leak at the bottle head? > Suction hose, pipette or adapter clogged? > Pump operating but no pumping speed? > Aspiration power setting is off (atmosphere, no vacuum)? 	<ul style="list-style-type: none"> ☑ Replace adapter with aspiration tube. ☑ Replace filter. ☑ Check seals, coupling and blind cap. Screw in or replace if necessary. ☑ Remove clog. ☑ Perform maintenance, replace diaphragm and valves if necessary. ☑ Set aspiration power to "+" (more vacuum).
<ul style="list-style-type: none"> ⊗ Pump switches frequently. 	<ul style="list-style-type: none"> > Leak in the system? > Filter clogged? > Diaphragms or valves of the pump defective? 	<ul style="list-style-type: none"> ☑ Check tubing, seals, coupling, and blind cap. Screw in or replace if necessary. ☑ Replace filter. ☑ Perform maintenance, replace diaphragm and valves if necessary.
<ul style="list-style-type: none"> ⊗ Pump is running continuously. 	<ul style="list-style-type: none"> > Leak in the system? > Too high underpressure for altitude / weather conditions? Vapor pressure of the media inside the collection bottle too high? > Other causes? 	<ul style="list-style-type: none"> ☑ Check tubing, seals, coupling, and blind cap. Screw in or replace if necessary. ☑ Reduce underpressure / aspiration power. ☑ Contact your local distributor.
<ul style="list-style-type: none"> ⊗ Pump is running continuously, all LEDs indicating the suction power are flashing. 	<ul style="list-style-type: none"> > Pressure sensor defective? 	<ul style="list-style-type: none"> ☑ Reduce suction power using key "-". If the error persists, contact your local distributor.
<ul style="list-style-type: none"> ⊗ Keys do not react. 	<ul style="list-style-type: none"> > Keys of the touch panel touched long-lasting? 	<ul style="list-style-type: none"> ☑ Wait 20 seconds and try again.

7.1 Replacing the fuse

DANGER

- > **Switch off the pump.**
- > **Disconnect the electrical power cord** before opening the terminal box.
After disconnecting from power, **wait 5 seconds** to allow the capacitors to discharge.
- > Identify and eliminate the cause of failure before switching on the equipment again.

- > Keep the snap-fit squeezed and pull the fuse holder out.
- > The fuse holder contains two fuses of the same type. Replace the defective fuse with a fuse of the same type (see "Technical data")
- > Insert the fuse holder into the housing of the pumping unit until it snaps into place.

8 Repair - Maintenance - Return

IMPORTANT

In case of service requests contact your local Starlab sales representative or Starlab service technician. Contact details can be found "Starlab Service Hotline" or the back page.

Return Form

Return code

Kindly enclose this form with every return!

Recipient	
STARLAB INTERNATIONAL GmbH Neuer Höltigbaum 38 D-22143 Hamburg	
<input type="checkbox"/> Warehouse	<input type="checkbox"/> Technical Service

Sender	
Company:	
Name:	
Address:	
Customer N°:	
Tel. / E-Mail:	

Reason for return (Please tick)		
<input type="checkbox"/> Wrong delivery	<input type="checkbox"/> Product Defect	<input type="checkbox"/> Wrong item ordered

Delivery note N°:	
-------------------	--

Details :

Returned Items			internal, filled in by the recipient	
Item N°	Lot N° / Serial N°	Quantity	Complete Quantity	Saleable Units
1.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
2.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
3.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
4.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
5.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
6.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
7.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:
8.			<input type="checkbox"/> yes <input type="checkbox"/> no	Units:

Return handling	
Exchange item <input type="checkbox"/>	Credit Note <input type="checkbox"/>

Number of boxes	
Total number of returned boxes:	

Was a decontamination of the returned goods necessary? Yes No

If yes, which was the nature of the contamination?

.....

Method of decontamination:

I confirm that the information on this form is accurate and the items detailed are safe to handle.

Date _____ Signature _____

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg, Germany
E-Mail: info@starlab.de
www.starlabgroup.com
Tel.: +49 (0)40 675 99 39 0
Fax: +49 (0)40 675 9939 20

Authorised Managing Director
Klaus Ambos

Registry Court
Amtsgericht [first instance cou
Hamburg
Register number: HRB 122020
VAT ID no. pursuant to
§ 27 a German
VAT Act DE 214 891 717

User Manual	3–22
Bedienungsanleitung	23–42
Mode d'emploi	43–61

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Starlab (UK), Ltd
Starlab S.r.l.