
# This is your **ErgoOne® FAST** User Manual.

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Latest and up-to-date version of the declaration of conformity is available at:  
[www.starlab.click/doc-ergoonefast](http://www.starlab.click/doc-ergoonefast)

<b>5</b>	<b>Operation</b>	<b>15</b>
5.1	Charging the rechargeable battery	15
5.1.1	Preserving the battery capacity	16
5.1.2	Long periods without operation – with charging stand	16
5.1.3	Long periods without operation – without charging stand	16
5.2	Inserting the pipette	17
5.3	Pipette types	18
5.4	Speed control	20
5.5	Aspirating liquid	21
5.6	Dispensing liquid	22
5.6.1	Flow-out	22
5.6.2	Blow-out	23
5.7	Using the wall mount	24
5.7.1	Mounting the wall mount	24
5.7.2	Removing the wall mount	24
5.8	Using the benchtop support	24
<b>6</b>	<b>Maintenance</b>	<b>25</b>
6.1	Disassembling the pipette clamp	25
6.2	Cleaning	25
6.2.1	Cleaning the Pipette Controller	25
6.2.2	Cleaning the pipette clamp	26
6.3	Mounting the pipette clamp	27
6.4	Checking the leak tightness	27
6.5	Replacing the rechargeable battery	27
<b>7</b>	<b>Troubleshooting</b>	<b>28</b>
7.1	General errors	28
<b>8</b>	<b>Technical data</b>	<b>29</b>
8.1	Weight/dimensions	29
8.2	Mains/power supply device	29
8.3	Rechargeable battery	29
8.4	Ambient conditions	29
<b>9</b>	<b>Transport, storage and disposal</b>	<b>30</b>
9.1	Decontamination before shipment	30
9.2	Transport	30
9.3	Storage	30
9.4	Disposal	31
<b>10</b>	<b>Ordering Information</b>	<b>34</b>

# 1 Operating instructions

## 1.1 Using this manual

- › Read this user manual completely before using the device for the first time. Also observe the instructions for use of the accessories.
- › This user manual is part of the product and must always be easily accessible.
- › Enclose this user manual when transferring the device to third parties.

## 1.2 Danger symbols and danger levels

### 1.2.1 Hazard icons

**Biohazard**

**Explosion**

**Cuts**

**Toxic substances**

**Hazard point**

**Material damage**

### 1.2.2 Degrees of danger

The safety instructions in this user manual indicate the following degrees of danger:

#### **WARNING**

May lead to severe injuries or death.

#### **CAUTION**

May lead to light to moderate injuries.

#### **NOTICE**

May lead to material damage.

# 1.3 Symbols used

Depiction	Meaning
>	You are requested to perform an action.
1. 2.	Perform these actions in the sequence described.
•	List.
<b>i</b>	References useful information.

# 1.4 Glossary

## E

Elapse time

Time required to empty a pipette from the upper part of the scale to the lowest scale graduation or to fully empty a measuring pipette or volumetric pipette.

## 2 Product description

### 2.1 Main illustration


Fig. 1: ErgoOne FAST®

- |   |                                      |    |   |
|---|--------------------------------------|----|---|
| 1 | Rechargeable battery status display  | 8  | Membrane filter                             |
| 2 | Aspirating button                    | 9  | Seal for filter adapter                     |
| 3 | Dispensing button                    | 10 | Filter adapter                              |
| 4 | Rechargeable battery compartment lid | 11 | Pressure compensation opening               |
| 5 | Connector socket                     | 12 | Pipette clamp                               |
| 6 | Aspirating cone                      | 13 | Benchtop support (P7166-6306; not included) |
| 7 | Pipette adapter                      |    |   |

## 2.2 Delivery package

Quantity	Description
1	ErgoOne® FAST
1	Wall mount
1	Sticky tape for wall mount
1	Membrane filter, 0.2 µm, non-sterile, PTFE
2	Membrane filter, 0.45 µm, non-sterile, PTFE (1 filter is built-in)
1	Seal for filter adapter
1	Filter adapter (Silicone)
1	Aspirating cone
1	Rechargeable battery compartment lid
1	Lithium Polymer rechargeable battery, 3.7 V
1	Power supply, 100-240V
1	User manual
1	Quickstart manual

## 2.3 Warranty

For warranty claims, please contact your local supplier. If the housing of the pipette controller is opened by unauthorized individuals, or the device is misused, no warranty claim may be made. The rechargeable battery and all other wear parts are excluded from the warranty.

## 2.4 Materials

**NOTICE! Aggressive substances may damage the ErgoOne® FAST and accessories.** Check the chemical resistance before using organic solvents or aggressive chemicals. Observe the cleaning instructions.

The pipette controller assemblies are composed of the following materials:

Component	Material
Housing, aspirating cone, aspiration button, dispensing button, membrane filter housing, wall holder	Polypropylene (PP)
Filter adapter	Polytetrafluoroethylene (PTFE)
Pipette adapter	Silicone
Filter membrane	Polytetrafluoroethylen (PTFE)
Seal for filter adapter	Hydrated acrylonitrile butadiene rubber (HNBR)
Tubes and valves	Polymer fluorine rubber (FKM), polybutylene terephthalate (PBT), polyphenylene sulfide (PPS), silicone
Battery status display	Cyclic olefin copolymer (COC)

The ErgoOne® FAST is resistant against UV-light. UV-light can therefore be used for decontamination of the pipette controller.

## 2.5 Features

The ErgoOne® FAST is a battery-driven pipette controller. You can use glass or plastic pipettes in a volume range of 0.1 ml to 100 ml.

A pump generates underpressure or overpressure to aspirate or dispense the liquid. The liquid can also be dispensed solely via the atmospheric pressure. The aspirating and dispensing speed is controlled by how far the control buttons are pressed in.

# 3 Safety

## 3.1 Intended use

The pipette controller is intended for dispensing liquids. In-vivo applications (in or on the human body) are not allowed. The pipette controller may only be operated by skilled personnel who have received the appropriate training. All users must have read the user manual carefully and must have become familiar with the device's mode of operation.

## 3.2 Warnings for intended use

**WARNING! Damage to health due to infectious liquids and pathogenic germs.**

- › When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets, and the manufacturer's application notes.
- › Wear personal protective equipment.
- › Follow the instructions regarding hygiene, cleaning and decontamination.
- › For complete instructions regarding the handling of germs or biological material in risk group II or higher, please refer to the "Laboratory Biosafety Manual" (source: World Health Organization, current edition of the Laboratory Biosafety Manual).

**WARNING! Risk of explosion from potentially explosive atmospheres and potentially explosive substances.**

- › Do not use the ErgoOne® FAST in potentially explosive atmospheres.
- › Do not operate the ErgoOne® FAST in rooms where work is conducted using potentially explosive substances.
- › Do not dispense any explosive, highly flammable (flash point < 21°C), extremely flammable (flash point < 0°C), or highly reactive substances with the ErgoOne® FAST.
- › Do not use the ErgoOne® FAST to dispense any substances that could create an explosive atmosphere.

**WARNING! Damage to health due to toxic, radioactive or aggressive chemicals.**

Wear personal protective equipment.

Observe the national regulations for handling these substances.

Observe the material safety data sheets and manufacturer's application notes.

**CAUTION! Poor safety due to incorrect accessories and spare parts.**

The use of accessories and spare parts other than those recommended by us may impair the safety, functioning and precision of the device. We cannot be held liable or accept any liability for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.

- › Only use accessories and original spare parts recommended by us.

**CAUTION! Danger to people due to grossly negligent use.**

- › Never direct the opening of a ErgoOne® FAST at yourself or other people.
- › Only initiate dispensing if it is safe to do so.
- › For all dispensing tasks, make sure that you do not place yourself or other people at risk.

**NOTICE! Damage to device due to penetration of liquids.**

- › The pipette controller itself must not come into contact with the liquid.
- › Only immerse the serological pipette into the liquid.
- › Do not put the pipette controller down when the pipette contains liquid.

**NOTICE! Damage to device from missing pipette.**

- › Use ErgoOne® FAST only when the pipette is inserted.

**Use only ErgoOne® FAST with inserted pipette.**

If you would like to dispense highly flammable or extremely flammable liquids, we recommend using a positive displacement system consisting of manual dispenser (e.g. RepeatOne®) and dispenser tips (e.g. TipOne® REPEAT). Check the chemical resistance and observe the safety notes before using the manual dispenser.

## 3.3 Information on product liability

In the following cases, the designated protection of the device may be compromised. Liability for any resulting property damage or personal injury is then transferred to the operator:


- The device is not used in accordance with the user manual.
- The device is used outside of its intended use.
- The device is used with accessories or consumables which are not recommended by us.
- The device is maintained or repaired by people not authorized by us.
- The user makes unauthorized changes to the device.

# 4 Installation



## 4.1 Preparing installation

- > Use the delivery package details to check that the delivery is complete.
- > Check all parts for any transport damage.
- > Keep the transport carton and the packing material for subsequent safe transport or storage.

## 4.2 Power supply assembly

- > Insert the appropriate power plug adapter into the opening of the power supply. If there is any doubt as to which power plug adapter is required, you should consult an electrician.




## 4.3 Removing the discharge protector

ErgoOne® FAST is delivered with a rechargeable battery inside.

### **WARNING! Injury caused by incorrect handling of the rechargeable battery.**

- > Only use rechargeable original batteries supplied by the ErgoOne® FAST manufacturer.
- > Never pierce, crush or throw the rechargeable battery.
- > Only use the rechargeable battery in the supplied device.
- > Do not touch a leaking rechargeable battery.
- > Do not use a damaged rechargeable battery.
- > Dispose of rechargeable batteries in accordance with the legal requirements.


Prior to first usage please proceed as follows:

1. Slide the battery compartment lid open.
2. Remove the battery and discharge protection.
3. Install the battery.
4. Close the battery compartment lid.


## 4.4 Battery status display during operation

### 4.4.1 Ready for operation


#### **Solid green light:**

The rechargeable battery is fully charged.

#### **Solid yellow light:**


The rechargeable battery is half-charged.

#### **Solid red light:**

The rechargeable battery is almost empty.


## 4.4.2 Empty battery

### **Flashing red light:**


The rechargeable battery needs to be charged.

## 4.4.3 Charge battery


### **Flashing red light:**

The rechargeable battery is <5% charged.

### **Flashing yellow light:**


The rechargeable battery is 5-60% charged.

### **Flashing green light:**

The rechargeable battery is 60-99% charged.

## 4.4.4 Battery is fully charged

The rechargeable battery is >99% charged.

If the power supply is connected when the battery is fully recharged, the battery status display lights up for approx. 30 seconds and no further charging will take place.

# 5 Operation

## 5.1 Charging the rechargeable battery

**WARNING! Incorrect or damaged power supplies can lead to severe personal injury or damage to the device.**

Incorrect or damaged power supplies can cause electric shock or cause the device to overheat, catch fire, melt, short circuit, or cause similar damage.

- › Use only supplied the power supply to charge the device. The correct power supply can be identified by the ‚ErgoOne FAST‘ name on the unit.
- › Do not use a damaged Power supply.

**WARNING! Personal injury caused by incorrect handling of the rechargeable battery.**

- › Do not disassemble or modify the rechargeable battery.
- › Never pierce, crush or throw the rechargeable battery.
- › Only use the battery in the ErgoOne® FAST.
- › Do not touch a leaking rechargeable battery.
- › Do not use a damaged rechargeable battery.
- › Dispose of a damaged rechargeable battery in accordance with the legal requirements.

**NOTICE! Loss of full battery charging capacity if charged incorrectly.**

- › The supplied rechargeable battery is not fully charged. The rechargeable battery only achieves its full capacity after several discharge/charge cycles.
- › Do not charge the rechargeable battery in a hot environment (> 60 °C).
- › Only the power supply provided with the device may be used to charge the rechargeable battery.

## Proceed as follows:


Fig. 2: Loading the rechargeable battery

1. Connect the power supply to the mains power socket.
2. Insert the charging plug of the power supply into the connector socket.

- i** The charging time depends on the charging state of the rechargeable battery. For a fully discharged rechargeable battery, this is about 3 hours.
- i** You can continue using the pipette controller during the charging process.

3 h

### 5.1.1 Preserving the battery capacity

The battery capacity can be preserved over its service life to a great extent.

### 5.1.2 Long periods without operation – with charging stand

- › The pipette controller can be stored, connected to the charging stand. The battery charge will be automatically monitored and recharged.

### 5.1.3 Long periods without operation – without charging stand

1. Completely charge the battery beforehand if the pipette controller is not to be used for a period longer than 4 weeks.
2. The battery should be completely recharged every 2 months.

## 5.2 Inserting the pipette

### **WARNING! Risk of cuts from shattered glass pipettes.**

Glass pipettes are fragile and may cause severe cuts if they break.

- › Do not use force to insert glass pipettes.
- › Wear your user protection (UP).
- › Use a towel to protect the hand used for insertion.

- › Pick up the pipette from above and carefully insert it in the aspirating cone until is positioned securely and air-tight.


Fig. 3: Inserting the pipette


## 5.3 Pipette types

Measuring pipettes and volumetric pipettes are divided into 3 classes: A, AS and B.

Class A and AS pipettes are more precise than Class B pipettes. They vary according to the elapse time, which depends on the nominal volume and the design. Class AS pipettes are quick-drain pipettes.

Volumetric pipettes have 1 or 2 marks and are adjusted to flow-out. Measuring pipettes have a scale and are divided into 4 types.


### Type 1 measuring pipette

- Adjusted to flow-out.
- Nominal volume is indicated by the lowest scale graduation.
- Aspiration up to the zero line at the upper end of the pipette.
- Dispensing to any scale graduation on the scale.


## Type 2 measuring pipette

- Adjusted to flow-out.
- Nominal volume is indicated by the top scale graduation.
- Aspiration to any scale graduation on the scale.
- Dispensing until the pipette is completely empty.

## Measuring pipette type 3 and type 4

### Type 3 measuring pipette

- Adjusted to flow-out.
- Nominal volume is indicated by the tip of the pipette.
- Aspiration up to the zero line at the upper end of the pipette.
- Dispensing to any scale graduation on the scale or until the pipette is completely empty.

### Type 4 measuring pipette

- Adjusted to blow-out.
- Nominal volume is indicated by the tip of the pipette.
- Aspiration up to the zero line at the upper end of the pipette.
- Dispensing to any scale graduation on the scale or until the pipette is completely empty.
- Dispensing the last drops using blow-out.
- Blow-out pipettes are only available in accuracy class B.

## 5.4 Speed control

The liquid aspiration speed is regulated continuously by controlling the pressure on the control buttons, i.e. how far in the button is pressed.

### Slowly aspirating or dispensing liquid


Fig. 4: Slow aspiration

To slowly aspirate or dispense liquid, press the corresponding control button lightly.

### Quickly aspirating or dispensing liquid


Fig. 5: Quick aspiration

To quickly aspirate or dispense liquid, press the corresponding control button firmly.

## 5.5 Aspirating liquid

The speed of aspirating liquid is adjusted steplessly. For this the control knob is pressed at different depths.


**NOTICE! Damage to device due to missing or damaged membrane filter.**

- > Do not use the pipette controller if the membrane filter is not inserted.
- > Replace the membrane filter if it is damaged.

Consider the type of pipette used.

1. Immerse the pipette into the liquid.
2. Slowly press the aspirating button and keep it pressed down.  
The further the aspirating button is pressed, the quicker the liquid will be aspirated.
3. Wipe the pipette on the tube inner wall and remove it.

Fig. 6: Aspirating liquid

## 5.6 Dispensing liquid


Consider the type of pipette used.

After liquid dispensing, hold class AS quick-drain pipettes on the tube inner wall for 5 seconds to allow the liquid to drain.

### 5.6.1 Flow-out

A valve will be opened during flow-out. The liquid drains from the pipette as a result of atmospheric pressure.


1. Hold the pipette vertically and place it on the tube inner wall.
2. Press the dispensing button lightly.

Fig. 7: Allowing liquid to flow-out

## 5.6.2 Blow-out

During blow out, the liquid will be dispensed using the pump.

1. Hold the pipette vertically and place it on the tube inner wall.
2. Press the dispensing button.

Fig. 8: Blowing out liquid

## 5.7 Using the wall mount

For storage, the pipette controller can be put in a wall mount. Please remove the serological pipette before using the wall mount.


Fig. 9: ErgoOne® FAST with wall mount

### 5.7.1 Mounting the wall mount

1. Clean the mounting location on the wall and allow it to dry.
2. Remove the protective foil.
3. Press the wall mount firmly against the wall. Allow the adhesive tape to dry for 24 hours.

### 5.7.2 Removing the wall mount

1. Rotate the wall mount to loosen the adhesive tape.
2. Remove the adhesive tape.

## 5.8 Using the benchtop support

The optional benchtop support can also be used for storage.


Fig. 10:  
Attaching the benchtop support

- › Insert the benchtop support in the groove.
- › To remove the benchtop support, press the sides of the benchtop support together.

# 6 Maintenance

## 6.1 Disassembling the pipette clamp

If liquid has entered the pipette clamp, the aspiration capacity may be decreased, or pipette clamp assemblies may be damaged. The pipette clamp must be disassembled in order to clean or replace the assemblies.

1. Turn the aspirating cone counterclockwise and remove it.
2. Remove the pipette adapter and membrane filter from the filter adapter.
3. Remove the membrane filter from the pipette adapter.
4. Use a narrow pointed object to pry the seal out of the filter adapter.

## 6.2 Cleaning

### 6.2.1 Cleaning the Pipette Controller

Special service is not required.

**NOTICE! Damage to the device due to autoclaving.**

- > Do not autoclave the pipette controller.


To clean contaminated surfaces, proceed as follows:

- > Wipe the housing using a damp cloth.
- > Disinfect surfaces using alcohol (ethanol, propanol) or alcohol-containing disinfectants.

## 6.2.2 Cleaning the pipette clamp

The pipette clamp assemblies can be replaced, cleaned or autoclaved as described below (121 °C, 1 bar overpressure for 20 min.)

### Aspirating cone

- Can be wiped using a damp cloth
- Can be disinfected with alcohol (ethanol, propanol) or alcohol-containing disinfectants.
- Repeatedly autoclavable
- Can be replaced

### Pipette adapter

- Can be rinsed with demineralized water
- Repeatedly autoclavable
- Can be replaced

### Membrane filter

- To be disposed of if contaminated
- Cannot be cleaned
- Can be autoclaved once
- Can be replaced

### Sealing



- Can be rinsed with demineralized water
- Repeatedly autoclavable
- Can be replaced

You can also use a membrane filter with a pore size of 0.2 µm.

## 6.3 Mounting the pipette clamp

1. Push the seal with the groove up into the filter adapter.
2. Push the wide opening of the membrane filter into the narrow opening of the pipette adapter.
3. Guide the aspirating cone over the pipette adapter and turn it until it engages.

## 6.4 Checking the leak tightness

1. Insert the pipette.
2. Fill the pipette with water.
3. Hold the pipette vertically.
4. Observe the pipette outlet for approx. 30 seconds.

Do not touch the pipette. Do not press the control buttons.

No water may be allowed to escape.

5. If water escapes, disassemble and carefully reassemble the pipette clamp.

## 6.5 Replacing the rechargeable battery

### Proceed as follows:

1. Slide the battery compartment lid open.
2. Remove the battery.
3. Insert a new battery.

# 7 Troubleshooting

## 7.1 General errors

Symptom/message	Cause	Remedy
Liquid drips out of the pipette.	<ul style="list-style-type: none"><li>• Pipette adapter and/or membrane filter inserted incorrectly.</li><li>• Pipette not inserted far enough.</li><li>• Pipette adapter damaged.</li><li>• Pipette damaged.</li><li>• Seal damaged.</li><li>• Liquid column too high.</li><li>• Speed too low.</li></ul>	<ul style="list-style-type: none"><li>› Remove the pipette adapter and membrane filter and reinsert them.</li><li>› Continue to carefully insert the pipette.</li><li>› Replace pipette adapter.</li><li>› Replace pipette.</li><li>› Replace seal.</li><li>› Slowly press the aspirating button down further.</li></ul>
Aspiration capacity reduced	<ul style="list-style-type: none"><li>• Membrane filter wetted.</li></ul>	<ul style="list-style-type: none"><li>› Replace membrane filter.</li></ul>
Pipette loose.	<ul style="list-style-type: none"><li>• Pipette adapter damaged.</li></ul>	<ul style="list-style-type: none"><li>› Replace pipette adapter</li></ul>
Bubble formation in the pipette during liquid aspiration.	<ul style="list-style-type: none"><li>• Speed too high.</li></ul>	<ul style="list-style-type: none"><li>› Apply less pressure to the aspirating button.</li></ul>
Rechargeable battery will not charge	<ul style="list-style-type: none"><li>• Rechargeable battery is charged.</li></ul>	<ul style="list-style-type: none"><li>› Disconnect the Power supply.</li><li>› Only recharge the battery if the status display is blinking.</li></ul>
No functional period despite fully charged battery.	<ul style="list-style-type: none"><li>• Battery is older than 3 years.</li></ul>	<ul style="list-style-type: none"><li>› Replace battery.</li></ul>
Shortened useful life and frequent loading necessary.	<ul style="list-style-type: none"><li>• Battery capacity is greatly reduced.</li><li>• Battery is older than 3 years.</li></ul>	<ul style="list-style-type: none"><li>› Replace battery.</li></ul>
Battery housing is deformed.	<ul style="list-style-type: none"><li>• Battery is older than 3 years.</li></ul>	<ul style="list-style-type: none"><li>› Replace battery.</li></ul>

## 8 TECHNICAL DATA

### 8.1 Weight/dimensions

---

<b>Weight</b>	134 g (0.295 lb) (without rechargeable battery, without pipette)
---------------	--

---

### 8.2 Mains/power supply device

---

<b>Input voltage</b>	100 V – 240 V AC, $\pm 10\%$
----------------------	------------------------------

---

<b>Frequency</b>	50 Hz – 60 Hz
------------------	---------------

---

<b>Output voltage</b>	5 V
-----------------------	-----

---

<b>Input current</b>	200 mA – 250 mA
----------------------	-----------------

---

<b>Output current</b>	1 A
-----------------------	-----

---

### 8.3 Rechargeable battery

---

Type	Lithium-Polymer
------	-----------------

---

Voltage	3.7 V
---------	-------

---

Capacity	1100 mAh
----------	----------

---

Charging time	~3 h
---------------	------

---

Weight	26 g (0.057 lb)
--------	-----------------

---

Number of dispensings	~2000 (with a 25 ml-pipette)
-----------------------	------------------------------

---

### 8.4 Ambient conditions

---

<b>Ambience</b>	For indoor use only.
-----------------	----------------------

---

<b>Ambient temperature</b>	5 °C – 40 °C
----------------------------	--------------

---

<b>Relative humidity</b>	10 % – 95 %, non-condensing
--------------------------	-----------------------------

---

<b>Atmospheric pressure</b>	79.5 kPa – 106 kPa
-----------------------------	--------------------

---

# 9 Transport, storage and disposal

## 9.1 Decontamination before shipment

Before sending the pipette controller to the authorized Technical Service for repairs you must decontaminate the pipette and fill out a Declaration of Decontamination. Please note the following:

## 9.2 Transport

- > The device may only be transported in its original packaging.

	<b>Air temperature</b>	<b>Relative humidity</b>	<b>Atmospheric pressure</b>
General transport	-25 °C – 60 °C	10 % – 95 %	30 kPa – 106 kPa
Air freight	-40 °C – 45 °C	10 % – 95 %	30 kPa – 106 kPa

## 9.3 Storage

### **NOTICE! Damage to device due to incorrect storage.**

- > Remove the rechargeable battery if you will not be using the device for an extended period (> 2 months).
- > Do not store the device while the pipette is inserted.
- > Select a secure storage location.
- > Do not expose the device to aggressive gases over an extended period.

	<b>Air temperature</b>	<b>Relative humidity</b>	<b>Atmospheric pressure</b>
In transport packaging	-25 °C – 55 °C	10 % – 95 %	70 kPa – 106 kPa
Without transport packaging	-5 °C – 45 °C	10 % – 95 %	70 kPa – 106 kPa

## 9.4 Disposal

In case the product is to be disposed of, the relevant legal regulations are to be observed.

### Information on the disposal of Electrical and Electronic Equipment in the United Kingdom:

Within the United Kingdom, the disposal of Electrical and Electronic Equipment (EEE) is regulated by national regulations based on Waste Electrical and Electronic Equipment Regulations 2013 (as amended) applicable for EEE.

According to this regulation, any devices put on the market after August 13, 2005, in the business-to-business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste.

They are marked with the following symbol to indicate this:

As the disposal regulations may differ from one country to another, please contact your supplier for more information.

### **WARNING! Risk of explosion and fire due to overheated accumulators and batteries.**

- > Do not heat accumulators and batteries to over 60 °C and do not throw them into fires.

### Disposing of accumulators and batteries

Do not dispose of accumulators and batteries as household waste.

Dispose of accumulators and batteries according to the locally applicable legal regulations.

Only dispose of completely discharged batteries and accumulators.

## **Disposal in the United Kingdom - Compliance with regards to UKCA regulatory requirements**

Electrical and electronic equipment (EEE) is regulated to reduce the amount of waste regarding EEE. EEE producers must help to protect the environment and human health.

Reduction is achieved through various measures in terms of recovery, reuse and recycling of products and components.

For the United Kingdom (UK) 'The Waste Electrical and Electronic Equipment Regulations 2013 (as amended)' is the underpinning UK legislation.

Exempted products from this regulations are:

- Items that protect the country's security such as arms, munitions, and items only for military use.
- A piece of equipment that is designed for and installed in another type of equipment – they can only function within that product, for example a built-in satellite navigation system installed into cars, boats or aeroplanes.
- Filament bulbs apart from LED filament bulbs which are not exempt.

Products that can be used for both military and civil purposes, such as laptops or keyboards, are EEE products.

Excluded products from this regulations are:

- Equipment designed to be sent into space.
- Large scale stationary industrial tools.
- Large scale fixed installations.
- Transport for persons or goods, excluding electric 2 wheeled vehicles which are not type-approved.
- Off-road mobile machinery only for professional use.
- Equipment designed only for research and development use and only available via business to business (B2B).
- Implantable medical devices.
- Medical devices that are expected to be infective at end-of-life.

The aforementioned regulations obliges producers of EEE to help protect natural resources and manage waste of EEE in the best way for people and the environment as well.

EEE is defined as follows:

- which is dependent on electric currents or electromagnetic fields to work properly.
- for generating, transferring and measuring these currents and fields.
- designed for use with a voltage rating 1,000 volts or less for alternating current, and 1,500 volts or less for direct current.

‘Dependent on electric currents or electromagnetic fields to work properly’ means that the equipment needs electric currents or electromagnetic fields (not petrol or gas) to fulfil its basic function. So when the electric current is off, the equipment cannot fulfil its basic function.

Where electrical energy is only used for support or control functions, the equipment is not covered by the regulations. Equipment that only needs a spark to start it (electronic ignition) and does not need electricity to fulfil its basic function includes:

- petrol lawn mowers,
- gas stoves.

Products that are not connected to a mains supply may still be EEE. They can be wind-up, battery-powered and solar-powered products. Where a product has several functions and only one needs an electrical current, the product may still be EEE.

### **Information on the disposal of Electrical and Electronic Equipment in the United Kingdom:**

Within the United Kingdom, the disposal of Electrical and Electronic Equipment (EEE) is regulated by national regulations based on Waste Electrical and Electronic Equipment Regulations 2013 (as amended) applicable for EEE.

According to this regulation, any devices put on the market after August 13, 2005, in the business-to-business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste.

They are marked with the following symbol to indicate this:

As the disposal regulations may differ from one country to another, please contact your supplier for more information.

# 10 Ordering information

<b>Art. No. / Art.-Nr.</b>	<b>Description</b>
S7166-0010	ErgoOne® FAST Pipette Controller
P7166-6402	Wall mount
P7166-6701	Sticky tape for wall mount
P7166-6306	Benchtop support (not included)
P7166-6605	Membrane filter, 0.2 µm, sterile, PTFE, pack of 5
P7166-6800	Membrane filter, 0.45 µm, sterile, PTFE, pack of 5
P7166-6904	Seal for filter adapter
P7166-6002	Pipette adapter (silicone)
P7166-6103	Aspirating cone
P7166-6200	Battery compartment lid
P7166-6501	Lithium Polymer Battery, 3.7 V
P7166-6109	Power supply, 100-240 V / 50-60Hz

<b>User Manual</b> .....	<b>3–34</b>
<b>Bedienungsanleitung</b> .....	<b>35–64</b>
<b>Mode d'emploi</b> .....	<b>65–94</b>
<b>Istruzioni d'uso</b> .....	<b>95–124</b>

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© June 2022. All rights reserved. Technical specifications are subject to change without notice. Our general terms and conditions apply. IMEOF-INT-INT-06-2022
**Starlab GmbH**
**Starlab S.r.l.**