# This is your **ErgoOne® E** Multi-Channel Pipettes user manual.

USER MANUAL · BEDIENUNGSANLEITUNG  
MODE D'EMPLOI · ISTRUZIONI PER L'USO · INSTRUCCIONES DE MANEJO

Passionate for science.

<b>User Manual</b>	<b>3–40</b>
<b>Bedienungsanleitung</b>	<b>41–78</b>
<b>Mode d'emploi</b>	<b>79–116</b>
<b>Istruzioni d'uso</b>	<b>117–154</b>
<b>Instrucciones de manejo</b>	<b>155–192</b>
<b>1. Introduction</b>	<b>5</b>
1.1 Scope of delivery	5
1.2 Terms of use	5
1.2.1 Hazard levels	5
1.2.2 Format	5
1.2.3 Operating manual symbols	5
<b>2. Safety regulations</b>	<b>6</b>
2.1 General safety regulations	6
2.2 Intended use	7
2.3 Limitations of use	7
2.4 Application restrictions	7
2.5 Operating exclusions	7
2.6 Rechargeable battery and power adapter specifications	7
<b>3. Functional and operational components</b>	<b>8</b>
<b>4. Operation</b>	<b>9</b>
4.1 First Steps	9
4.2 Set volume	10
4.3 Set the aspiration and dispensing speed	11
4.4 Pipetting	12
4.5 Direct trigger blow-out	13
<b>5. Pipetting programs</b>	<b>14</b>
5.1 PIP mode	14
5.2 PIPmix mode	15
5.3 revPIP mode	16
5.4 Electrophoresis (GEL) mode	18
5.5 DISP mode	19
<b>6. Checking the volume</b>	<b>21</b>
<b>7. Accuracy table</b>	<b>22</b>
<b>8. Adjustment – Easy Calibration</b>	<b>23</b>
8.1 Adjustment	23
8.2 Restore factory setting	24

<b>9.</b>	<b>Disinfection/autoclaving</b>	<b>24</b>
9.1	UV sterilization	24
9.2	Autoclaving	25
9.3	Reference run (rEF)	25
<b>10.</b>	<b>Maintenance</b>	<b>25</b>
10.1	Separate the pipetting unit from the handle	26
10.2	Disassembly of ErgoOne® E Multi-Channel pipettes up to 300 µl	26
10.3	Disassembly of ErgoOne® E Multi-Channel pipettes 1250 µl	29
10.4	Charge and replace rechargeable battery	33
10.5	Rechargeable battery regeneration function	34
<b>11.</b>	<b>Troubleshooting</b>	<b>35</b>
<b>12.</b>	<b>Product markings</b>	<b>36</b>
<b>13.</b>	<b>Order Information</b>	<b>37</b>
13.1	Ordering Information	37
13.2	Spares	37
13.2.1	ErgoOne® E Multi-channel pipettes to 300 µl	37
13.2.2	ErgoOne® E Multi-Channel pipette 1250 µl	38
13.3	Additional accessories	38
<b>14.</b>	<b>Repairs</b>	<b>39</b>
14.1	Sending for repair	39
<b>15.</b>	<b>Calibration service</b>	<b>39</b>
<b>16.</b>	<b>Warranty and guarantee</b>	<b>39</b>
<b>17.</b>	<b>Disposal</b>	<b>40</b>

# 1. Introduction

## 1.1 Scope of delivery

The ErgoOne® E Multi-Channel pipette, one battery pack, one reagent reservoir, grease, operating manual.

Nominal volume	Scope of delivery
Up to 300 µl	Plug-in power supply unit with battery charging cable, 1 instrument stand, 1 set of O-rings made of FKM/silicone, 1 TipOne® pipette tip rack
1250 µl	Universal power adapter with battery charging cable, 1 shelf mount, 1 set of O-rings made of silicone, 2 TipOne® pipette tip racks

## 1.2 Terms of use

- › Carefully read the operating manual before using the device for the first time.
- › The operating manual is part of the device and must be kept in an easily accessible place.
- › Be sure to include the operating manual if you transfer possession of this device to a third party.

### 1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

### 1.2.2 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

### 1.2.3 Operating manual symbols

Icon	Meaning	Icon	Meaning
	Hazardous, may cause moderate injuries.		Explosion hazard

## 2. Safety regulations

### 2.1 General safety regulations

#### Please read carefully!

The instrument ErgoOne® E Multi-Channel pipette can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual before operation.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.  
When working with infectious or hazardous samples, standard laboratory rules and precautions must be adhered to.
3. Observe all specifications provided by the reagent manufacturer.
4. Do not operate the instrument in potentially explosive atmospheres and do not pipette highly flammable media.
5. Use the instrument only for pipetting liquids within the defined limitations and restrictions of use. Comply with the operating exclusions; see Operating exclusions.  
If in doubt, contact the manufacturer or supplier.
6. Always perform work in a way that does not endanger the user or any other person. Avoid splashes. Use only suitable vessels.
7. Avoid touching the tip opening when working with aggressive media.
8. Never use force.
9. Use only original spare parts. Do not attempt to make any technical modifications.  
Do not disassemble the instrument any further than is described in the operating manual.
10. Always check that the instrument is in proper working condition before use. If there are any signs of the instrument malfunctioning (e.g., sluggish piston, leakage), stop pipetting immediately, and follow the instructions in the chapter Troubleshooting.  
Contact the manufacturer, if necessary.
11. The original rechargeable battery may not be exchanged for batteries or rechargeable batteries from other manufacturers.
12. Use only the original power adapter to charge the nickel-metal hydride battery.
13. Protect the power adapter from moisture and use it only with this instrument.
14. Dispose of batteries according to the applicable regulations.

#### WARNING!

#### Potential risk of explosion due to damaged battery

In extreme cases, improper handling of the instrument or the battery (short circuit, mechanical damage, overheating, etc.) may cause the battery to explode.

### 2.2 Intended use

ErgoOne® E Multi-Channel is a microprocessor-controlled, battery-operated piston-stroke pipette based on the air cushion principle for pipetting aqueous solutions of medium density and viscosity. If the instrument is handled correctly, the sample to be dispensed comes into contact only with the tip and not with the ErgoOne® E Multi-Channel.

## 2.3 Limitations of use

This instrument is intended for pipetting samples, within the following limitations:

- › Operating temperature of instrument and reagent should be between +15 °C and +40 °C (59 °F to 104 °F) (other temperatures upon request)
- › Vapor pressure up to 500 mbar
- › Viscosity: 260 mPas

For viscous media, the speed must be adjusted if necessary.

## 2.4 Application restrictions

Viscous and wetting liquids may compromise volumetric accuracy. Volumetric accuracy may also be affected when pipetting liquids whose temperature deviates from the ambient temperature by more than  $\pm 1$  °C/ $\pm 1.8$  °F.

## 2.5 Operating exclusions

The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:

**Liquids with very high steam pressure**

**Liquids that corrode the following materials:**

- › Fluoroelastomer rubber (FKM),
- › Polycarbonate (PC)
- › Polycarbonate/polybutylene terephthalate (PC/PBT)
- › Polyetherimide (PEI)
- › Polypropylene (PP)
- › Polyvinylidene fluoride (PVDF)
- › Silicone

### NOTICE!

#### Autoclaving

The hand grip is not autoclavable!

a) Heat and steam during autoclaving will damage the electronics.

## 2.6 Rechargeable battery and power adapter specifications


#### Rechargeable battery

Nickel metal hydride rechargeable battery with 3 cylindrical single cells size AAA, 3.6 V, 700 mAh

#### Power adapter

Output voltage 6.5 V DC, 200 mA

### 3. Functional and operational components

\*) The instrument is switched on by pressing the 'Enter' key. By subsequently pressing the pipetting key, the instrument is ready for pipetting.

The ErgoOne® E Multi-Channel pipette switches itself off automatically 10 minutes after the last operation (auto power off). The ErgoOne® E Multi-Channel pipette fits ergonomically in your hand. For absolutely effortless operation of the function keys, the hand position can be further optimized by means of the finger grip, which is height-adjustable via a screw.

## 4. Operation


### 4.1 First Steps

#### 1. Insert rechargeable battery

- a. Open the cover of the battery compartment.
- b. Insert the rechargeable battery into the compartment. Make sure that the plug of the rechargeable battery is firmly inserted into the socket in the instrument.
- c. Replace the cover of the battery compartment and close it.


#### 2. Activate instrument

The ErgoOne® E Multi-Channel pipette automatically requests a reference run immediately after the battery is inserted. After pressing the pipetting key, the reference run is carried out, and the instrument is ready for pipetting.

Rechargeable  
battery capacity  
indicator

Pipetting

The display shows the factory-set default pipetting mode (PIP) and the respective nominal volume. The aspiration and dispensing speed are set to maximum at the factory.

Arrow icon  
'Aspirate'

Volume  
display


The simple volume and speed setting is described on the following pages.

## 4.2 Set volume


The volume is factory-set to the respective nominal volume of the ErgoOne® E Multi-Channel pipette and can be easily and quickly changed individually.


- a. Pressing one of the arrow keys directly selects a volume.  
⇒ 'VOL' flashes.

- b. Pressing the arrow key (-) decreases the volume. Pressing and holding the arrow key causes the volume to change quickly.  
⇒ 'VOL' continues to flash.

- c. Pressing the arrow key (+) increases the volume. Pressing and holding the arrow key causes the volume to change quickly.  
⇒ 'VOL' continues to flash.

- d. To confirm the volume selection, press the 'Enter' key.  
⇒ The display now shows the newly set volume. Here, for example, the display of the PIP mode set by default.


### NOTICE!

By pressing the menu key, any setting procedure can be canceled. The display then jumps to the next setting option or back to the initial display.

### 4.3 Set the aspiration and dispensing speed

Aspiration and dispensing speed are separately adjustable. When the menu is called up, the last speed set in each case is displayed. Five speed levels are available in each case.


#### Set aspiration speed


1x

- a. Pressing the menu key once briefly takes you to the aspiration speed menu.  
⇒ 'Speed' flashes.

- b. Press the arrow keys (+/-) to select the speed level (e.g. level 5).  
⇒ 'Speed' continues to flash.

- c. Press 'Enter' key.  
⇒ The display returns to the basic state of the respective set mode. Here, for example the display of the standard PIP mode.


#### Set dispensing speed

2x

- a. Pressing the menu key twice briefly takes you to the dispensing speed menu.  
⇒ 'Speed' flashes.

- b. Press the arrow keys (+/-) to select the speed level (e.g. level 2).  
⇒ 'Speed' continues to flash.


- c. Press 'Enter' key.  
The display returns to the basic state of the respective set mode. Here, for example the display of the standard PIP mode.

## 4.4 Pipetting

The volume is factory-set to the respective nominal volume of the ErgoOne® E Multi-Channel pipette and can be easily and quickly changed individually, see Set volume.

### NOTICE!

- › Perfect analysis results can only be achieved by using quality tips. We recommend the TipOne® tips which are delivered with the pipette or other compatible pipette tips, see TipFinder on Starlab's website [www.starlabgroup.com](http://www.starlabgroup.com).
- › If other pipette tips are used, check whether these tips fit the ErgoOne® E Multi-Channel pipette and have the required quality for the desired application.



- Insert tip vertically:**  
Use the correct tips, in accordance with the volume range or color code!  
Make sure that the tips are firmly in place and leak tight. Pipette tips are disposable products!
- Align pipetting unit:** The pipetting unit can turn freely in both directions.
- Aspirate liquid:** Hold the instrument vertically and immerse the tip in the liquid. The liquid is aspirated by pressing the pipetting key. The arrow in the display points upwards (aspiration).

Leave the tip immersed in the liquid for a few seconds, so that the set volume is aspirated completely. This is especially important when pipetting viscous media and when using pipettes with large volumes.

Volume range	Immersion depth	Wait time
0.5 – 100 µl	2 – 3 mm	1 s
100 – 300 µl	2 – 4 mm	1 s
> 1000 µl	3 – 6 mm	3 s

- d. Dispense liquid: When the liquid aspiration is finished, the arrow in the display points downwards (dispensing). Place the pipette tip against the vessel wall. Hold the pipette at an angle of 30–45°. By pressing the pipetting key again, the liquid is dispensed completely with automatic blow-out. While doing this, wipe the pipette tip against the vessel wall.


- e. Ejecting a tip: Hold the pipette Pipetting unit over a suitable disposal container and press the tip ejection key.

**NOTICE!**


ISO 8655 requires that pipette tips are pre-wetted once before the actual pipetting procedure.

**4.5 Direct trigger blow-out**

The blow-out can also be triggered directly at any time if necessary.

- a. Call up the blow-out function: Press 'Enter' key. The display shows 'blo' for blow-out.

- b. Trigger blow-out: Press the pipetting key once to trigger the blow-out. The display will then jump back to the set pipetting mode (start position).

**NOTICE!**

During blow-out, the piston moves completely downwards. Ensure that any residual liquid is discharged safely. **Keeping the pipetting key pressed keeps the piston down and thus prevents accidental aspiration of liquid. Releasing causes the piston to return to the start position.**

## 5. Pipetting programs

Activity	Designation	Info
<b>Normal pipetting</b>	PIP mode, see PIP mode	Standard program. A previously entered volume is aspirated and dispensed again
<b>Pipetting for electrophoresis</b>	GEL mode, see Electrophoresis (GEL) mode	Program for loading electrophoresis gels. A previously defined sample volume is aspirated at a high, variable speed and slowly released again.
<b>Mixing samples</b>	PIPMix mode, see PIPmix mode	Program for mixing liquids. Sample is aspirated and dispensed repeatedly.
<b>Reverse pipetting</b>	revPIP mode, see revPIP mode	Program especially for pipetting liquids with high viscosity, high vapor pressure, or foaming media.
<b>Dispensing</b>	DISP mode, see DISP mode	Program for dispensing liquids. A volume taken up is dispensed again in partial steps.

### NOTICE!


#### GEL mode

The GEL mode is not available for ErgoOne® E 1000 µl and 5000 µl Multi-Channel pipettes.


#### 5.1 PIP mode

The default program – a previously entered volume is aspirated and dispensed again.


Volume and speed setting, see Set volume and Set the aspiration and dispensing speed.

- a. Call up menu selection: Pressing the menu key three times takes you to the program selection.  
⇒ 'Mode' flashes.

- b. Set PIP mode: Use one of the arrow keys to scroll through the modes until 'PIP' appears.  
⇒ 'Mode' continues to flash.


- c. Confirm PIP mode: Press 'Enter' key.
- ⇒ The display now shows 'blo' for blow-out.


- d. Preparing for pipetting: Press the pipetting key once to move the piston to its start position.
- ⇒ The arrow in the display points upwards (aspiration).

- e. Aspirate liquid: Press the pipetting key once to aspirate the liquid.


- f. Dispense liquid: The liquid is dispensed by pressing the pipetting key once.
- ⇒ The arrow in the display points downwards (dispensing).


- g. Trigger blow-out? You do not have to do anything. When pipetting in PIP mode, the blow-out is performed automatically.

## 5.2 PIPmix mode

Program for mixing liquids. Sample is aspirated and dispensed repeatedly. Volume and speed setting, see Set volume and Set the aspiration and dispensing speed.

- a. Call up menu selection: Pressing the menu key three times takes you to the program selection.
- ⇒ 'Mode' flashes.

- b. Set PIPmix mode: Use one of the arrow keys to scroll through the modes until 'PIPmix' appears.
- ⇒ 'Mode' continues to flash.


- c. Confirm PIPmix mode: Press 'Enter' key.
- ⇒ The display now shows 'blo' for blow-out.


- d. Preparing for pipetting: Press the pipetting key once to move the piston to its start position.
- ⇒ The arrow in the display points upwards (aspiration).

- e. Aspirate liquid: Press the pipetting key once to aspirate the liquid.

- f. Dispense liquid in PIPmix mode: Pressing and holding the pipetting key causes the liquid to be alternately dispensed and aspirated. The display alternately shows the arrow symbol for aspiration or dispensing as well as the number of cycles.


- g. End pipetting: Pressing the pipetting key once dispenses the liquid and triggers the blow-out. After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).

**NOTICE!**


The display shows a maximum of 19 cycles.

**5.3 revPIP mode**


Program especially for pipetting liquids with high viscosity, high vapor pressure, or foaming media. Volume and speed setting, see Set volume and Set the aspiration and dispensing speed.

- a. Call up menu selection: Pressing the menu key three times takes you to the program selection.
- ⇒ 'Mode' flashes.

- b.** set revPIP mode: Use one of the arrow keys to scroll through the modes until 'revPIP' appears.  
 ⇨ 'Mode' continues to flash.

- c.** Confirm revPIP mode: Press 'Enter' key.  
 ⇨ The display now shows 'blo' for blow-out.

- d.** Preparing for pipetting: Press the pipetting key once to move the piston to its start position.  
 ⇨ The arrow in the display points upwards (aspiration).

- e.** Aspirate liquid: Press the pipetting key once to aspirate the liquid.


**NOTICE!**


When aspirating the liquid, slightly more volume is aspirated than set.

- f.** Dispense liquid in revPIP mode: Press the pipetting key once to dispense. On the display, the arrow points down (dispensing). The set volume is now dispensed and some liquid remains in the tip.

- g.** Re-aspirate liquid in revPIP mode: Pressing the pipetting key again now resumes the set volume. (Pressing the pipetting key again releases the volume)

- h.** Trigger blow-out: After the last pipetting, press 'Enter' key.  
 ⇨ The display now once again shows 'blo' for blow-out.


- i. End pipetting: By pressing the pipetting key once, the blow-out is triggered, and the residual liquid is dispensed.
- ⇒ After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).

### 5.4 Electrophoresis (GEL) mode

Program for loading electrophoresis gels. A previously defined sample volume is aspirated at a high, variable speed and slowly released again.

Volume and speed setting, see Set volume and Set the aspiration and dispensing speed.

- a. Call up menu selection: Pressing the menu key three times takes you to the program selection.
- ⇒ 'Mode' flashes.

- b. Set GEL mode: Use one of the arrow keys to scroll through the modes until 'GEL' appears.
- 'Mode' continues to flash.

- c. Confirm GEL mode: Press 'Enter' key. The display now shows 'blo' for blow-out.


- d. Preparing for pipetting: Press the pipetting key once to move the piston to its start position.
- ⇒ The arrow in the display points upwards (aspiration).

- e. Aspirate liquid: Press the pipetting key once to aspirate the liquid.

- f. In order aspirate more liquid than set (up to max. 110% of the nominal volume), keep the pipetting key pressed during the aspiration process until the desired volume has been aspirated.
- ⇒ A rhombus appears in the display.

Hold pressed


**g.** Dispense liquid in GEL mode: To dispense, briefly press the pipetting key once. A rhombus appears in the display. The aspirated volume is slowly released again.

**h.** The dispensing of the sample can be interrupted by pressing the pipetting key again.

⇒ The display shows the volume of the liquid dispensed.

**i.** Trigger blow-out: After the last pipetting, press 'Enter' key.

⇒ The display now once again shows 'blo' for blow-out.

**j.** End pipetting: By pressing the pipetting key once, the blow-out is triggered, and the residual liquid is dispensed.

⇒ After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).


**NOTICE!**

The GEL mode requires very slow dispensing speeds in order to prevent sample turbulence. In order to ensure optimum dispensing, the dispensing speed is set at the factory. It is considerably slower than the adjustable level 1 and cannot be selected individually.


**5.5 DISP mode**

Program for dispensing an aspirated liquid in partial steps. Slightly more fluid is aspirated than is mathematically necessary.


Volume and speed setting, see Set volume and Set the aspiration and dispensing speed.

**a.** Call up menu selection: Pressing the menu key three times takes you to the program selection.  
⇒ 'Mode' flashes.

**b.** Set DISP mode: Use one of the arrow keys to scroll through the modes until 'DISP' appears.  
⇒ 'Mode' continues to flash.

- c. Confirm DISP mode: Press 'Enter' key.
- ⇒ The display now shows 'blo' for blow-out.

- d. Preparing for pipetting: Press the pipetting key once to move the piston to its start position.
- ⇒ The arrow in the display points upwards (aspiration).

- e. Set partial volume: Pressing the arrow key (+/-) sets the volume. Pressing and holding the arrow key causes the volume to change quickly.
- ⇒ 'VOL' flashes.

- f. Confirm partial volume: Press 'Enter' key. The display shows the newly set partial volume.
- ⇒ 'Steps' flashes. The maximum possible number of steps is displayed.


- g. Set the number of steps: Pressing the arrow key (+/-) sets the number of steps.
- ⇒ 'Steps' continues to flash.


- h. Confirm number of steps: Press 'Enter' key.
- ⇒ The display shows the set number of steps.

- i. Aspirate liquid: Press the pipetting key once to aspirate the liquid.


- j. Dispense liquid: Each time the pipetting key is pressed, a dispensing step is performed. The arrow in the display points downwards (dispensing). The 'Step' display shows the number of steps remaining.

- k. Trigger blow-out: After the last dispensing step, press 'Enter' key.  
⇒ The display now once again shows 'blo' for blow-out.
- l. End dispensing: By pressing the pipetting key once, the blow-out is triggered, and the residual liquid is dispensed.
- m. After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).

## 6. Checking the volume

We recommend testing the device every 3 to 12 months, depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655, Part 6.

### 1. Setting the nominal volume

- a. Set the maximum specified instrument volume (for procedure, see Pipetting).

### 2. Conditioning the pipette

- b. Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.

### 3. Performing the test

#### NOTICE!

In accordance with DIN EN ISO 8655-2, a tip change is recommended after each individual measurement. An exception to this rule can be made, according to DAkkS guideline DKD-R8-1.

- a. Aspirate the test liquid and pipette into the weighing vessel.

#### NOTICE!

Each individual channel must be inspected separately.

- b. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- c. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- d. At least 10 pipetting series and weighings in 3 volume ranges (100 %, 50 %, 10 %) are recommended.

## Calculation (for nominal volume)

$x_i$  = weighing results

$n$  = Number of weighings

$V_0$  = nominal volume

$Z$  = Correction factor (e.g. 1.0029  $\mu\text{l}/\text{mg}$  at 20 °C, 1013 hPa)

### Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

### Mean volume:

$$\bar{V} = \bar{x} \cdot Z$$

### Accuracy\*:

$$A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

### Coefficient of variation\*:

$$CV\% = \frac{100 s}{\bar{V}}$$

### Standard deviation\*:

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

\*) = Calculation for accuracy (A%) coefficient of variation (CV%): A% and CV% are calculated using the formulas of statistical quality control.


## NOTICE!

Test Instructions (SOPs) are available for download from [www.starlabgroup.com](http://www.starlabgroup.com).

## 7. Accuracy table

Volume range [ $\mu\text{l}$ ]	Partial volume [ $\mu\text{l}$ ]	A* $\leq \pm$ %	CV* $\leq$ %	Substeps [ $\mu\text{l}$ ]	Recommended tip type [ $\mu\text{l}$ ]
0.5 - 10	10	1.2	0.8	0.01	0.5 - 20
	5	2.0	1.5		
	1	8.0	4.0		
1 - 20	20	1.0	0.5	0.02	0.5 - 20
	10	2.0	1.0		
	2	8.0	3.0		
5 - 100	100	0.8	0.25	0.1	2 - 200
	50	1.6	0.4		
	10	4.0	1.5		
15 - 300	300	0.8	0.25	0.5	5 - 300
	150	1.4	0.4		
	30	4.0	1.2		
50 - 1250	1250	1.0	0.25	1.0	50 - 1250
	625	1.2	0.4		
	125	5.0	1.2		

\*A = Accuracy, CV = Coefficient of Variation



Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

## 8. Adjustment – Easy Calibration


### 8.1 Adjustment

The nominal volume or volume to be tested is set, standard mode Pipetting (PIP), e.g. 200.0 µl (for procedure, see PIP mode). Ex: Volume according to volume test 301.3 µl.

- a.** Bring up CAL mode: Press and hold the MENU key (> 3 sec) until CAL appears.
- ⇒ The display reads ,off'.
  - ⇒ 'CAL' blinks.

- b.** Activate CAL mode: Press one of the arrow buttons to activate the CAL mode.
- ⇒ The display changes from ,off' to ,on'.
  - ⇒ 'CAL' continues to blink.

- c.** Confirm CAL mode: Press the Enter key.
- ⇒ The display now shows the set pipetting volume.
  - ⇒ 'CAL' blinks.

- d.** Setting the volume: Use the arrow keys (+/-) to set the volume, which was previously determined and tested.
- ⇒ 'CAL' blinks.


- e.** Confirm volume: Press the Enter key.
- ⇒ The display shows the tested and corrected volume.
  - ⇒ The CAL symbol is continuously displayed to confirm that an adjustment has been made.

## 8.2 Restore factory setting

The CAL symbol permanently shown in the display indicates that an adjustment has been made.

- a. Enter CAL mode: Pressing and holding (> 3 s) the menu key invokes the CAL mode.
  - ⇒ The display shows 'on'.
  - ⇒ 'CAL' flashes.

- b. Switch off CAL mode: Pressing one of the arrow keys deactivates the CAL mode.
  - ⇒ The display changes from 'on' to 'off'.
  - ⇒ 'CAL' continues to flash.

- c. Restore factory condition: Press 'Enter' key.
  - ⇒ The permanently displayed CAL icon has disappeared.
  - ⇒ The instrument has been returned to the original factory condition.

### NOTICE!

For the ErgoOne® E Multi-Channel pipettes, a volume offset is applied during adjustment (i.e. the volume changes by the same amount over the entire volume range of the pipette). It is recommended to perform the adjustment at 50% of the nominal volume.

### NOTICE!

The instrument is permanently adjusted for aqueous solutions but can also be adjusted for solutions of different density, viscosity, and temperature. The ErgoOne® E Multi-Channel pipettes can be adjusted in any mode (except GEL mode).

## 9. Disinfection/autoclaving

### 9.1 UV sterilization

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

## 9.2 Autoclaving

The highlighted part of ErgoOne® E Multi-Channel pipettes can be autoclaved at 121°C (250 °F), 2 bar, and with a holding time of at least 15 minutes according to DIN EN 285.

- a. Eject the pipette tips.
- b. Separate the pipetting unit from the handle, see Maintenance.
- c. Autoclave the complete pipetting unit without further disassembly.
- d. Allow the pipetting unit to cool down completely and dry.
- e. Screw the pipetting unit back into the handle.
- f. Carry out reference run (rEF).

### NOTICE!

The effectiveness of autoclaving must be verified by the user. Maximum safety is achieved through vacuum sterilization. We recommend the use of sterilization bags.

If the pipette is autoclaved frequently, the piston and seal should be greased with the grease supplied in order to ensure proper movement.

## 9.3 Reference run (rEF)

After each change of the pipette shaft, a manual reference run must be carried out. The reference run is used for the secure coupling of the piston.

- a. Call up rEF mode: Pressing the 'Menu' and 'Enter' keys simultaneously activates the rEF mode.
  - ⇒ The display shows 'rEF'.
- b. Carry out reference run: Pressing the pipetting key once triggers the reference run.
  - ⇒ A clear functional noise can be heard.
  - ⇒ After the reference run, the display automatically switches back to the previously set program.

## 10. Maintenance

In order to ensure proper functioning, the ErgoOne® E Multi-Channel pipettes should be serviced at regular intervals and cleaned as necessary.


### What must be checked?

- a. Check the pipette shafts, pistons, and seals for damage and contamination.
- b. Check the device for leaks. Aspirate a sample and hold the instrument vertically for approx. 10 s. If drops form at the pipette tips, see Troubleshooting.

## 10.1 Separate the pipetting unit from the handle

Example pipetting unit up to 100 µl

Example pipetting unit up to 1250 µl

- Eject the pipette tips.
- To separate, pull the pipetting unit from the handle downward as far as possible. **Only then**, turn it clockwise. After one revolution, the unit no longer needs to be pulled downward while turning.
- If the pipetting unit is unscrewed, it must be pulled down again in order to disconnect the magnetic coupling.

### NOTICE!

- › When assembling, the pipetting unit must be screwed counter-clockwise onto the handle so that it clicks into place.
- › When assembling, do not pull the pipetting unit downwards.
- › Improper handling can lead to damage.

## 10.2 Disassembly ErgoOne® E Multi-Channel pipettes up to 300 µl


The three main components of the pipetting unit can be easily separated and disassembled for maintenance, cleaning or replacing parts. The procedure is explained in detail on the following pages.

Changing the O-rings on the individual shafts is described in detail in the instructions included with the replacement part.

### Cleaning instructions


- Clean individual shafts and pistons (only these parts) with soap solution or isopropanol. Then rinse with distilled water.
- Allow the parts to completely dry and cool. Liquid residues in the shafts result in deviations in accuracy.
- Re-grease pistons with a very thin coat of the grease supplied. For the central guide axis (Z), use only the prescribed fluorine static grease.
- After assembling the unit, carry out a reference run (rEF).

## Removing shafts and seals for cleaning or replacement


**b.**

- a. Unscrew pipetting unit: Separate the pipetting unit from the handle.
- b. Remove the pipetting housing: Turn both closures of the pipetting housing cover by 90°, and pull off the pipetting housing.

**c.**

- c. Unscrew the shaft: Place the installation tool onto a single shaft and unscrew the shaft.


**d.**

- d. Remove seal: Push the piston unit all the way down. Remove and inspect the seal and clean or change as necessary.

### NOTICE!

After removing the shaft, the seal is located either in the shaft or on the piston.

**e.**


- e. Mount the seal: If required, re-grease the piston with the silicone grease supplied. Slide the seal onto the piston with the flat side upwards. Tighten the cleaned or new shaft using the installation tool.
- f. Assemble pipetting unit: Reassemble the pipetting unit and mount it on the handle. Check the instrument for leaks, free movement, and precision.

### NOTICE!


The pipetting unit must be screwed onto the handle in a counter-clockwise direction until it clicks into place. Do not pull the pipetting unit downwards.

- g. Carry out reference run: Carry out reference run (rEF).

## Remove piston for cleaning or replacement


- a. Unscrew pipetting unit: Separate the pipetting unit from the handle.
- b. Remove the pipetting housing: Turn both closures of the pipetting housing cover by 90°, and pull off the pipetting housing.


- c. Remove screws: Remove both outer Phillips-head screws on the shaft unit.

### NOTICE!


The central guide axis (Z) must not be removed.


- d. Separate the piston and shaft unit: Pull apart and separate the piston and shaft unit.

- e. Remove piston: Place the installation tool onto the piston nut and unscrew it. Remove piston nut, and pull out piston.

- f.** Mount piston: Insert cleaned or new piston. Tighten the piston nut again using the installation tool.

- g.** Assemble the piston and shaft unit: Loosen the shafts by half a turn. Then place the shaft unit on the piston unit and fasten it. Then screw on the shafts.
- h.** Assemble pipetting unit: Reassemble the pipetting unit and mount it on the handle. Check the instrument for leaks, free movement, and precision.

**NOTICE!**

The pipetting unit must be screwed onto the handle in a counter-clockwise direction until it clicks into place. Do not pull the pipetting unit downwards.


- i.** Carry out reference run: Carry out reference run (rEF).

### 10.3 Disassembly of ErgoOne® E Multi-Channel pipettes 1250 µl

The three main components of the pipetting unit can be easily separated and disassembled for maintenance, cleaning or replacing parts. The procedure is explained in detail on the following pages.

Changing the O-rings on the individual shafts is described in detail in the instructions included with the replacement part.

## Main components of the pipetting unit

### A

Complete **piston shaft system** with central guide axis (Z) connected to the pipetting housing cover [II] (channel identification can be read here).

### B

**Piston–shaft units** inserted in the frame [I] of the pipetting unit. The units consist of the piston (with seal) (B') and the stem (with O-ring) (B''), which can be disassembled for cleaning, lubrication, and replacement.


### C

Pipetting housing connected to the pipetting housing cover [II] by two sliding catches.

## Cleaning instructions

- Clean individual shafts and pistons (only these parts) with soap solution or isopropanol. Then rinse with distilled water.
- Allow the parts to completely dry and cool. Liquid residues in the shafts result in deviations in accuracy.
- Re-grease pistons with a very thin coat of the grease supplied. For the central guide axis (Z), use only the prescribed fluorine static grease.
- After assembling the unit, carry out a reference run (REF).


## Removing shafts and pistons for cleaning or replacement

**b.**

**c.**

**d.**

**e.**


**f.**

- a. Unscrew pipetting unit: Separate the pipetting unit from the handle.
- b. Remove the pipetting housing: Pull the sliding catches of the pipetting housing cover out to the side as far as they will go and remove the pipetting housing.
- c. Move the piston to the lower position: Push the plunger down until the pistons are in the lowest position.
- d. Remove the shaft attachment: Pull out the shaft attachment from the side.
- e. Gently push out the piston shaft unit directly at the pipetting housing.
- f. Remove the piston shaft unit. Only the entire piston–shaft unit can be replaced.


- g.** Pull the piston with the seal out of the shaft.

- h.** Dip the brush into the piston grease supplied and wipe it well against the vessel wall. Hold brush against seal, and turn piston with seal 1–2 turns. Insert the piston with seal back into the shaft.


**NOTICE!**

Apply piston grease in a thin layer to the seal with a brush. Used only the grease supplied (order no.: G9999-1012). Do not use silicone grease or fluorstatic grease!

- i.** Assemble piston and shaft. Replace the piston–shaft unit. To do this, push the shaft into the holder, and carefully press the piston into the original position. The shafts must then be aligned again.

Proceed as described with all 4 or 6 piston shaft units on this page. To disassemble/replace the remaining 4 or 6 channels, turn the pipetting unit.

- j.** Reassemble the pipetting unit. To do this, slide the stem attachment over the stems until it clicks into place. Observe the direction from channel 8/12 to channel 1 (see marking on the pipetting housing cover). Then fix the pipetting housing again and mount the pipetting unit on the handle section. Check the instrument for leaks, free movement, and precision

**NOTICE!**

The pipetting unit must be screwed onto the handle in a counter-clockwise direction until it clicks into place. Do not pull the pipetting unit downwards.

- k.** Carry out reference run: Carry out reference run (rEF).

## 10.4 Charge and replace rechargeable battery

A fully charged rechargeable battery allows approx. 8 h of continuous pipetting (over 4000 pipetting cycles) of samples of water-like viscosity and density.

### NOTICE!

- › Before charging, ensure that the power adapter is suitable for the voltage present in the laboratory.
- › The instrument must not be charged in an explosive environment.
- › The rechargeable battery can be charged only inside the ErgoOne® E Multi-Channel pipettes.

### Recharge rechargeable battery

- a. Insert the charging cable plug of the power adapter into the socket provided for this purpose at the top of the .
  - ⇒ The charging process starts automatically.
  - ⇒ During the charging process, the bars of the rechargeable battery capacity indicator continually run from the bottom to the top.
  - ⇒ Once the bars of the display have stopped, this means that the rechargeable battery is fully charged.

### Pipetting during the charging process?

During charging, you can continue to work with the ErgoOne® E Multi-Channel pipettes. When the rechargeable battery is fully discharged, it takes a few minutes to reach a certain minimum charge capacity. This is necessary for safe operation of the instrument. The last settings made are stored in the EEPROM of the instrument. In the case of complete discharge or when changing the battery, these settings remain saved!

### Replace rechargeable battery

- a. Open the cover of the battery compartment, remove the rechargeable battery, and pull the plug out of the socket.

- b. Insert the plug of the new rechargeable battery into the socket and insert the new rechargeable battery.

- c. Replace the cover of the battery compartment and close it.

Remove the rechargeable battery from the instrument during longer breaks in operation.

## Battery indicator after reinserting a battery pack

After inserting a rechargeable battery, the full capacity indicator appears in the display with a flashing frame (the instrument does not yet recognize the charge status at first). After 3.5 h charging time – safe complete charging of the battery – the frame stops flashing.

### NOTICE!

Always charge for at least 3.5 h after inserting a battery!  
The full charge capacity is reached after several charge/discharge cycles!

## 10.5 Rechargeable battery regeneration function

### Refresh function

To extend the service life and increase the performance of the batteries, the ErgoOne® E Multi-Channel pipette has a regeneration function (refresh function). This function allows the batteries to be fully discharged and recharged under program control. To optimize the performance of the batteries, the refresh function should be used from time to time.

### Carry out refresh function

- a. Insert the charging cable plug (connection) of the power adapter into the socket provided for this purpose at the top of the ErgoOne® E Multi-Channel pipette.
- b. Press the lower arrow key for longer than 3 s. During discharging, the capacity bars of the battery display run continuously from top to bottom.
- c. After discharging (up to 3 h), the charging process (3.5 h) is started automatically. During charging, the capacity bars of the battery display run continuously from bottom to top.

### Canceling the refresh function

Pressing any key terminates the program. The instrument automatically switches back to the standard pipetting mode (PIP) and to the nominal volume, and the normal loading process is automatically started, see Charge and replace rechargeable battery. Unplugging the power supply also terminates the program. Cancellation of the refresh function must not be performed at the end of the discharge cycle.

## 11. Troubleshooting

Fault	Display	Cause	Corrective action
Instrument does not respond	ERR 1	Rechargeable battery empty or defective	Charge the battery for at least 5 min without actuation. Then continue working with the charging cable only until the battery is recharged. Replace the rechargeable battery if necessary
		Electronic components defective	Send instrument for repair
Instrument does not respond	ERR 2	Electronic components defective	Send instrument for repair
Instrument does not respond	ERR 3	Unforeseen program error	Error confirmation by pressing the 'Enter' key, instrument is reinitialized
Instrument does not respond	ERR 4	No rechargeable battery in the instrument	Insert rechargeable battery
		Rechargeable battery defective	Replace rechargeable battery
		Electronic components defective	Send instrument for repair
Tip dripping/ instrument leaking or volume error	—	Unsuitable tip	Use only high-quality tips
		Tip not seated tightly	Press tip on more firmly/other interchangeable clip
		Piston, shaft, or seal dirty or damaged	Clean instrument/replace seal, grease piston
No indication in the display	—	Electrostatic discharge	Remove rechargeable battery and insert again
		Electronic components defective	Send instrument for repair
Aspiration is not possible	—	Motor has no connection to the pipetting unit.	Carry out reference run (rEF), see Reference run (rEF)

## 12. Product markings

Symbol or number	Meaning
	With this mark, we confirm that the product complies with the requirements set out in the EC Directives and has been subjected to the specified testing procedures.
	UKCA: United Kingdom Conformity Assessed With this mark, we confirm that the product meets the requirements specified in the UK Designated Standards .
	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
XXZXXXXX	Serial number
	Follow the instructions listed on the instrument, the accessory parts and in the operating manual.
	The instrument or rechargeable battery should be disposed of properly.
	China RoHS (EFUP) EFUP defines the time period in years, in which the hazardous materials found in the electrical and electronic equipment should not leak or mutate under normal operating conditions. When used under normal conditions, such electrical or electronic products do not lead to severe environmental pollution, severe personal injuries or damage to the user's property.
	This electrical instrument may not be disposed of in household waste.

## 13. Order Information

### 13.1 Ordering Information

**ErgoOne® E-8 Channel incl. AC adapter** for Continental Europe + UK (100 - 240 V/50 - 60 Hz)

Volume	0.5 - 10 µl	1 - 20 µl	5 - 100 µl	15 - 300 µl	50 - 1250 µl
Article No.	Article No.	Article No.	Article No.	Article No.	Article No.
G9008-0010	G9008-0020	G9008-0100	G9008-0300	G9008-1250	


**ErgoOne® E-12 Channel incl. AC adapter** for Continental Europe + UK (100 - 240 V/50 - 60 Hz)

Volume	0.5 - 10 µl	1 - 20 µl	5 - 100 µl	15 - 300 µl	50 - 1250 µl
Article No.	Article No.	Article No.	Article No.	Article No.	Article No.
G9012-0010	G9012-0020	G9012-0100	G9012-0300	G9012-1250	

### 13.2 Spares

#### 13.2.1 ErgoOne® E Multi-channel pipettes to 300 µl

The appearance and dimensions of the spare parts correspond to the respective nominal volume.

**A** Piston-shaft unit

**B** Nose cone


**C** O-ring

**D** Seal

Volume	A	B*	C	D
0.5 - 10 µl	P7141-0101	P4031-0101	P4021-0101	P7111-0111
1 - 20 µl	P7141-0201	P4031-0201	P4021-0101	P7111-0211
5 - 100 µl	P7141-1001	P7131-1001	P7121-1301	P7111-1001
15 - 300 µl	P2041-3001	P2031-3101	P7121-1301	P7111-3001

\* incl. seal, O-ring and imounting wrench

### 13.2.2 ErgoOne® E Multi-Channel pipette 1250 µl

**E** Piston-shaft unit

**F** O-ring

Volume	E	F
50–1250 µl	P7141-1211	P7121-1201

### 13.3 Additional accessories

Description	Article No.
Stand for one ErgoOne® E Multi-Channel pipette up to 300 µl	G9999-1004
Shelf Mount	G9999-1009
Wall Mount	G9999-1010
Universal stand for one ErgoOne® E Multi-Channel pipette - all volumes	G9999-1011
Replacement battery	G9999-1005
AC adapters for Continental Europe + UK (100 - 240 V/50 - 60 Hz)	G9999-1000
Silicone grease variants up to 300 µl	G9999-1007
Fluorstatic grease	G9999-1008
Piston grease for 1250 µl variant	G9999-1012

## 14. Repairs

### 14.1 Sending for repair

**NOTICE!**

Transporting of hazardous materials without a permit is a violation of federal law.

#### **Clean the instrument thoroughly and decontaminate!**

- › When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- › Only send the device without a battery installed.
- › Shipment is at the risk and the cost of the sender.

#### **Contact addresses:**

##### **Starlab International GmbH**

Neuer Höltigbaum 38

22143 Hamburg

T: +49 (0)40 675 99 39 0

F. +49 (0)40 675 99 39 20

info@starlab.de

**all other contacts: [www.starlabgroup.com](http://www.starlabgroup.com)**

## 15. Calibration service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3 - 12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

Just send in the instruments to be calibrated, accompanied by an indication of which kind of calibration you require. Your instruments will be returned within a few days together with a Certificate of Calibration. For further information, please contact Starlab Pipette Calibration Service or your Starlab partner.

## 16. Warranty and guarantee

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

ErgoOne® E pipettes are guaranteed for 12 months from the delivery date, providing that the product is used according to the operating instructions. The guarantee is not valid in the event of damage or defects caused by improper use. The guarantee covers repairs to or replacement of the equipment. Our “General terms and conditions” apply. Starlab reserves the right to make changes to specifications without prior notice in order to implement short-term innovations.

## 17. Disposal

This symbol means that at the end of their service life, batteries/accumulators and electronic devices must be disposed of separately from household waste (unsorted municipal waste).

Electronic devices must be disposed of in accordance with Directive 2012/19/ EU of the European Parliament and of the Council from July 04, 2012 on waste from electrical and electronic equipment and in compliance with national disposal regulations.

Both batteries and accumulators (rechargeable batteries) contain materials that can be damaging to the environment and human health. Therefore, they must be properly disposed of in accordance with Directive 2006/66/EC of the European Parliament and of the Council from September 06, 2006 on batteries and accumulators and in compliance with national disposal regulations. Only dispose of fully discharged batteries and accumulators.

<b>User Manual</b>	<b>3–40</b>
<b>Bedienungsanleitung</b>	<b>41–78</b>
<b>Mode d'emploi</b>	<b>79–116</b>
<b>Istruzioni d'uso</b>	<b>117–154</b>
<b>Instrucciones de manejo</b>	<b>155–192</b>

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**Starlab (UK), Ltd**
**Starlab S.r.l.**