Mini Dry Bath Incubator

Instruction Manual

N2400-4020

Mini Dry Bath Incubator

CE

www.starlabgroup.com

STAR
LAB

ENGLISH.....3
DEUTSCH 10
FRANCAIS 18
ITALIANO 25

DECLARATION OF CONFORMITY

STARLAB certifies that

N2400-4020 Mini Dry Bath Incubator

*is manufactured and inspected in accordance with
the following International and European standards and
is in compliance with council directive 98/79/EC.*

RoHS 2011-65-EU

Harmonized Standards:

EN 61326-1:2006; EN 61326-2-2:2006

CISPR 11:2003+A1:2004 + A2:2006, Class B

EN 61000-3-2:2006+A1:2009 and EN 61000-3-3:2008

IEC 61000-4-2:2008, IEC 61000-4-3:2010

IEC 61000-4-4:2011, IEC 61000-4-5:2005

IEC 61000-4-6:2008, IEC 61000-4-11:2004

EN 61010-1:2010

As an ISO 9001 and ISO 13485 certified company,
our working procedures and processes are checked and audited regularly
and the quality of our products is continuously monitored.

STARLAB International GmbH

This unit must be disposed of in an environmentally-friendly manner.

Mini Dry Bath Incubator

Cat. No. N2400-4020

1. WARNING!

This equipment has been tested and found to comply with the limits for the CE regulation. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy, and if not installed and used in accordance with the instruction manual may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case, the user will be required to correct the interference at their expense. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. It is strongly recommended that only those who have read this manual operates the unit.

Please read and follow the instructions in this manual:

1. Do not make any alterations to the unit as these could result in personal and/or laboratory hazards, as well as invalidating the warranty.
2. Use a properly grounded electrical outlet with the correct voltage and current.
3. Disconnect from the power supply before maintenance and servicing. Servicing and maintenance should be performed only by qualified personnel.
4. Never use the unit without the safety cover on the bottom of the unit correctly in position.
5. Do not use the unit if there is any sign of damage, and replace any damaged parts before use.
6. Do not use in the presence of flammable or combustible material as fire or explosion may result as the unit contains components which may ignite such materials.
7. Use appropriate materials and operate correctly to avoid possible hazards of explosion, implosion or release of toxic or flammable gases arising from overheated materials.
8. Do not use the lid when using the unit in water bath mode.
9. Always use the supplied block lifter to remove hot blocks, and wear appropriate protection to avoid burning your hands.

ATTENTION: Hot surface!

2. Safety information

As with any electrical device, safety precautions should be taken. Before connecting to the mains check that the supply voltage is within the range stated on the rating label and that the device is earthed. Place the unit on a level, dry surface in a safe location. The unit MUST NOT touch anything in the surrounding area. Always follow the relevant safety precautions for all chemicals / dangerous materials and Hot Surface!

2.1 Environmental conditions

Ensure the instrument is installed and operated strictly in the following conditions:

1. Indoor use only.
2. ≤ 95 % RH.
3. 75 kPa – 106 kPa.
4. Altitude must not exceed 2000 meters.
5. Ambient to 40 °C operating temperature.
6. Pollution degree: 2.
7. Mains supply voltage fluctuations up to ± 10 % of the normal voltage.

2.2 Avoiding electrical shock

The Mini Dry Bath Incubator has been designed for use with shielded wires to minimise any potential shock hazard to the user. DO NOT use the unit with unshielded wires.

1. In the event of accidental spills into the instrument, the unit must be allowed to dry for AT LEAST 2 HOURS and restored to NORMAL CONDITION before reuse.
2. NEVER connect or disconnect the wire lead (9) when the power switch (8) is on.
3. WAIT at least 5 seconds after stopping a run before handling output leads or connected apparatus.
4. ALWAYS make sure that your hands, the work area and the unit are clean and dry before touching or operating the unit.
5. ONLY connect the power cord to a properly grounded AC outlet.

2.3 Avoiding damage to the unit

1. Do not attempt to operate the unit if it is damaged in any way.
2. Protect the unit from damage, corrosive agents and extreme temperatures (direct sunlight etc).
3. For proper ventilation, the unit should be situated and operated with at least 10 cm of clear space behind and 5 cm of clear space on either side.
4. Do not operate the unit outside of the environmental conditions listed above.
5. Before using any cleaning or decontamination method (except those recommended in this manual), you should check with STARLAB or your STARLAB distributor that the proposed method will not damage the equipment.

2.4 Equipment operation

Follow the guidelines below to ensure safe operation of the unit:

1. The displayed temperature should be checked with an external temperature probe at least once a day to ensure it is not overheating.
2. DO NOT use dangerous chemicals or other materials which may cause hazard of explosion or damage to the unit.
3. Do not apply lids to the tubes heated inside the Mini Dry Bath Incubator in order to prevent the possible hazard of explosion and damage. To prevent caps from popping open it is recommended to use cap locks, like STARLAB's E1405-9000 for 0.5 ml tubes or I1415-1508 for 1.5/2.0 ml tubes.
4. A temporary conductivity caused by condensation might occur even though the unit is rated Pollution Degree 2 in accordance with IEC 664.

2.5 Symbols

The symbols used on the Mini Dry Bath Incubator are explained below.

Danger Electricity: indicates a potential shock hazard may exist.

ATTENTION: Hot surface!

3. Product overview

The Mini Dry Bath Incubator is suitable for a wide variety of applications and provides excellent temperature control to deliver accurate and reliable experimental results, again and again. The single block dry bath uses interchangeable heating blocks for a variety of applications, including restriction enzyme digestion, denaturing DNA, Blood Urea Nitrogen, melting agar, coagulation studies, hybridization, and Hot Start thermocycle reactions. Due to its moulded aluminium alloy chamber, the unit can also be used as a mini water bath. The unit incorporates a PID controller for easy temperature selection, rapid heating and temperature stability. The temperature may be set in 0.1 °C increments from 5 °C above ambient to 100 °C. The timer may be set from 1 to 9999 minutes.

Features

- Microprocessor controller with digital display
- User temperature calibration
- Leak-proof, moulded aluminium alloy heating chamber
- A variety of blocks available

4. Controls and features


Fig 1. Front View


Fig 2. Rear View

5. Operating instructions

The Mini Dry Bath Incubator is a pre-installed instrument and is ready for use.

1. Locate the unit on a sturdy, level surface in a safe, dry place away from laboratory traffic.
2. Ensure that the **power switch (8)** is OFF, then attach the power cord into socket (9) at the rear of the unit and plug into a grounded AC outlet of appropriate voltage.
3. Insert an aluminium block or the appropriate water volume in the unit.
Note: Do not use the lid in water bath mode!
4. Turn the **Power Switch (8)** to ON.
5. Perform a temperature calibration procedure when using the instrument for the first time (see 6).
6. Press the **UP (1)** or **DOWN (2)** button to set the desired temperature.
7. To set the timer for the heating process, press **MODE (4)** button, and then press **UP (1)** or **DOWN (2)** button to adjust to required time. The unit will stop heating and an alarm will sound when time is up.
8. Press the **RUN/STOP (3)** button to start heating.

9. To reset the timer during heating, simply press **RUN/STOP (3)** to deactivate heating.
10. Press **RUN/STOP (3)** again to start the heating, after setting the timer.

6. Temperature calibration

The Mini Dry Bath Incubator with optional blocks has been calibrated as a set. However, different block types or water sources may have different influences and may cause different results. For optimum accurate temperature control, or when swapping to different block types or water, the unit should be calibrated in accordance with the procedure outlined below.

1. Insert a calibrated laboratory thermometer into the thermometer holding port located on the block or into the middle of the chamber if being used as a water bath.
2. Press and hold the **RUN/STOP (3)** button while switching the **Power Switch (8)** ON.
The unit will omit a beep and the LED display will start flashing here. 
This indicates the unit is in Calibration Mode.
3. Press the **UP (1)** or **DOWN (2)** to adjust the display value to the temperature you want to control accurately, and then press the **RUN/STOP (3)** button.
4. As the unit reaches the target temperature the timer will start to count down.
Wait for approximately 30 mins, when you should hear a "beeping" sound.
5. Adjust the display value according to the value shown on the thermometer. Now press the **RUN/STOP (3)** button and the screen will display 
which indicates it is saving the calibrated value.
6. The calibration procedure is now finished. Wait for a few more minutes for the microprocessor to automatically adjust the displayed temperature to the same temperature the thermometer measured.

7. Power recover function

1. If the unit is turned off during the heating stage and then turned back on again, an alarm will sound and the display will countdown from 10 to 0, then it will recover back to the last state. The timer will also resume to the last state.
2. Press the **RUN/STOP (3)** button during the countdown period to terminate the alarm and the timer. The unit is then ready for operating again.

8. Troubleshooting

Many operating problems may be solved by following the instructions in this manual. Some suggestions for troubleshooting are given below. Should these suggestions not resolve the problem, please contact STARLAB or your STARLAB distributor for assistance, providing full information on the problem.

PROBLEM	RECOMMENDATIONS
LED light does not light up	Check the FUSE
	Ensure that the AC power switch (8) is ON
	Check that the power cord is plugged securely in the back of the unit and also plugged into a grounded AC outlet of the appropriate voltage

8.1 Operation protection mechanism

PROBLEM	CAUSE	SOLUTION
 Alarms with short buzz sound	Overheating: when the temperature reaches 5 °C above the set value	System will automatically return to normal
 Alarms with long buzz sound	Overheating: when the temperature reaches above 110 °C	Restart the unit if AL1 and Er1 occur simultaneously, restart the unit (as Er1 takes priority). If the problem persists contact STARLAB for assistance
 Alarms with long buzz sound	Unit temperature sensor is abnormal	Restart the unit. If the problem persists contact STARLAB for assistance

9. Specifications

Controller	Digital Microprocessor Controller
Display	LED
Operating Temperature	Ambient to 40 °C
Temperature Control Range	<i>Dry Bath</i> : 5 °C above ambient to 100 °C. <i>Water Bath</i> : 5 °C above ambient to 90 °C
Temperature Increase	0.1 °C
Temperature Calibration	Yes
Temperature Uniformity	± 0.2 °C at 37 °C
Temperature Accuracy	± 0.25 °C at 37 °C
Timer	1 to 9999 minutes, continuous
Safety Device	Leak-proof heating chamber. Over-heat alarm protection
Heating Chamber Material	PTFE-coated aluminium
Block Material	Aluminium alloy
Block Dimensions (WxDxH)	71 x 47 x 32 mm (71 x 47 x 75 mm, N2400-4025 and N2400-4026)
Water Bath Volume (approx.)	85 ml
Heating Power	50 W
Rated Voltages	100-240 V~, 2 A 50/60 Hz
Unit Weight and Dimensions	0.6 kg (without blocks) / 125 x 153 x 97 mm (WxDxH)

Features and specifications subject to change without notice.

10. Maintenance

The Mini Dry Bath Incubator may be cleaned with a moist cloth containing a mild soap solution. The heating chamber and metal blocks are constructed of aluminium alloy and may be cleaned with any commercial aluminium cleaners available on the market.

The heating surface contains a PTFE coating. Please avoid contact with sharp objects.

11. Ordering information

N2400-4020	Mini Dry Bath Incubator*
N2400-4021	Mini Dry Bath Metal Block for 32 x 0.2 ml PCR Tubes (ø 6.35 mm)
N2400-4022	Mini Dry Bath Metal Block for 12 x 0.5 ml Microcentrifuge Tubes (ø 8.0 mm)
N2400-4023	Mini Dry Bath Metal Block for 12 x 1.5 ml Microcentrifuge Tubes (ø 10.8 mm)
N2400-4024	Mini Dry Bath Metal Block for 12 x 1.5/2.0 ml Microcentrifuge Tubes (ø 10.9 mm)
N2400-4025	Mini Dry Bath Metal Block for 6 x 15 ml Centrifuge Tubes (ø 17.3 mm)
N2400-4026	Mini Dry Bath Metal Block for 2 x 50 ml Centrifuge Tubes (ø 29.2 mm)

*Metal blocks are not included and must be purchased separately.

Please see our website www.starlabgroup.com for information on our wide range of consumables including microcentrifuge tubes, centrifuge tubes and PCR tubes.

12. Warranty

The Mini Dry Bath Incubator is under warranty against defects in materials and workmanship under normal service, for two years from the shipping date to the purchaser. This warranty excludes damages resulting from transport, misuse, carelessness, or neglect. STARLAB's liability under the warranty is limited to the receipt of reasonable proof from the customer that the defect is embraced within the terms of the warranty. All claims made under this warranty must be presented to STARLAB or a STARLAB Distributor within two years following the date of delivery of the product.

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STARLAB GmbH
STARLAB (UK), Ltd
STARLAB S.r.l.