Orbital Shaker

Instruction Manual (v.1)

N2400-8030

Orbital Shaker

(Shown with N2400-8011 Universal Platform with Clamping Rolls)

CE

www.starlabgroup.com

S T A R
L A B

ENGLISH 3
DEUTSCH10
FRANCAIS 17
ITALIANO24

EC DECLARATION OF CONFORMITY

STARLAB certifies that

N2400-8030 Orbital Shaker

has been manufactured and inspected in accordance

with the following European standards:

2006/95/EC Low Voltage Directive

2004/108/EC Electromagnetic Compatibility Directive

Harmonized Standards:

EN 61010-1:2010

EN 61326-1:2006

As an ISO 9001 and ISO 13485 certified company,
our working procedures and processes are checked and audited regularly
and the quality of our products is continuously monitored.

STARLAB International GmbH

ORBITAL SHAKER

Cat. No. N2400-8030

1. Safety instructions

Before operation

- Read the Instruction Manual carefully before use.
- Ensure that only trained staff use the instrument. The instrument may only be opened by a qualified and trained technician.
- Make sure the socket is earthed before use.
- Set the instrument on a stable, clean, dry, non-slip and fireproof surface.
- Do not operate in hazardous areas.
- Safe operation is only guaranteed with the original accessories from STARLAB. You will find the complete list in chapter 8.
- Accessories must be securely attached.
- Before each use, check both the instrument and the accessories for damage. Do not use damaged components.
- The operation voltage noted on the label must match.
- Always disconnect the plug before fitting accessories.
- When working, wear personal protection to avoid the risk from splashing of liquids and possible glass breakage which can result from mechanical vibrations.

During operation

- Do not touch any running parts or attempt to move the instrument during operation.
- Set the "speed-control knob" to the lowest position first and then gradually increase the speed.
- Reduce the speed if the instrument does not run smoothly.
- Firmly secure the accessories and vessels in place to avoid damage or risk.
- Preparation of certain samples may lead to dangerous flammability. Only process samples that will not react dangerously.
- To ensure safety, follow the instructions for use for accessories. Please switch off the power before assembly of accessories and ensure both the instrument and the accessories are intact before switching on again.

2. Product description

The instrument is designed for mixing liquids. Normal use: The instrument is designed for use in professional laboratories. Do not use it in residential areas.

3. Inspection

Unpack the equipment carefully and check for any damages. If the instrument is faulty you should contact STARLAB or your STARLAB distributor promptly.

The packing includes:

- One Orbital Shaker
- One EU-power supply and one UK-power supply
- One instruction manual

A platform is not included and must be purchased separately. See chapter 8.1.

Please check all items are present.

4. Ambient conditions

Power supply: 100 – 240 V 50/60 Hz

Ambient temperature: 4 – 40° C

Relative humidity: < 80 %

Place the instrument on a safe and stable surface.

5. Operation

Installation

An initial trial operation should be performed as follows:

Ensure the instrument is standing on a stable, even surface. If the mounting surface is not even, you can improve the stability of the instrument by adjusting the feet. To do this, turn the appropriate foot downward with the wrench until the instrument is standing securely on the surface.

- Make sure the required operating voltage and power supply voltage match.
- Ensure the socket is securely earthed.
- Plug in the power cable and power ON the instrument.
- LCD displays the safe rotary speed limit.

- LCD displays operating mode.
- Turn the speed control knob to set the required shaking speed.
- Press the speed control knob and the shaking function is switched ON.
- Press the speed control knob again, and the shaking function is switched OFF.
- If the operations above are normal, the instrument is ready to operate following the instructions in the "operating" chapter. If the operations are not normal, please contact STARLAB or your STARLAB distributor.

5.1 Control elements


Image shows Shaker with Universal Attachment with Clamping Rolls (N2400-8011) not included.

Speed control knob	Sets the rated shaking speed. Shaking is switched ON or OFF by pressing the knob
Time knob	Sets the working time. Shaking is switched ON or OFF by pressing the knob
LCD	Displays the status of instrument and any time/speed set.
Power switch	Switches the instrument ON or OFF.

5.2 Operating

Set time

The rated shaking time can be set by turning the time. After the instrument is powered OFF and restarted, the set time is erased to zero and it switches into offline operating mode.

If both a target time (max. 19h 59min) and the speed are set, users can activate the instrument by pressing the time knob and the count-down starts.

If:

A: Pressing the time knob, the speed and the time will stop. Pressing the time knob again restarts the shake function for the pre-set time.

B: Pressing the speed knob, the speed and the time will stop. Pressing the speed knob again will restart the shake function but not the time. There will be flashing on the right LCD (time). Pressing the time knob again will restart the time function at the pre-set time.

Note: The current set time can be varied at any time.

Setting the speed

The rated speed and upper speed limit can be set by turning the speed control knob. The continuous mode can be switched ON by pressing the speed control knob without the time setting, and switched OFF by pressing the speed control knob again.

Adjust the speed control knob slowly in order to keep the instrument running smoothly.

Additional settings

Mode A

Operating mode A is the default factory setting. After turning ON the instrument with the main switch, the shaking and timing functions are turned OFF. The LCD displays the set time and speed values. They are adopted or varied when the corresponding functions are turned ON.

After turning ON the instrument, the following will be shown on the display.

The left LCD panel displays "SAFE" and the right LCD panel displays the upper speed limit (--- rpm) which can be set by pressing and turning the speed control knob at the same time.

The left LCD displays "StA", the right LCD displays the running mode "A" or "B" for about 2 seconds.

The set or stored speed value is shown in the left display, and the set or stored time value is shown in the right display. Turn the speed control knob and time knob to set the speed and time respectively.

After pressing the speed control knob or time knob, the instrument begins to run at the set speed. The actual speed and time remaining are shown on the LCD panels. If the remaining time has reached zero, the shaker stops its motion. By pressing the speed control knob alone, the instrument will work in continuous mode (with the time function OFF).

Mode B

The upper speed limit that has been set is taken over from operating mode A and cannot be modified. Operating mode B is not the default setting.

After turning ON the instrument, the following process will be shown on the display.

- The left LCD panel displays "SAFE" and the right LCD panel displays the upper speed limit (---rpm) which cannot be changed.
- The LCD displays run modes "StA" for about 2 seconds.
- The set speed value is shown in the left display, and the set time value is shown in the right display. Then set the motor rotary speed and time.
- After pressing the speed knob or time knob, the instrument begins to run at the set speed. The actual speed and time remaining are shown in the LCD. If the time has reached zero, the shaker stops its motion.

Switching modes

Switch the operating mode as follows:

- Turn OFF the instrument with the main switch.
- Hold down both the speed control and the time knob together and turn ON the instrument at the mains switch. After about 5 seconds, you can let go of the knobs.
- Switch between operating Modes A and B in order.

Supported load

- In order to ensure safe operation, the shaker must only be operated with a maximum load of 7.5 kg.
- Make sure the shaking platform is always clean and level. The individual shaking vessels must be placed evenly and fastened securely in the middle of the shaking platform whether either single or multiple vessels are being used.

Troubleshooting

- Instrument does not work when switched ON:
 - check the power cable is fitted securely and the power is switched on.
 - check that the fuse is not faulty.
- The speed cannot reach set value:
 - check that the platform is not overloaded.
- The motor does not start by pressing the speed knob and time knob:
 - check that the time is not set to zero.

Interface

- The shaker can be equipped with a special adapter and cable which can be used to link it to an external device (Standard 9-pin interface to connect with the PC).
- The cable between the shaker and the computer uses EIA Standard RS232C communication line, corresponding to the DIN 66020 interface.
- Transmission method: Asynchronous signal transmission.
- Mode of transmission: Fully Duplex.
- 1 start bit; 8 character bits; 1 stop bit.
- Transmission speed: 9600 bit/s.

6. Maintenance and cleaning

- Proper maintenance and cleaning will keep the instrument in good working order and will lengthen its life time.
- Be careful not spray cleaning fluid into the instrument. Unplug the power cord when cleaning.
- Please use only Isopropyl alcohol and water containing tenside. Before using other methods or cleaning agents the user must ascertain with the manufacturer that this will not destroy the instrument.
- Wear protective gloves during cleaning of the instrument.
- When not in use for long periods, switch off the instrument and store it in a dry, clean and stable place at room temperature.

Warranty

This instrument is under warranty to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty claim expires if the instrument is opened by untrained or unqualified staff. Our "General terms and conditions" apply.

Service

When help needed, please contact STARLAB or your STARLAB distributor.

See www.starlabgroup.com for your local contact.

Please provide the customer care representative with the following information:

- Serial Number (on the rear panel)
- Description of problem (i.e., hardware or software)
- Methods and procedures adopted to resolve the problems
- Your contact information

7. Technical Data

Parameter	Data
Voltage	100 – 240 V
Frequency	50/60 Hz
Power	30 W
Shaking movement	Orbital
Orbital diameter	10 mm
Max. shaking weight (with attachment)	7.5 kg
Motor type	External rotor / brushless motor
Motor rating input	28 W
Motor rating output	15 W
Speed range	100 – 500 rpm
Speed display	LCD
Timer	Yes
Timer display	LCD
Time setting range	1 – 1199 min
Operation type	Timer / Continuous
Dimensions (WxDxH)	370 x 420 x 100 mm
Weight	13.5 kg
Permissible ambient temperature	4 – 40 °C
Permissible relative humidity	< 80 %
Protection class acc. to DIN EN60529	IP21
RS232 interface	Yes

8. Ordering information

Cat. No.	Description
N2400-8030	Orbital Shaker

8.1 Accessories

Cat. No.	Description
N2400-8011	Universal Attachment with clamping rolls
N2400-8021	Universal/ Replacement attachment with bungee cords

Please visit our website www.starlabgroup.com, for details on our comprehensive range of plastic laboratory consumables.

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STARLAB GmbH
STARLAB (UK), Ltd
STARLAB S.r.l.