













Metal block information.

	Art. No	Description	Well shape	Well diameter	Well depth	Block dimensions (WxDxH)
	N2400-4021	Metal block for 32 x 0.2 ml PCR tubes		6.35 mm	19.0 mm	71 x 47 x 32 mm
	N2400-4022	Metal block for 12 x 0.5 ml microcentrifuge tubes		8 mm	25.0 mm	71 x 47 x 32 mm
	N2400-4023	Metal block for 12 x 1.5 ml microcentrifuge tubes		10.8 mm	28.5 mm	71 x 47 x 32 mm
	N2400-4024	Metal block for 12 x 1.5/2.0 ml microcentrifuge tubes		11 mm	30 mm	71 x 47 x 32 mm
	N2400-4025	Metal block for 6 x 15 ml centrifuge tubes		17.3 mm	70 mm	71 x 47 x 75 mm
	N2400-4026	Metal block for 2 x 50 ml centrifuge tubes		29.2 mm	72 mm	71 x 47 x 75 mm

Blocks are supplied with a screw-in handle for insertion/removal.

 www.starlab.click/mini-dry-bath



N2400-4020 Mini Dry Bath.

