


Metal block information.

	Art. No	Description	Well shape	Well diameter	Well depth	Block dimensions (WxDxH)
	N2400-4100	Metal block for 64 x 0.2 ml tubes		6.35 mm	17.0 mm	80 x 105 x 50 mm
	N2400-4101	Metal block for 20 x 1.5 ml or 2.0 ml tubes		11.00 mm	33.0 mm	80 x 105 x 50 mm
	N2400-4102	Metal block for 12 x 15 ml centrifuge tubes		17.30 mm	45.0 mm	80 x 105 x 50 mm
	N2400-4103	Metal block for 4 x 50 ml centrifuge tubes		29.00 mm	45.0 mm	80 x 105 x 50 mm
	N2400-4104	Metal block for 20 x 1.5 ml tubes (V-bottoms)		11.00 mm	33.4 mm	80 x 105 x 50 mm
	N2400-4105	Metal block for 20 x 0.5 ml tubes		7.90 mm	25.0 mm	80 x 105 x 50 mm
	N2400-4106	Metal block for 20 Tubes up to 13 mm diameter		13.70 mm	36.5 mm	80 x 105 x 50 mm
	N2400-4107	Metal block for 96-well PCR plate (Single) ¹⁾		6.28 mm	15.0 mm	87 x 128 x 70 mm
	N2400-4108	Metal block for 96-well PCR plate (Dual) ²⁾		6.28 mm	15.0 mm	105 x 158 x 50 mm

Blocks are supplied with a screw-in handle for insertion/removal. ¹⁾ For use with N2400-4001 only. ²⁾ For use with N2400-4002 only.

 www.starlab.click/drybaths



N2400-4002 Dry Bath System (Dual Block), shown with blocks N2400-4105 and N2400-4103.



N2400-4001 Dry Bath System (Single Block).