

Freezer Racking Guide


► UPRIGHT FREEZERS



1. Stainless Steel Racking with Side Access
2. Aluminium Racking with Side Access
3. Stainless Steel Racking with Pull Out Shelves

The freezer racking listed on this flyer is suitable for use with any standard storage box up to 136 x 136 x 55 mm.

How to measure upright freezer racking



See catalogue or website for the external dimensions of the freezer racking listed on this guide

Eppendorf (New Brunswick)

EPPENDORF (NEW BRUNSWICK)					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
F570h CryoCube	1. Stainless Steel, Side Access	I2602-4400	25	16 (4x4)	400
	2. Aluminium, Side Access	I2612-4400	25	16 (4x4)	400
	3. Stainless Steel, Pull Out Shelves	I2632-4400	25	16 (4x4)	400
F740h comp3 CryoCube	1. Stainless Steel, Side Access	I2602-4700	18	28 (4x7)	504
	2. Aluminium, Side Access	I2612-4700	18	28 (4x7)	504
	3. Stainless Steel, Pull Out Shelves	I2632-4700	18	28 (4x7)	504
U360VIP	1. Stainless Steel, Side Access	I2602-4700	9	28 (4x7)	252
	2. Aluminium, Side Access	I2612-4700	9	28 (4x7)	252
	3. Stainless Steel, Pull Out Shelves	I2632-4700	9	28 (4x7)	252
U410	1. Stainless Steel, Side Access	I2602-4400	15	16 (4x4)	240
	2. Aluminium, Side Access	I2612-4400	15	16 (4x4)	240
	3. Stainless Steel, Pull Out Shelves	I2632-4400	15	16 (4x4)	240

EPPENDORF (NEW BRUNSWICK)					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
U535VIP	1. Stainless Steel, Side Access	I2602-4700	12	28 (4x7)	336
	2. Aluminium, Side Access	I2612-4700	12	28 (4x7)	336
	3. Stainless Steel, Pull Out Shelves	I2632-4700	12	28 (4x7)	336
U570	1. Stainless Steel, Side Access	I2602-4400	25	16 (4x4)	400
	2. Aluminium, Side Access	I2612-4400	25	16 (4x4)	400
	3. Stainless Steel, Pull Out Shelves	I2632-4400	25	16 (4x4)	400
U725VIP	1. Stainless Steel, Side Access	I2602-4700	18	28 (4x7)	504
	2. Aluminium, Side Access	I2612-4700	18	28 (4x7)	504
	3. Stainless Steel, Pull Out Shelves	I2632-4700	12	28 (4x7)	336

Freezer racking listed on this guide for both upright and chest freezers, is suitable for use with any standard storage box up to 136 x 136 x 55 mm, including the following storage boxes supplied by STARLAB:

- 50 mm high cardboard storage boxes
- StarStore 130 x 130mm range of storage boxes
- Cryogenic storage boxes (except 25- & 81-Tall)
- Expanding work rack
- StarRack 64



UPRIGHT FREEZERS

PHCBI (Panasonic / Sanyo)

PHCBI (PANASONIC / SANYO)						
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer	
MDF-DU302VX	1. Stainless Steel, Side Access	I2602-4500	12	20 (4x5)	240	
	2. Aluminium, Side Access	I2612-4500	12	20 (4x5)	240	
	3. Stainless Steel, Pull Out Shelves	I2632-4500	24	20 (4x5)	240	
MDF-DU502VH-PE MDF-DU502VX	1. Stainless Steel, Side Access	I2602-4600	16	24 (4x6)	384	
	2. Aluminium, Side Access	I2612-4600	16	24 (4x6)	384	
	3. Stainless Steel, Pull Out Shelves	I2632-4600	16	24 (4x6)	384	
MDF-DU702VH-PE MDF-DU702VX	1. Stainless Steel, Side Access	I2602-4600	24	24 (4x6)	576	
	2. Aluminium, Side Access	I2612-4600	24	24 (4x6)	576	
	3. Stainless Steel, Pull Out Shelves	I2632-4600	24	24 (4x6)	576	
MDF-DU900V-PE	1. Stainless Steel, Side Access	I2602-4500	14	20 (4x5)	280	
		I2602-4600	14	24 (4x6)	336	
				28 Total	44 Total	616 Total
	2. Aluminium, Side Access	I2612-4500	14	20 (4x5)	280	
		I2612-4600	14	24 (4x6)	336	
				28 Total	44 Total	616 Total
3. Stainless Steel, Pull Out Shelves	I2632-4500	14	20 (4x5)	280		
	I2632-4600	14	24 (4x6)	336		
			28 Total	44 Total	616 Total	
MDF-U33V	1. Stainless Steel, Side Access	I2602-4400	6	16 (4x4)	96	
		I2602-4500	6	20 (4x5)	120	
				12 Total	36 Total	216 Total
	2. Aluminium, Side Access	I2612-4400	6	16 (4x4)	96	
		I2612-4500	6	20 (4x5)	120	
				12 Total	36 Total	216 Total
3. Stainless Steel, Pull Out Shelves	I2632-4400	6	16 (4x4)	96		
	I2632-4500	6	20 (4x5)	120		
			12 Total	36 Total	216 Total	
MDF-U3386S	1. Stainless Steel, Side Access	I2602-4400	6	16 (4x4)	96	
		I2602-5400	6	20 (5x4)	120	
				12 Total	36 Total	216 Total
	2. Aluminium, Side Access	I2612-4400	6	16 (4x4)	96	
		I2612-5400	6	20 (5x4)	120	
				12 Total	36 Total	216 Total
3. Stainless Steel, Pull Out Shelves	I2632-4400	6	16 (4x4)	96		
	I2632-5400	6	20 (5x4)	120		
			12 Total	36 Total	216 Total	

PHCBI (PANASONIC / SANYO)						
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer	
MDF-U4186S	1. Stainless Steel, Side Access	I2602-3400	8	12 (3x4)	96	
		I2602-3500	8	15 (3x5)	120	
				16 Total	27 Total	216 Total
	2. Aluminium, Side Access	I2612-3400	8	12 (3x4)	96	
		I2612-3500	8	15 (3x5)	120	
				16 Total	27 Total	216 Total
3. Stainless Steel, Pull Out Shelves	I2632-3400	8	12 (3x4)	96		
	I2632-3500	8	15 (3x5)	120		
			16 Total	27 Total	216 Total	
MDF-U53V MDF-U55V MDF-U500VX	1. Stainless Steel, Side Access	I2602-4500	8	20 (4x5)	160	
		I2602-4600	8	24 (4x6)	192	
				16 Total	44 Total	352 Total
	2. Aluminium, Side Access	I2612-4500	8	20 (4x5)	160	
		I2612-4600	8	24 (4x6)	192	
				16 Total	44 Total	352 Total
3. Stainless Steel, Pull Out Shelves	I2632-4500	8	20 (4x5)	160		
	I2632-4600	8	24 (4x6)	192		
			16 Total	44 Total	352 Total	
MDF-U5386S	1. Stainless Steel, Side Access	I2602-4400	16	20 (4x5)	320	
	2. Aluminium, Side Access	I2612-4400	16	20 (4x5)	320	
	3. Stainless Steel, Pull Out Shelves	I2632-4400	16	20 (4x5)	320	
MDF-U7386S	1. Stainless Steel, Side Access	I2602-4500	24	20 (4x5)	480	
	2. Aluminium, Side Access	I2612-4500	24	20 (4x5)	480	
	3. Stainless Steel, Pull Out Shelves	I2632-4500	24	20 (4x5)	480	
MDF-U76V MDF-U700VX	1. Stainless Steel, Side Access	I2602-4500	12	20 (4x5)	240	
		I2602-4600	12	24 (4x6)	288	
				24 Total	44 Total	528 Total
	2. Aluminium, Side Access	I2612-4500	12	20 (4x5)	240	
		I2612-4600	12	24 (4x6)	288	
				24 Total	44 Total	528 Total
3. Stainless Steel, Pull Out Shelves	I2632-4500	12	20 (4x5)	240		
	I2632-4600	12	24 (4x6)	288		
			24 Total	44 Total	528 Total	

Haier

HAIER					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
DW-86L 288	1. Stainless Steel, Side Access	I2602-3500	12	15 (3x5)	180
	2. Aluminium, Side Access	I2612-3500	12	15 (3x5)	180
	3. Stainless Steel, Pull Out Shelves	I2632-3500	12	15 (3x5)	180
DW-86L 338	1. Stainless Steel, Side Access	I2602-4400	6	16 (4x4)	96
		I2602-4500	6	20 (4x5)	120
		12 Total	36 Total	216 Total	
	2. Aluminium, Side Access	I2612-4400	6	16 (4x4)	96
		I2612-4500	6	20 (4x5)	120
		12 Total	36 Total	216 Total	
	3. Stainless Steel, Pull Out Shelves	I2632-4400	6	16 (4x4)	96
		I2632-4500	6	20 (4x5)	120
		12 Total	36 Total	216 Total	
DW-86L 388	1. Stainless Steel, Side Access	I2602-3500	16	15 (3x5)	240
	2. Aluminium, Side Access	I2612-3500	16	15 (3x5)	240
	3. Stainless Steel, Pull Out Shelves	I2632-3500	16	15 (3x5)	240
DW-86L 486 DW-86K 728J	1. Stainless Steel, Side Access	I2602-4500	16	20 (4x5)	320
	2. Aluminium, Side Access	I2612-4500	16	20 (4x5)	320
	3. Stainless Steel, Pull Out Shelves	I2632-4500	16	20 (4x5)	320
DW-86L 628	1. Stainless Steel, Side Access	I2602-4500	20	20 (4x5)	400
	2. Aluminium, Side Access	I2612-4500	20	20 (4x5)	400
	3. Stainless Steel, Pull Out Shelves	I2632-4500	20	20 (4x5)	400
DW-86L 728	1. Stainless Steel, Side Access	I2602-4500	24	20 (4x5)	480
	2. Aluminium, Side Access	I2612-4500	24	20 (4x5)	480
	3. Stainless Steel, Pull Out Shelves	I2632-4500	24	20 (4x5)	480

HAIER					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
DW-86L 490	1. Stainless Steel, Side Access	I2602-4400	8	16 (4x4)	128
		I2602-4500	8	20 (4x5)	160
		16 Total	36 Total	288 Total	
	2. Aluminium, Side Access	I2612-4400	8	16 (4x4)	128
		I2612-4500	8	20 (4x5)	160
		16 Total	36 Total	288 Total	
	3. Stainless Steel, Pull Out Shelves	I2632-4400	8	16 (4x4)	128
		I2632-4500	8	20 (4x5)	160
		16 Total	36 Total	288 Total	
DW-86L 578J	1. Stainless Steel, Side Access	I2602-5500	16	25 (5x5)	400
	2. Aluminium, Side Access	I2612-5500	16	25 (5x5)	400
	3. Stainless Steel, Pull Out Shelves	I2632-5500	16	25 (5x5)	400
DW-86L 729	1. Stainless Steel, Side Access	I2602-5500	20	25 (5x5)	500
	2. Aluminium, Side Access	I2612-5500	20	25 (5x5)	500
	3. Stainless Steel, Pull Out Shelves	I2632-5500	20	25 (5x5)	500
DW-86L 828J	1. Stainless Steel, Side Access	I2602-5500	24	25 (5x5)	600
	2. Aluminium, Side Access	I2612-5500	24	25 (5x5)	600
	3. Stainless Steel, Pull Out Shelves	I2632-5500	24	25 (5x5)	600

Thermo

THERMO					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
TSX600V TLE60086V V/D	1. Stainless Steel, Side Access	I2602-5500	24	25 (5x5)	600
	2. Aluminium, Side Access	I2612-5500	24	25 (5x5)	600
	3. Stainless Steel, Pull Out Shelves	I2632-5500	24	25 (5x5)	600
TSX700	1. Stainless Steel, Side Access	I2602-5500	28	25 (5x5)	700
	2. Aluminium, Side Access	I2612-5500	28	25 (5x5)	700
	3. Stainless Steel, Pull Out Shelves	I2632-5500	28	25 (5x5)	700

Freezer Racking Guide

▶ CHEST FREEZERS



4. Stainless Steel Racking

5. Stainless Steel Racking with Spring Lock

6. Aluminium Racking

The freezer racking listed on this flyer is suitable for use with any standard storage box up to 136 x 136 x 55 mm.

How to measure chest freezer racking



See catalogue or website for the external dimensions of the freezer racking listed on this guide

PHCBI (Panasonic / Sanyo)

PHCBI (PANASONIC / SANYO)					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
MDF-192 MDF-193	4. Stainless Steel	I2602-7000	9	7	63
	5. Stainless Steel, Spring Lock	I2613-7000	9	7	63
	6. Aluminium	I2612-7000	9	7	63
MDF-393 MDF-394	4. Stainless Steel	I2602-9000	21	9	189
	5. Stainless Steel, Spring Lock	I2613-9000	21	9	189
	6. Aluminium	I2612-9000	21	9	189
MDF-594	4. Stainless Steel	I2602-1300	24	13	312
	5. Stainless Steel, Spring Lock	I2613-1300	24	13	312
	6. Aluminium	I2612-1300	24	13	312
MDF-794	4. Stainless Steel	I2602-1300	36	13	468
	5. Stainless Steel, Spring Lock	I2613-1300	36	13	468
	6. Aluminium	I2612-1300	36	13	468
MDF-C8V1	4. Stainless Steel	I2602-7000	6	7	42
	5. Stainless Steel, Spring Lock	I2613-7000	6	7	42
	6. Aluminium	I2612-7000	6	7	42

Eppendorf (New Brunswick)

EPPENDORF (NEW BRUNSWICK)					
Freezer Model	Freezer Racking Type	Freezer Racking Cat. No.	Racks per Freezer	Boxes per Rack	Boxes per Freezer
C340	4. Stainless Steel	I2602-1300	13	20	260
	5. Stainless Steel, Spring Lock	I2612-1300	13	20	260
	6. Aluminium	I2632-1300	13	20	260
C585VIP	4. Stainless Steel	I2602-1300	13	32	416
	5. Stainless Steel, Spring Lock	I2612-1300	13	32	416
	6. Aluminium	I2632-1300	13	32	416
C660	4. Stainless Steel	I2602-1300	13	40	520
	5. Stainless Steel, Spring Lock	I2612-1300	13	40	520
	6. Aluminium	I2632-1300	13	40	520
C760VIP	4. Stainless Steel	I2602-1300	13	40	520
	5. Stainless Steel, Spring Lock	I2612-1300	13	40	520
	6. Aluminium	I2632-1300	13	40	520

Copyright February 2019. All rights reserved. Technical specifications are subject to change without notice. Our general terms and conditions apply.