

Get green with Starlab.

PASSIONATE FOR SUSTAINABILITY IN SCIENCE.



Passionate for science.



A company with a mission.

OUR GREEN MISSION

When you take one step forward in the lab, you should not have to take two steps back where environmental protection is concerned. Only when scientific progress is also sustainable, will it unfold long-term benefits for all of us in terms of climate issues.


That is why we at Starlab have been thinking about sustainability for some time, with long-term corporate planning and durable products. We have already succeeded in saving plastic with our TipOne® pipette tip system. However, the truth is that we all still have a long way to go towards a fully sustainable laboratory.

For you and for our planet, our mission is to continually look for intelligent, climate-friendly products and processes. Together with you, we will strive to achieve our goal of becoming a thoroughly green company. Our promise along the way: to continue to take steps wherever we can – whether large or small – to reduce our impact on the environment.

Sustainability is at the heart of our business culture, and we hope this will help us attract the best talent. Through this mission, Starlab will incentivise and reward ecologically sustainable behaviour in everyday business tasks, and we will actively support our employees' engagement with sustainability initiatives in the local or wider communities.

“We need to take responsibility. It requires self-reflection, engagement and the will to make a change for a sustainable future. In our business we need to constantly look for and provide better solutions that direct us to take the next step for being a role model in our industry.”

Klaus Ambos
CEO Starlab International



A certified company.

Taking care of the environment and the prevention of causing negative impacts on the environment are two of the most important challenges facing businesses today. STARLAB is ISO 14001 certified for Environmental Management. At Starlab, we care enough to reduce our environmental footprint having a detailed list of Environmental Aspects with associated Environmental Impacts which we constantly work on reducing.



HERE ARE JUST A FEW OF OUR INTERNAL, SUSTAINABILITY INITIATIVES ACROSS THE STARLAB GROUP

- **OFFICES:** The paperless office ... we know a 100% paperless office is an unrealistic goal, but we do regularly review and audit our internal processes and digitise our administrative processes to reduce paper and printing to the absolute minimum.
- **WAREHOUSE:** We reuse cardboard boxes where possible. Plastic padding has been replaced with 100% biodegradable paper packing. All waste (cardboard, pallets etc) is recycled. Warehouse lighting replaced with energy-saving lamps.
- **PURCHASING:** Product is purchased for the whole Starlab group for more efficient transportation.
- **MARKETING:** Reduced printed material and increased digital solutions. Online shop provides paperless ordering process.
- **SALES TEAM:** Promotes TipOne® refills to customers as the plastic-saving and economical alternative to racked tips. Provide free empty racks for refill users. Help customers rationalise product inventory for efficient ordering and storage of products.
- **SERVICING LABORATORY:** Calibration certificates provided in pdf format. Car sharing or use of public transport to on-site clinics where possible. Reusable boxes provided for sending your pipettes to our labs for servicing.
- **PRODUCT MANAGEMENT:** Reviews and streamlines product portfolio. Audits packaging on existing products to reduce and improve where possible.
- **R&D DEPARTMENT:** Development of greener products that result in the use of less raw materials without compromising quality.

TipOne® – one system, one greener solution.

SAVING PLASTIC SINCE 1998

TipOne® was designed to be a plastic saving pipette tip system. The system is reusable, refillable and recyclable, with every part of the system being interchangeable with the rest. It was the first tip system on the market that provided an easy to use refill system.

The system was redesigned in 2016 with two aims – to update its look, and to save more plastic. Through our expertise in manufacturing, the new system now takes less plastic to manufacture, without compromise to its quality. Compared to the original system, up to 59 % less plastic (by weight, per pack), is required to make the clear refill shells, and a 6 % reduction of plastic achieved per rack. The introduction of a new, sterile refill system provides a saving of up to 63 % less plastic (per pack) compared to the equivalent, racked product.

On top of that, our four, new-generation of quality pipette tips saved up to 40 % plastic, compared to the equivalent older tip of the same volume.



REDUCE AND REUSE!

When it comes to plastic waste in the lab, it's obviously good to recycle what you can, but two other important considerations are to REDUCE the amount of plastic coming into your lab in the first place, and second, to REUSE whatever you can.

This is where TipOne® comes into its own – when using our refill systems, you do both! You reduce the amount of new plastic coming into your lab AND you are reusing TipOne® racks at the same time. Here's the amount of plastic saved with every pack of refills compared to pre-racked tips:

PLASTIC SAVED PER PACK OF REFILLS vs PER PACK OF RACKED TIPS

Tip volume	Per pack (960 tips)
10 µl & 10/20 µl XL	466 g
200 µl	546 g
300 µl	546 g
1000 µl	737 g
1250 µl XL	807 g

PLASTIC SAVED PER PACK OF STERILE REFILLS vs PER PACK OF STERILE RACKED TIPS

Tip volume	Per pack (960 tips)
10 µl & 10/20 µl XL	370 g
200 µl	370 g
300 µl	370 g
1000 µl	475 g
1250 µl XL	475 g

TipOne® IS A CONVINCING, SUSTAINABLE PIPETTE TIP SYSTEM

100%



GREEN ENERGY

TipOne® is manufactured in a factory powered entirely by wind.

UP TO 40%



PIPETTE TIPS

Our newest generation of TipOne® tips save up to 40% plastic compared to their equivalent previous version.

100%



CARDBOARD BOXES

Pack and case cartons for TipOne® are made from 100% recycled cardboard.

100



AUTOClave CYCLES

TipOne® racks are designed to be reused. We have tested them; they can withstand at least 100 autoclave cycles (with wafers inserted) without any loss of shape or quality.

UP TO 63%



PLASTIC SAVED

when using a sterile refill compared to using the racked equivalent.

ZERO



OTHER MATERIAL

TipOne® racks, refills and stack racks are made entirely from recyclable polypropylene (PP) and contain no paper labels, so recycling is quick and easy.

66%



SPACE SAVING

The 1000 µl TipOne® refill redesign provided a 66% space saving compared to the original, which saves both packaging and transportation resources.

UN LIMITED




UPCYCLING OPTIONS

Reuse your TipOne® racks for whatever creative storage solution you can come up with!

TipOne® recycling collections.

TipOne® IS THE ORIGINAL REUSABLE, REFILLABLE AND RECYCLABLE TIP SYSTEM



REUSE

TipOne® racks are designed for reuse with the refill systems. Our durable racks (with tip wafers inserted) can withstand at least 100 autoclave cycles without loss of shape or strength.

REDUCE

Save empty TipOne® racks to reuse with the refills. Aim to reuse the racks as many times possible. Why not liaise with other people in your lab or building, and set up a pool of empty racks for sharing to ensure even less waste.

REFILL

TipOne® refills are an economical alternative to pre-racked tips. They save on plastic waster, are quick and easy to use, and do not compromise the tip's cleanliness when racking up.

RECYCLE

All your polypropylene (PP) TipOne® packaging: the racks, the coloured tip wafers, the clear refill bottom shells and lids, and the Stack Rack outer layers and the lid, can be recycled.



THE STARLAB RECYCLING COLLECTION SERVICE

The TipOne® recycling collection service is available to all TipOne® customers in Germany and Austria. It's easy to use, and you can organise the regularity of collections to suit your needs.

For more information go to www.starlab.click/recycling-DE

The TipOne® recycling process.

FOR CUSTOMERS IN GERMANY AND AUSTRIA

REGISTER

You must register before you can use the service. Click here:

REGISTRATION FORM

or go to: www.starlab.click/recycling-DE

- Contact your Sales Representative or Customer Services if you need help
- Email the completed form to: info@starlab.de
- We will set up your details on our Recycling Collection Portal
- You will receive an email with your personal log-in details. (This may take up to 5 working days)

1

BUY A CASE OF TIPS

Once you're registered, you can start collecting for recycling!

- First, you must buy your racked or refill tips by the case, but don't worry – it's both cheaper and greener!
- Why do you need to do this? So you can reuse the outer case box to store and send back your recycling. Clever eh?




You can easily spot our case Article Nos. as they all end with -C

2

WHAT CAN I RECYCLE?

All TipOne® polypropylene (PP) packaging



Tips cannot be recycled. All recycling must be contamination free. TipOne® plastic only please!

When you have filled your TipOne® case with recycling, you can arrange collection

3

ARRANGE COLLECTION

Collections are arranged via our:

RECYCLING COLLECTION PORTAL

<https://login.remondis-recycling.de/versandlabel>

- Log-in using your personal log-in details
- To print a returns label, click on the "Download Shipping Label" icon
- Confirm the recycling is free from contamination
- Confirm Data Protection Regulations
- Print a separate label for each box to be shipped

4

WAVE GOODBYE

Prepare your recycling for collection:

- Seal the box and contents securely
- Affix the returns label
- Leave the box in the same area where parcels are usually collected by couriers
- Your recycling will be collected by UPS or DHL Express, and taken to our recycling partner for processing

5

REPEAT

Repeat steps 1 to 4

- And continue to Reduce, Reuse, Refill and Recycle with Starlab TipOne®


MORE INFO

www.starlab.click/recycling-DE


What happens to my recycling?




You collect your
TipOne®
polypropylene (PP)
packaging for recycling,
and arrange a collection



The recycling is taken to
the recycling company



The first step in the
recycling process is the
crushing of the TipOne®
plastic. The shredded
material is then processed
into plastic granulate



The plastic granulate
is then processed further



And made into
new plastic products

Tips on **greening** your lab.



REUSE TipOne® RACKS FOR MORE THAN JUST TIP REFILLS

You can reuse empty TipOne® racks to store tubes, caps and other small items to keep your bench or hood tidy and organised. The hinged lid and how the racks stack one on top of the other makes the racks ideal storage boxes.

THINK BEFORE YOU AUTOCLAVE!

Many products are certified. Nine times out ten they probably don't need to be autoclaved.

Think of the energy you can save.



FREEZER MANAGEMENT!

Set ultra low freezers to -70°C rather than -80°C; this won't compromise the sample's integrity (always test run with unimportant samples first!) but will save a lot of energy. Buy efficient freezers when replacing old freezers. Defrost regularly to avoid using more energy when loaded with ice. Discard old samples so freezer space is not wasted, and optimise freezer space with sample storage boxes for efficient energy use.

BUY IN BULK!

Buy by the case (Starlab Article Nos ending with -C). It's greener, cheaper, and saves both transport emissions and packing materials.



SWITCH IT OFF!

Equipment not being used should be switched off. Even on standby, equipment can significantly increase the energy used.

CLOSE IT!

Close the fume hood sash when not in use. Even for fixed flow hoods, this can reduce the amount of heated or cooled air escaping the room.

MORE INFO

www.starlab.click/newsroom

See our online **NEWSROOM** for more tips on greening your lab and to read our articles on sustainability.

Our greenest products.

REDUCE



TipOne® refills.
(non-sterile/sterile) significantly reduce single-use plastic waste. Ensures the reuse of tip racks. Space-saving which reduces packaging needs and transport emissions.



ErgoOne® Pipettes.
Greased system reduces the waste and cost associated with annual seal and o-ring replacement. Stainless steel tip cone on micro models prevents breakage and further part needs.



StarGuard® Gloves.
Comfort & Sensitive packed in boxes of 250 & 200 gloves. Up to 40% more gloves in the same size boxes as 100 gloves to reduce packaging needs and transport emissions.



GuardOne® Workstations.
Designed with the environment in mind. Minimal packing volume and weight reduces transport needs. Latest fan technology and energy-saving LED light provides one of the lowest power consumptions on the market. Long-life UV lamp.

MORE INFO

www.starlab.click/sustainability

REUSE



TipOne® Racks.
Reuse with the refill systems over and over again. The sturdy racks can withstand at least 100 autoclave cycles. Great storing small lab items too!



Pipette Servicing Mailing Boxes.
Reuse again and again to mail-in your pipettes for servicing. Reduces waste.

RECYCLE



Centrifuge Tubes.
Racked tubes in recyclable cardboard trays. No hard to recycle polystyrene!



TipOne® tips.
All polypropylene (PP) packaging components of the system are recyclable.

The Starlab green team.

THE STARLAB GREEN TEAM

Green Teams exist in each country within the Starlab Group (UK, Germany, France and Italy). The GREEN TEAM members are passionate about sustainability and green issues; their role is to research, identify, share and galvanise ideas for improving Starlab's sustainability ratings, be that for both our products, our internal processes or initiatives with the local or wider communities.

If you would like to discuss ways of making your lab greener, or have a suggestion on how we can improve or help you reach your own sustainable goals, please let us know. We are always happy to talk “all things sustainability”, and are open to ideas and suggestions or ways in which Starlab could partner with you so together, we can build a greener future in the lab.

Please contact your area Account Manager, or you can contact the **Starlab Green Officer** at: vanessa.whight@starlab.de



● **STARLAB LOCATIONS:**
France, Italy, United Kingdom
and Germany.

**Starlab International GmbH**

Neuer Höltingbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlabgroup.com

Starlab (UK), Ltd

5 Tanners Drive
Milton Keynes MK14 5BU
United Kingdom

T: +44 (0)1908 283800
F: +44 (0)1908 283802
infoline@starlab.co.uk
www.starlab.co.uk

Starlab GmbH

Neuer Höltingbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlab.de

Starlab FRANCE SARL

30 Rue Jean Rostand
91400 Orsay
France

T: +33 (0)1 60 13 71 70
F: +33 (0)1 69 41 48 65
info@starlab.fr
www.starlab.fr

Starlab S.r.l.

Via Pinturicchio 1
20133 Milano
Italy

T: +39 (0)2 7020 1040
F: +39 (0)2 7020 1033
info@starlab.it
www.starlab.it



THE UNIVERSAL PIPETTE TIP SYSTEM: ONE SYSTEM, ONE **GREENER** SOLUTION

TipOne® is a plastic-saving, eco-friendly pipette tip system that is designed to be reusable, refillable and recyclable.



Scan QR code for more
information & ordering

See our online Newsroom for more tips on greening your lab and to read our articles on sustainability.

NEWSROOM

www.starlab.click/newsroom

For more info on Starlab's sustainability initiatives:
www.starlab.click/sustainability