

X-Clear Advanced Polyolefin StarSeal (qPCR)



Single coated tape consisting of a 0.05mm clear polyolefin film coated on one side with a clear pressure sensitive silicone-based adhesive. The tape is supplied on a white polyester release liner.

Features

- /// Suitable for use in bio-analytical research
- /// For use on polypropylene, polystyrene or polycarbonate plates
- /// Compatible with aqueous solutions and organic solvents
- /// Compatible with real-time thermal cyclers
- /// Suitable for use between -70°C and 110°C
- /// Minimal evaporation of microplate contents over use conditions
- /// Superior optical properties; best for visual microscopic analyses
- /// Minimal cross-well contamination
- /// Low autofluorescence
- /// Easy to handle delivery system, low tack to touch
- /// Certified RNase, DNase, DNA and Endotoxin free

Shelf Life

Product as supplied in original packaging will maintain stated test properties for a period of two years from date stamped on shipping container when stored at a temperature between 10–27°C (50–80°F) and a relative humidity between 40–60%.

Ordering Information

Cat. No.	Description	Pack Size
E2796-9795	X-Clear Advanced Polyolefin StarSeal (qPCR), Self-Adhesive	100

STARLAB International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlabgroup.com

STARLAB GmbH
Neuer Höltigbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlab.de

STARLAB (UK), Ltd
4 Tanners Drive
Milton Keynes MK14 5NA
United Kingdom

T: +44 (0)1908 283800
F: +44 (0)1908 283802
info@starlab.co.uk
www.starlab.co.uk

STARLAB FRANCE SARL
30 rue Jean Rostand
91400 ORSAY
France

T: +33 (0)1 60 1371 70
F: +33 (0)1 69 4148 65
info@starlab-france.com
www.starlab-france.com

STARLAB S.r.l.
Via Pinturicchio 1
20133 Milano
Italy

T: +39 02 7020 1040
F: +39 02 7020 1033
info@starlabitalia.it
www.starlabitalia.it

www.starlabgroup.com

