

Freezing racking guide.

FIND THE RIGHT RACK FOR YOUR FREEZER

UPRIGHT FREEZERS

EPENDORF (NEW BRUNSWICK)

- F570h CryoCube
- F740h comp3 CryoCube
- U360VIP
- U410
- U535VIP
- U570
- U725VIP

PHCBI (PANASONIC / SANYO)

- MDF-DU302VX
- MDF-DU502VH-PE
- MDF-DU502VX
- MDF-DU702VH-PE
- MDF-DU702VX
- MDF-DU900V-PE
- MDF-U33V
- MDF-U3386S
- MDF-U4186S
- MDF-U53V
- MDF-U55V
- MDF-U500VX
- MDF-U5386S
- MDF-U700VX
- MDF-U7386S
- MDF-U76V

CHEST FREEZERS

EPENDORF (NEW BRUNSWICK)

- C340
- C585VIP
- C660
- C760VIP

PHCBI (PANASONIC / SANYO)

- MDF-C8V1
- MDF-192
- MDF-193
- MDF-393
- MDF-394
- MDF-594
- MDF-794



UPRIGHT FREEZERS



1 Stainless Steel Racking



2 Aluminium Racking with Side Access



3 Stainless Steel Racking with Pull Out Shelves

The freezer racking listed is suitable for use with any standard storage box up to 136 x 136 x 55 mm.



EPENDORF (NEW BRUNSWICK)


	F570h CryoCube			F740h comp3 CryoCube			U360VIP			U410		
	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S.Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S. Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S. Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S. Steel, Pull Out Shelves
Freezer racking Cat. No.	I2602-4400	I2612-4400	I2632-4400	I2602-4700	I2612-4700	I2632-4700	I2602-4700	I2612-4700	I2632-4700	I2602-4400	I2612-4400	I2632-4400
Racks per freezer	25	25	25	18	18	18	9	9	9	15	15	15
Boxes per rack	16 (4x4)	16 (4x4)	16 (4x4)	28 (4x7)	28 (4x7)	28 (4x7)	28 (4x7)	28 (4x7)	28 (4x7)	16 (4x4)	16 (4x4)	16 (4x4)
Boxes per freezer	400	400	400	504	504	504	252	252	252	240	240	240

EPPENDORF (NEW BRUNSWICK)

	U535VIP			U570			U725VIP		
	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S.Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S.Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S. Steel, Pull Out Shelves
Freezer racking Cat. No.	I2602-4400	I2612-4400	I2632-4400	I2602-4400	I2612-4400	I2632-4400	I2602-4700	I2612-4700	I2632-4700
Racks per freezer	12	12	12	25	25	25	18	18	12
Boxes per rack	28 (4x7)	28 (4x7)	28 (4x7)	16 (4x4)	16 (4x4)	16 (4x4)	28 (4x7)	28 (4x7)	28 (4x7)
Boxes per freezer	336	336	336	400	400	400	504	504	336



HOW TO MEASURE UPRIGHT FREEZER RACKING



See website for the external dimensions of the freezer racking listed on this guide.

PHCBI (PANASONIC / SANYO)


	MDF-DU302VX			MDF-DU502VH-PE MDF-DU502VX			MDF-DU702VH-PE MDF-DU702VX			MDF-DU900V-PE					
	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S.Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S. Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Aluminium, Side Access	3 S.Steel, Pull Out Shelves	1 Stainless Steel, Side Access		2 Aluminium, Side Access		3 Stainless Steel, Pull Out Shelves	
Freezer racking Cat. No.	I2602-4500	I2612-4500	I2632-4500	I2602-4600	I2612-4600	I2632-4600	I2602-4600	I2612-4600	I2632-4600	I2602-4500	I2602-4600	I2612-4500	I2612-4600	I2632-4500	I2632-4600
Racks per freezer	12	12	24	16	16	16	24	24	24	14	14	14	14	14	14
										28 Total		28 Total		28 Total	
Boxes per rack	20 (4x5)	20 (4x5)	20 (4x5)	24 (4x6)	24 (4x6)	24 (4x6)	24 (4x6)	24 (4x6)	24 (4x6)	20 (4x5)	24 (4x6)	20 (4x5)	24 (4x6)	20 (4x5)	24 (4x6)
										44 Total		44 Total		44 Total	
Boxes per freezer	240	240	240	384	384	384	576	576	576	280	336	280	336	280	336
										616 Total		616 Total		616 Total	

UPRIGHT FREEZERS

PHCBI (PANASONIC / SANYO)

	MDF-U33V						MDF-U3386S					
	1 Stainless Steel, Side Access		2 Aluminium, Side Access		3 Stainless Steel, Pull Out Shelves		1 Stainless Steel, Side Access		2 Aluminium, Side Access		3 Stainless Steel, Pull Out Shelves	
Freezer racking Cat. No.	I2602-4400	I2602-4500	I2612-4400	I2612-4500	I2632-4400	I2632-4500	I2602-4400	I2602-5400	I2612-4400	I2612-5400	I2632-4400	I2632-5400
Racks per freezer	6	6	6	6	6	6	6	6	6	6	6	6
	12 Total		12 Total		12 Total		12 Total		12 Total		12 Total	
Boxes per rack	16 (4x4)	20 (4x5)	16 (4x4)	20 (4x5)	16 (4x4)	20 (4x5)	16 (4x4)	20 (4x5)	16 (4x4)	20 (4x5)	16 (4x4)	20 (4x5)
	36 Total		36 Total		36 Total		36 Total		36 Total		36 Total	
Boxes per freezer	96	120	96	120	96	120	96	120	96	120	96	120
	216 Total		216 Total		216 Total		216 Total		216 Total		216 Total	

	MDF-U4186S						MDF-U53V MDF-U55V MDF-U500VX					
	1 Stainless Steel, Side Access		2 Aluminium, Side Access		3 Stainless Steel, Pull Out Shelves		1 Stainless Steel, Side Access		2 Aluminium, Side Access		3 Stainless Steel, Pull Out Shelves	
Freezer racking Cat. No.	I2602-3400	I2602-3500	I2612-3400	I2612-3500	I2632-3400	I2632-3500	I2602-4500	I2602-4600	I2612-4500	I2612-4600	I2632-4500	I2632-4600
Racks per freezer	8	8	8	8	8	8	8	8	8	8	8	8
	16 Total		16 Total		16 Total		16 Total		16 Total		16 Total	
Boxes per rack	12 (3x4)	15 (3x5)	12 (3x4)	15 (3x5)	12 (3x4)	15 (3x5)	20 (4x5)	24 (4x6)	20 (4x5)	24 (4x6)	20 (4x5)	24 (4x6)
	27 Total		27 Total		27 Total		44 Total		44 Total		44 Total	
Boxes per freezer	96	120	96	120	96	120	160	192	160	192	160	192
	216 Total		216 Total		216 Total		352 Total		352 Total		352 Total	



UPRIGHT FREEZERS

PHCBI (PANASONIC / SANYO)

	MDF-U5386S			MDF-U7386S			MDF-U76V MDF-U700VX					
	1 S. Steel, Side Access	2 Alu, Side Access	3 S.Steel, Pull Out Shelves	1 S. Steel, Side Access	2 Alu, Side Access	3 S. Steel, Pull Out Shelves	1 Stainless Steel, Side Access		2 Aluminium, Side Access		3 Stainless Steel, Pull Out Shelves	
Freezer racking Cat. No.	I2602-4400	I2612-4400	I2632-4400	I2602-4500	I2612-4500	I2632-4500	I2602-4500	I2602-4600	I2612-4500	I2612-4600	I2632-4500	I2632-4600
Racks per freezer	16	16	16	24	24	24	12	12	12	12	12	12
							24 Total		24 Total		24 Total	
Boxes per rack	20 (4x5)	20 (4x5)	20 (4x5)	20 (4x5)	20 (4x5)	20 (4x5)	20 (4x5)	24 (4x6)	20 (4x5)	24 (4x6)	20 (4x5)	24 (4x6)
							44 Total		44 Total		44 Total	
Boxes per freezer	320	320	320	480	480	480	240	288	240	288	240	288
							528 Total		528 Total		528 Total	



CHEST FREEZERS



4. Stainless Steel Racking



5. Stainless Steel Racking with Spring Lock



6. Aluminium Racking

The freezer racking listed on this flyer is suitable for use with any standard storage box up to 136 x 136 x 55 mm.



EPPENDORF (NEW BRUNSWICK)

	C340			C585VIP			C660			C760VIP		
	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium
Freezer racking Cat. No.	I2602-1300	I2613-1300	I2612-1300	I2602-1300	I2613-1300	I2612-1300	I2602-1300	I2613-1300	I2612-1300	I2602-1300	I2613-1300	I2612-1300
Racks per freezer	13	13	13	13	13	13	13	13	13	13	13	13
Boxes per rack	20	20	20	32	32	32	40	40	40	40	40	40
Boxes per freezer	260	260	260	416	416	416	520	520	520	520	520	520

CHEST FREEZERS

PHCBI (PANASONIC / SANYO)

	MDF-C8V1			MDF-192 MDF-193		
	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium
Freezer racking Cat. No.	I2602-7000	I2613-7000	I2612-7000	I2602-7000	I2613-7000	I2612-7000
Racks per freezer	6	6	6	9	9	9
Boxes per rack	7	7	7	7	7	7
Boxes per freezer	42	42	42	63	63	63

	MDF-393 MDF-394			MDF-594			MDF-794		
	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium	4. Stainless Steel	5. S. Steel Spring Lock	6. Aluminium
Freezer racking Cat. No.	I2602-9000	I2613-9000	I2612-9000	I2602-1300	I2613-1300	I2612-1300	I2602-1300	I2613-1300	I2612-1300
Racks per freezer	21	21	21	24	24	24	36	36	36
Boxes per rack	9	9	9	13	13	13	13	13	13
Boxes per freezer	189	189	189	312	312	312	468	468	468

HOW TO MEASURE CHEST FREEZER RACKING



See website for the external dimensions of the freezer racking listed on this guide