

5.0 ml Pipette Tip for ErgoOne® E Pipettes.



starlab.click/largevolumetip

SPECIALLY DESIGNED FOR OPTIMAL TIP FIT ON ERGOONE® E PIPETTES

Ideal for delicate pipetting out of tall and narrow vessels, these new 5.0 ml tips are provided in an easy-to-handle round rack.

- Special slim design: 160 mm long and approx. 9.6 mm in diameter, these tips are perfect for pipetting from tall and narrow vessels
- Available in an easy-to-handle round rack or in bags
- Compatible with ErgoOne® E, Brand® and Finnpiquette®
- Racked tips are certified RNase, DNA and ATP free
- All components made of recyclable polypropylene

Article No.	Description	Pack Size
E1070-2596	5.0 ml Pipette Tip, Loose	200
E1070-2601	5.0 ml Pipette Tip, Loose	5 x 200
E1070-2607	5.0 ml Pipette Tip, Racked	28-Tip Rack