

Rigid PCR Plate – Cyclers Compatibility.



Continued next page



	96-Well Semi-Skirted	96-Well Skirted, Low-Profile	384-Well Skirted
APPLIED BIOSYSTEMS® (ABI)			
3100, 3130XL 3500, 3500XL 3730, 3730XL			●
GeneAmp: 2700, 2720, 9600, 9700	●		
MiniAmp MiniAmp Plus	●		
ProFlex	●		
Proflex Multiblock system			●
QuantStudio: 3, 5, 6, 7, 12 K Real-Time PCR Systems	●		●
SimpliAmp	●		
Veriti 384-well			●
Veriti 96-well (with 0.2 ml block)	●		
ViiA7 7900 HT FAST			●
AGILENT®			
AriaMx qPCR System		●	
Mx3000P Mx3005P Mx4000	●		
Robocycler Gradient	●	●	
Surecycler 8800	●		●

	96-Well Semi-Skirted	96-Well Skirted, Low-Profile	384-Well Skirted
ANALYTIK JENA®			
Flexcycler2	●	●	●
T1 Thermocycler, T Gradient, T One, T Advanced, T Professional series	●	●	●
qTOWER3 (G, touch) Topical	●		●
BIOER® / BULLDOG®			
GeneExplorer	●		
GeneTouch 384			●
GeneTouch, GeneTouch Plus	●		
LifeECO LifeTouch LifePro	●		
BIOMETRA®			
Biometra T Advanced			●
Biometra T Advanced: 96, 96 G, 96 S, 96 SG	●	●	
Biometra T: One, Professional series	●	●	
Biometra Uno	●	●	
BIO-RAD® / MJ RESEARCH®			
C1000 C1000 Touch S1000	●	●	●
CFX Connect CFX96 CFX96 Touch		●	

	96-Well Semi-Skirted	96-Well Skirted, Low-Profile	384-Well Skirted
BIO-RAD® / MJ RESEARCH® continued			
CFX96 Touch Deep Well	●	●	
CFX384 CFX384 Touch			●
Chromo4	●	●	
DNA Engine Opticon Opticon 2		●	
DNA Engine DNA Engine Tetrad DNA Engine Dyad Dyad Disciple Tetrad 2		●	●
iCycler	●		
iQ, iQ-4, iQ-5 myiQ PCR System	●		
PTC200, 220, 221, 225, 240	●		●
MyCycler	●		
PTC 200		●	●
T100	●		
CORBETT® / QIAGEN®			
Palm Cycler	●		
Palm Cycler 384			●

Starlab recognises that the trademarks and names in this table belong to their respective owners.


Rigid PCR Plate – Cyclers Compatibility.

	96-Well Semi-Skirted	96-Well Skirted, Low-Profile	384-Well Skirted
EPPENDORF®			
MasterCycler ep 384 / Pro 384			●
MasterCycler: ep, ep gradient, ep realplex Pro, Pro S	●	●	
Nexus, Nexus gradient, Nexus SX1, Nexus GSX1	●	●	
GE / MEGABACE®			
MegaBACE 4000			●
MegaBACE 5000 MegaBACE 1000 Mark II		●	
ROCHE®			
Lightcycler 480 (384 well)			●
TECHNE® / COLE-PARMER®			
Flexigene	●	●	●
Genius Genius Quad	●	●	●
PCRmax Alpha Cycler 1, 2, 4			●
Prime PrimeG Prime Elite Prime Elite Satellite	●		●
Quantica		●	●
TC412, TC512		●	●
TC-PLUS	●	●	
Touchgene Touchgene Gradient	●	●	

	96-Well Semi-Skirted	96-Well Skirted, Low-Profile	384-Well Skirted
OTHER BRANDS & MODELS			
Apollo ATC-201	●		
BIONEER AllInOneCycler	●		
DLab TC1000-G	●	●	
G-Storm GS1, GS2, GS4	●	●	
HiMedia Prima-96	●		
Idaho Technology LightScanner		●	●
peqLab peqSTAR 96X	●	●	
Sensoquest Labcycler	●		●
Takara Dice touch, Gradient	●		
ThermoFisher Hybaid MBS, Px2/PxE	●		
VWR Collection UNO96	●	●	
Starlab recognises that the trademarks and names in this table belong to their respective owners.			

Rigid PCR Plates

Discover our Rigid-Frame PCR Plates.



starlab.click/pcrplates

