

Product Change Notice

1.2 ml Dilution Tube Racks with individual and strip tubes⁶ New molded alphanumeric


SUMMARY OF CHANGES

The 1.2 ml Dilution Tube Racks have been modified to incorporate the alphanumeric indexing directly during the moulding process (previously the indexing was pad printed).

- Indexing is now permanently molded into the rack surface (no printed ink).
- Production is simplified, improving efficiency and lead times.
- Removal of solvent-based inks reduces VOC emissions and improves the environmental profile.
- These changes are limited to the indexing and have no effect on FIT, FORM or FUNCTION of the racks or tubes.

| Article No. | Description |
|-------------|---|
| E1750-0400 | 1.2ml Reaction Tubes in Rack (96 Individual Tubes) |
| E1750-0410 | 1.2ml Reaction Tubes in Rack (96 Individual Tubes) (Sterile) |
| E1760-8400 | 1.2ml Reaction Tubes in Rack (12 x 8-Strip Tubes) |
| E1760-8410 | 1.2ml Reaction Tubes in Rack (12 x 8-Strip Tubes) (Sterile) |
| E1770-9400 | 1.2ml Reaction Tubes in Rack (8 x 12-Strip Tubes) |
| E1770-9410 | 1.2 ml Reaction Tubes in Rack (8 x 12- Strip Tubes) (Sterile) |

Passionate for science.



Effective from lot: 2209X (and all subsequent lots)

