


TipOne® Robotic Tips for Tecan®.


Starlab's TipOne® robotic tips require no instrument hardware or software adjustment. These high-quality robotic tips are seated in white hanging racks for use with Tecan® / Packard® robots with Varispan 4 or 8 head probes. The racks are supplied in individual blister packs to ensure integrity of product.

CERTIFIED









Compatible with the following robots with Varispan 4 or 8 head probes:







- Tecan® Freedom EVO®
- Tecan® Genesis
- Tecan® MiniPrep
- Packard® MultiProbe® II
- Packard® MultiProbe® IIHT/EX



Also compatible with:

Qiagen® QIAgility® (CAS-1200) with single-channel head and Aurora Biomed® Systems with 8-channel head.

Natural Tips	Article No.
 50 µl Tip	R1181-5700 (24 x 96-Tip Racks)
 50 µl Filter Tip	R1128-5810 (24 x 96-Tip Racks)
 200 µl Tip	R1181-0700 (24 x 96-Tip Racks)
 200 µl Filter Tip	R1128-0810 (24 x 96-Tip Racks)
 1000 µl Tip	R1181-2700 (24 x 96-Tip Racks)
 1000 µl Filter Tip	R1128-2810 (24 x 96-Tip Racks)

Conductive Tips	Article No.
 50 µl Conductive Tip	R1181-5709 (24 x 96-Tip Racks)
 50 µl Conductive Filter Tip	R1128-5819 (24 x 96-Tip Racks)
 200 µl Conductive Tip	R1181-0709 (24 x 96-Tip Racks)
 200 µl Conductive Filter Tip	R1128-0819 (24 x 96-Tip Racks)
 1000 µl Conductive Tip	R1181-2709 (24 x 96-Tip Racks)
 1000 µl Conductive Filter Tip	R1128-2819 (24 x 96-Tip Racks)

 www.starlab.click/robotictips

Passionate for science.

