

Cell scraper with pivoting blade.

300 mm SCRAPER WITH PIVOTING BLADE

Make the process of scraping off and collecting cells both easier and effective! The freely rotating blade allows you to twist it in the desired direction by applying slight pressure on the floor of the container you are scraping. The long 300 mm handle is made from flexible ABS for easy handling inside flasks. Individually wrapped and gamma sterilised to SAL 10^{-6} to ensure they meet the needs of even the most sensitive cell culture applications.

PRODUCT FEATURES

- Freely rotating 20 mm blade
- Blade made from clear polyethylene
- Flexible ABS handle for easy handling in flasks
- Gamma sterilised (SAL 10^{-6})
- RNase and DNase free
- Non-cytotoxic and non-pyrogenic
- Individually wrapped
- Manufactured in a Class 100,000 clean environment under ISO13485 & ISO9001 quality management systems



DESCRIPTION

Cell Scraper, 300 mm,
20 mm Pivoting Blade
(Sterile)

PACK SIZE

150 scrapers

ARTICLE NO.

CC7600-0300



www.starlab.click/cellscrapers