

# Stainless steel freezer racking.

IT'S COOL TO BE ORGANISED

MADE FROM  
**85%**  
RECYCLED  
MATERIAL



Passionate for science.



# Stainless steel freezer racking.

QUALITY RACKING UNITS TO OPTIMISE YOUR FREEZER CAPACITY

For use with all standard storage boxes up to **136 x 136 x 53 mm**

Racking with pull out shelves also available for use with our popular **StarStore 100** boxes.

Folding handle on front for easy removal from the freezer



## UPRIGHT FREEZER RACKS, SIDE ACCESS

Side access shelving system for optimising space in your upright freezer

Grippers on the front are designed for ease of use even with specialist gloves



## UPRIGHT FREEZER RACKS, PULL OUT SHELVES

Pull out shelves provide the convenience of access to your samples without having to remove the rack from the freezer



Folding handle on top for easy removal from the freezer



Locking rod keeps sample boxes securely in the rack

## CHEST FREEZER TOWER RACKS

Tower racks for optimising storage space and easy access to your samples from your chest freezer

Storage boxes not included.

## BE COOL, GET ORGANISED!

Starlab's new range of freezer racking provides a lightweight but sturdy solution to keep your samples safe. Manufactured to the highest specifications, laser cutting technology results in precision racking with no sharp edges.

Keeping your freezer tidy and optimising storage space is an efficient use of energy. Our racks are made with 85% recycled stainless steel, in a factory powered 100% by sustainable energy.

So get your samples organised with Starlab's new range of freezer racking. To help, we have even stamped the article number on the rack to make reordering more racks easier. Oh, and we have very short lead times too. Now that is organised!



### SMOOTH

Laser cutting technology results in high quality units with no sharp edges




### BRANDED

Article No. etched onto every rack for easy reordering

## SUSTAINABILITY JUST GOT COOLER!

- Freezer racking helps optimise capacity by keeping freezers tidy and maximising storage space for efficient energy use
- Made from 85% recycled stainless steel
- Manufactured in a factory powered by 100% sustainable energy
- Made in Germany. Short journey from factory to customer reduces carbon footprint

#EcolutionMovement



# Stainless steel racking for upright freezers.

## 1 RACKS WITH SIDE ACCESS



[starlab.click/sideracks](https://starlab.click/sideracks)

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No
9-Place	3 deep x 3 high	139 x 426 x 169 mm	E6000-0033
12-Place	3 deep x 4 high	139 x 426 x 224 mm	E6000-0034
15-Place	3 deep x 5 high	139 x 426 x 297 mm	E6000-0035
12-Place	4 deep x 3 high	139 x 565 x 169 mm	E6000-0043
16-Place	4 deep x 4 high	139 x 565 x 224 mm	E6000-0044
20-Place	4 deep x 5 high	139 x 565 x 279 mm	E6000-0045
24-Place	4 deep x 6 high	139 x 565 x 334 mm	E6000-0046

Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.



For use with all storage boxes up to **136 x 136 x 53 mm.**



## 2 RACKS WITH PULL OUT SHELVES




[starlab.click/shelfracks](https://starlab.click/shelfracks)

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No
9-Place	3 deep x 3 high	140 x 450 x 173 mm	E6001-0033
12-Place	3 deep x 4 high	140 x 450 x 231 mm	E6001-0034
15-Place	3 deep x 5 high	140 x 450 x 288 mm	E6001-0035
12-Place	4 deep x 3 high	140 x 569 x 174 mm	E6001-0043
16-Place	4 deep x 4 high	140 x 569 x 231 mm	E6001-0044
20-Place	4 deep x 5 high	140 x 569 x 281 mm	E6001-0045
24-Place	4 deep x 6 high	140 x 569 x 337 mm	E6001-0046
28-Place	4 deep x 7 high	140 x 569 x 403 mm	E6001-0047

Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.

## HOW TO MEASURE RACKING FOR UPRIGHT FREEZERS



Storage boxes are sold separately and are not included with any of the freezer racking.

## 3 RACKS WITH PULL OUT SHELVES

For use with **StarStore 100** storage boxes



Capacity	Format	Ext. Dimensions (WxDxH)	Art. No
9-Place	3 deep x 3 high	147 x 500 x 179 mm	E6001-0133
12-Place	3 deep x 4 high	147 x 500 x 239 mm	E6001-0134
16-Place	4 deep x 4 high	147 x 655 x 239 mm	E6001-0144

Racking suitable for use with StarStore 100 storage boxes (WxDxH) 140 x 155 x 55 mm.

[starlab.click/starstore100racks](https://starlab.click/starstore100racks)

# Racking guide: upright freezers.

## 1 SIDE ACCESS FOR STANDARD STORAGE BOXES up to 136 x 136 x 53 mm



### EPENDORF (NEW BRUNSWICK)

F570h CryoCube	F740h CryoCube	U360VIP	U535VIP	U570	U725VIP
----------------	----------------	---------	---------	------	---------

### PHCBI (PANASONIC / SANYO) - MDF RANGE

DU302VX	DU502VH-PH	DU502VX	DU702VH-PE	DU702VX	DU900V-PE	U33V	U4186S	U500VX	U53V	U5386S	U55V	U700VX	U7386S	U76V
---------	------------	---------	------------	---------	-----------	------	--------	--------	------	--------	------	--------	--------	------

### HAIER

DW86L288	DW86L338(J)	DW86L388(A/J)	DW86L486	DW86L490(J)	DW86L628	DW86L728(J)
----------	-------------	---------------	----------	-------------	----------	-------------

### THERMO

ExF60086
----------

E6000-0033	9-place (3d x 3h)	139 x 426 x 169 mm													
E6000-0034	12-place (3d x 4h)	139 x 426 x 224 mm													
E6000-0035	15-place (3d x 5h)	139 x 426 x 297 mm													
E6000-0043	12-place (4d x 3h)	139 x 565 x 169 mm													
E6000-0044	16-place (4d x 4h)	139 x 565 x 224 mm	25												
E6000-0045	20-place (4d x 5h)	139 x 565 x 279 mm													
E6000-0046	24-place (4d x 6h)	139 x 565 x 334 mm													

## 2 PULL OUT SHELVES FOR STANDARD STORAGE BOXES up to 136 x 136 x 53 mm



E6001-0033	9-place (3d x 3h)	140 x 450 x 173 mm													
E6001-0034	12-place (3d x 4h)	140 x 450 x 231 mm													
E6001-0035	15-place (3d x 5h)	140 x 450 x 288 mm													
E6001-0043	12-place (4d x 3h)	140 x 569 x 174 mm													
E6001-0044	16-place (4d x 4h)	140 x 569 x 231 mm	25												
E6001-0045	20-place (4d x 5h)	140 x 569 x 281 mm													
E6001-0046	24-place (4d x 6h)	140 x 569 x 337 mm													
E6001-0047	28-place (4d x 7h)	140 x 569 x 403 mm		18	9	12									18

## 3 PULL OUT SHELVES FOR STARSTORE 100 STORAGE BOXES 140 x 155 x 55 mm



E6001-0133	9-place (3d x 3h)	147 x 500 x 179 mm													
E6001-0134	12-place (3d x 4h)	147 x 500 x 239 mm													
E6001-0144	16-place (4d x 4h)	147 x 655 x 239 mm													

Table shows the maximum number of racks that fit into the model of freezer.

\*These freezers hold racks of different sizes on different shelves; table shows maximum number of each rack size. The maximum number of racks per freezer is the total of each added together.

Version 1: more info to be added as soon as it becomes available.

# Stainless steel racking for chest freezers.




## 4 TOWER RACKS

For use with standard storage boxes up to **136 x 136 x 53 mm**

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No
9-Place	9 high	141 x 148 x 501 mm	<b>E6010-0009</b>
10-Place	10 high	141 x 148 x 556 mm	<b>E6010-0010</b>
11-Place	11 high	141 x 148 x 611 mm	<b>E6010-0011</b>
12-Place	12 high	141 x 148 x 666 mm	<b>E6010-0012</b>
13-Place	13 high	141 x 148 x 721 mm	<b>E6010-0013</b>

Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.



**HOW TO MEASURE RACKING FOR CHEST FREEZERS**

 [starlab.click/chesttracks](http://starlab.click/chesttracks)

## Racking guide: chest freezers.

### 4 RACKING FOR STANDARD STORAGE BOXES up to 136 x 136 x 53 mm



		EPPENDORF (NEW BRUNSWICK)				PHCBI (PANASONIC / SANYO) - MDF RANGE				ARTICO	THERMO
		C340	C585VIP	C660	C760VIP	C2156-VAN	394	594	794	ULTF 420	CXF ULT2090
<b>E6010-0009</b>	9-high	141 x 148 x 501 mm					21				
<b>E6010-0010</b>	10-high	141 x 148 x 556 mm				14					
<b>E6010-0011</b>	11-high	141 x 148 x 611 mm								24	
<b>E6010-0012</b>	12-high	141 x 148 x 666 mm									33
<b>E6010-0013</b>	13-high	141 x 148 x 721 mm	13	13	13	13		24	36		

Table shows the maximum number of racks that fit into the model of freezer.



# Storage boxes for freezer racking.

FOR USE WITH FREEZER RACKING THAT HOLDS STANDARD BOXES UP TO 136 x 136 x 53 mm



## StarStore 81

81-place polypropylene storage box with transparent lift-off lid. Printed reference grid on base. Natural or neon colours.

WxDxH: 130 x 130 x 45 mm



## StarStore 49

49-place polypropylene storage box with transparent lift-off lid. Natural colour only.

WxDxH: 130 x 130 x 45 mm



## StarStore

Polypropylene storage box with no dividers. Transparent lift-off lid. Natural colour only.

WxDxH: 130 x 130 x 45mm



## Cardboard storage

Moisture resistant cardboard boxes with a choice of grid insert configurations. Boxes also available without grids.

WxDxH: 133 x 133 x 53 mm



## StarRack 64

64-place polypropylene double-sided rack with transparent lift-off lid. Neon colours.

WxDxH: 130 x 130 x 50 mm

[starlab.click/ppstorage](http://starlab.click/ppstorage)

[starlab.click/cardboardstorage](http://starlab.click/cardboardstorage)


[starlab.click/starrack](http://starlab.click/starrack)



## 81-place coloured cryo boxes

Polycarbonate storage boxes with transparent lift-off lid. 9x9 fixed, grid dividers. Printed reference grid on base.


WxDxH: 133 x 133 x 49 mm



## 100-place coloured cryo boxes

Polycarbonate storage boxes with transparent lift-off lid. 10x10 fixed, prong dividers. Printed reference grid on base.

WxDxH: 133 x 133 x 49 mm



## 81- and 100-place cryo boxes

Polycarbonate storage boxes. Transparent lift-off lid with printed reference grid. 9x9 box, fixed grid dividers; 10x10 fixed pronged dividers. White base.

WxDxH: 130 x 130 x 55 mm



## StarStore 81Z

Zylar resin storage box with transparent lift-off lid. 9x9 fixed, grid dividers. Natural colour only.

WxDxH: 130 x 130 x 55 mm

[starlab.click/cryogenicstorage](http://starlab.click/cryogenicstorage)



## StarStore 100

100-place polypropylene storage box with fixed 10x10 dividers. Hinged, lockable lid. Positioning lugs on lid for safe stacking. Natural or neon colours.

WxDxH: 140 x 155 x 55 mm

[starlab.click/starstore100](http://starlab.click/starstore100)

For use with the freezer racks that hold StarStore 100 storage boxes only.



**Starlab International GmbH**

Neuer Höltigbaum 38  
22143 Hamburg  
Germany

T: +49 (0)40 675 99 39 0  
F: +49 (0)40 675 99 39 20  
info@starlab.de  
www.starlabgroup.com

**Starlab GmbH**

Neuer Höltigbaum 38  
22143 Hamburg  
Germany

T: +49 (0)40 675 99 39 0  
F: +49 (0)40 675 99 39 20  
info@starlab.de  
www.starlab.de

**Starlab (UK), Ltd**

5 Tanners Drive  
Milton Keynes MK14 5BU  
United Kingdom

T: +44 (0)1908 283800  
F: +44 (0)1908 283802  
infoline@starlab.co.uk  
www.starlab.co.uk

**Starlab France SARL**

30 Rue Jean Rostand  
91400 Orsay  
France

T: +33 (0)1 60 13 71 70  
F: +33 (0)1 69 41 48 65  
info@starlab.fr  
www.starlab.fr

**Starlab S.r.l.**

Via Bracco 6  
20159 Milano  
Italy

T: +39 (0)2 7020 1040  
F: +39 (0)2 7020 1033  
info@starlab.it  
www.starlab.it



**IT'S COOL TO BE ORGANISED.**

Starlab's new range of freezer racking provides a lightweight but sturdy solution to keep your samples safe. Manufactured to the highest specifications, laser cutting technology results in precision racking with no sharp edges.

Keeping your freezer tidy and optimising storage space is an efficient use of energy. Our racks are made with 85% recycled stainless steel, in a factory powered 100% by sustainable energy.

Get organised, with Starlab's new range of freezer racking!



[www.starlab.click/samplestorage](http://www.starlab.click/samplestorage)