



Rotor-Gene® Style 4-Strip Tubes and Caps, 0.1ml



Product Highlights

- 0.1 ml 4-strip tubes and caps
- For use with the Rotor-Gene® Q Real-Time PCR Cycler

Product Applications

Realttime PCR (quantitative PCR = qPCR)

–

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	I1402-0400
Sterile	No
Autoclavable	No
Volume	0.1 ml
Pack Size	250 Pcs. (1 Box × 250 Pcs.)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes
For use with (SL)	Rotor-Gene® Q Real-Time PCR Cycler

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Rotor-Gene® Style 4-Strip Tubes and Caps, 0.1ml

0.1 ml 4-strip tubes and caps are for use with the Rotor-Gene® Q Real-Time PCR Cycler

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk