



TPE Sealing Mat for 48-Well Deepwell Plates



Product Highlights

- Autoclavable
- Certified RNase, DNase, DNA and PCR Inhibitor free
- Endotoxin free

Product Applications

PCR

Sample storage and preparation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E2848-0001
Pack Size	20 Pcs. (1 Pack × 20 Pcs.)
Free of endotoxins	Yes
Pyrogen free	Yes
DNA free	Yes
PCR inhibitor free	Yes
RNase free	Yes
DNase free	Yes
Easy to Pierce with Needle or Metal Probe	yes
Easy to Pierce with Pipette Tip	yes
Number of wells	48 wells
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about TPE Sealing Mat for 48-Well Deepwell Plates

Sealing mats for deepwell and microplates. Pierceable mats reseal after use. Chemical-resistant mats can be used with DMSO and other solvents.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk