



## TipOne® Repeat Dispenser Tips



### Product Highlights

- Compatible with a wider range of dispensers
- Available in eight volumes – new 0.1 ml tip!
- High accuracy and precision
- Lower dispensing forces
- ISO 8655 manufactured
- Certified RNase, DNase, DNA and Pyrogen free

## Product Applications

High-precision repetitive dispensing regardless of viscosity, vapour pressure and density of the liquid

Safe dispensing by avoidance of aerosol contaminations

Ideal for the dispensing of radioactive, infectious or toxic substances

Quick dispensing of long series

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|              |                |
|--------------|----------------|
| Art. No.     | See variations |
| Autoclavable | No             |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about TipOne® Repeat Dispenser Tips

TipOne® Repeat [Dispenser Tips](#) are available in volumes 0.1 ml to 50 ml. These high quality [dispenser](#) tips have a polypropylene cyclinder and a polyethylene piston. They are certified RNase, DNase, DNA and Pryogen free.

- › Tip volume recognition with Multipette® Plus
- › Up to 100 [dispenser](#) steps with Multipette® Plus
- › Extended piston length of most tips minimises dead air volume, improves reagent delivery and provides easier access to narrow vessels
- › High accuracy and precision
- › Lower dispensing forces
- › Available non-sterile bulk, or sterile individually wrapped
- › Colour coded adapters supplied with 25 ml and 50 ml tips
- › ISO 8655 manufactured

### TipOne® Repeat Compatibility













- › STARLAB: RepeatOne® [Dispenser](#)
- › Brand: Handystep® Manual, HandyStep® S and HandyStep® electronic
- › Denville: Faststepper®
- › Eppendorf: Multipette® 4780, Multipette® Plus\*, Multipette® stream and Multipette® Xstream\*\*
- › Gilson: Pipetman® and DistriMan®
- › Rainin: AutoRep® E, AutoRep® S and AutoRep® M
- › Ritter: Ripette®, Ripette® pro and Ripette® genX
- › Sartorius: Biohit® Mechanical Stepper, Biohit® eLine® Pro, Biohit® and eLine® Lite
- › Scilogex: Stepmate®
- › Thermo Fisher: Fisherbrand® Repeat Stepper

\* Works only with old sensor (prior to service).

\*\* For use with models manufactured prior to Dec 2013 only.

















## All Variations

| PRODUCT NAME   | PACKAGING SIZE   | ART. NO.   |
|--|--|------------|
| <br><a href="#">0.1 ml TipOne® Repeat Dispenser Tips (sterile)</a><br><b>Volume:</b> 0.1 ml<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a>     | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-0010 |
| <br><a href="#">0.5 ml TipOne® Repeat Dispenser Tips (sterile)</a><br><b>Volume:</b> 0.5 ml<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a>     | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-0050 |
| <br><a href="#">1.0 ml TipOne® Repeat Dispenser Tips (sterile)</a><br><b>Volume:</b> 1.0 ml<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a>   | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-0100 |
| <br><a href="#">2.5 ml TipOne® Repeat Dispenser Tips (sterile)</a><br><b>Volume:</b> 2.5 ml<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a> | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-0250 |
| <br><a href="#">5.0 ml TipOne® Repeat Dispenser Tips (sterile)</a><br><b>Volume:</b> 5.0 ml<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a> | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-0500 |
| <br><a href="#">10 ml TipOne® Repeat Dispenser Tips (sterile)</a><br><b>Volume:</b> 10 ml<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a>   | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-1000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)









| PRODUCT NAME  | PACKAGING SIZE   | ART. NO.   |
|---|--|------------|
|  <p><a href="#">25 ml TipOne® Repeat Dispenser Tips (sterile)</a><br/> <b>Volume:</b> 25 ml<br/> <b>Sterile:</b> Yes</p> <p> <a href="#">Create PDF Data Sheet</a></p>          | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-2500 |
|  <p><a href="#">50 ml TipOne® Repeat Dispenser Tips (sterile)</a><br/> <b>Volume:</b> 50 ml<br/> <b>Sterile:</b> Yes</p> <p> <a href="#">Create PDF Data Sheet</a></p>          | 100 Individually wrapped sachets<br>(1 Box × 100 Individually wrapped sachets) | S4761-5000 |
|  <p><a href="#">0.1 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br/> <b>Volume:</b> 0.1 ml<br/> <b>Sterile:</b> No</p> <p> <a href="#">Create PDF Data Sheet</a></p>     | 100 Tips (1 Bag × 100 Tips)  | S4751-0010 |
|  <p><a href="#">0.5 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br/> <b>Volume:</b> 0.5 ml<br/> <b>Sterile:</b> No</p> <p> <a href="#">Create PDF Data Sheet</a></p> | 100 Tips (1 Bag × 100 Tips)  | S4751-0050 |
|  <p><a href="#">1.0 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br/> <b>Volume:</b> 1.0 ml<br/> <b>Sterile:</b> No</p> <p> <a href="#">Create PDF Data Sheet</a></p> | 100 Tips (1 Bag × 100 Tips)  | S4751-0100 |
|  <p><a href="#">2.5 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br/> <b>Volume:</b> 2.5 ml<br/> <b>Sterile:</b> No</p> <p> <a href="#">Create PDF Data Sheet</a></p> | 100 Tips (1 Bag × 100 Tips)  | S4751-0250 |
|  <p><a href="#">5.0 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br/> <b>Volume:</b> 5.0 ml<br/> <b>Sterile:</b> No</p> <p> <a href="#">Create PDF Data Sheet</a></p> | 100 Tips (1 Bag × 100 Tips)  | S4751-0500 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



| PRODUCT NAME  | PACKAGING SIZE              | ART. NO.   |
|---|-----------------------------|------------|
| <br><a href="#">10 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br><b>Volume:</b> 10 ml<br><b>Sterile:</b> No<br><br> <a href="#">Create PDF Data Sheet</a> | 100 Tips (1 Bag × 100 Tips) | S4751-1000 |
| <br><a href="#">25 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br><b>Volume:</b> 25 ml<br><b>Sterile:</b> No<br><br> <a href="#">Create PDF Data Sheet</a> | 100 Tips (1 Bag × 100 Tips) | S4751-2500 |
| <br><a href="#">50 ml TipOne® Repeat Dispenser Tips (non-sterile)</a><br><b>Volume:</b> 50 ml<br><b>Sterile:</b> No<br><br> <a href="#">Create PDF Data Sheet</a> | 100 Tips (1 Bag × 100 Tips) | S4751-5000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)