



200 µl Wide Orifice Tips



Product Highlights

- Large orifice tip ideal for pipetting macromolecules and viscous samples
- Compatible with a wide range of pipettes
- Certified

Product Applications

Next generation sequencing. Cell culture. Plant sciences.

Eliminate cell fragmentation caused by mechanical shearing

Maintaining the integrity of high molecular weight DNA like gDNA or rDNA

Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup, etc.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	See variations
Filter tip	Non-filter
Volume	200 µl
Tip Type	wide orifice
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk











More informations about 200 µl Wide Orifice Tips

- › Protect delicate samples from shear force and reduce flow resistance
- › Eliminate cell fragmentation caused by mechanical shearing
- › Maintaining the integrity of high molecular weight DNA like gDNA or rDNA
- › Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup, etc.



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl Graduated Tip, Wide Orifice, Loose Sterile: No Graduations: Yes Packaging Unit: Bag  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	E1011-8000
 200 µl Graduated Tip / Wide Orifice, Racked Sterile: No Graduations: Yes Packaging Unit: Racks  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	E1011-8400
 200 µl Graduated Tip, Wide Orifice, Racked (Sterile) Sterile: Yes Graduations: Yes Packaging Unit: Racks  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	E1011-8410
 200 µl Tip, Wide Orifice, Loose, non-sterile Sterile: No Graduations: No Packaging Unit: Bag  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	I1011-8000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk