



StarChill MCT



Product Highlights

- 24-well chiller rack with lid for use with 0.5, 1.5 and 2.0 ml microcentrifuge tubes
- Temperature colour change: 7 °C
- Estimated time at 4 °C: + 3 hours*
- StarChill Rack must be frozen for a minimum of 8 hours at -20°C before use
- Freeze the StarChill Rack upside down in order for the gel to freeze evenly

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|------------------------|-------------------------|
| Art. No. | See variations |
| Dimensions (W × D × H) | 146 mm × 105 mm × 70 mm |
| Dimensions (W × D) | 146 mm × 105 mm |
| Autoclavable | No |
| Number of wells | 96 wells |
| Rack Type | Chiller Rack |
| Color | Various |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about StarChill MCT





StarChill high performance [chiller racks](#)

Maintain a fridge environment on the benchtop! The centre section of StarChill racks changes colour when frozen. It then gradually turns back to its original colour as it warms up, thus giving an instant, visual-indication of the temperature.

*Time will vary, depending on ambient temperature. All StarChill racks must be frozen for a minimum of 8 hours at -20 °C before use. Not Autoclavable.



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|------------|
|  StarChill MCT, Purple/Pink Color: Pink (warm)/purple (cold)  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | E2396-5634 |
|  StarChill MCT, Green/Yellow Color: Yellow (warm)/green (cold)  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | E2396-5636 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk