



300 µl graduated TipOne® Filter Tip, Natural (sterile)



Product Highlights

- Graduated at 100 μl, 200 μl and 300 μl
- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal tip fit
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- All components, excluding tips and shrink wrap, are recyclable

Product Applications

DNA applications (e.g. PCR)
RNA applications (e.g. Gene expression analysis)
Protein applications (e.g. Antibody Research)
Cell Culture applications (e.g. Media)
Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







General Data

Art. No.	See variations
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	1.0 - 300 µl
Volume	300 μl
Tip Type	standard
Material	Polypropylene
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







More informations about 300 μ l graduated TipOne® Filter Tip, Natural (sterile)

Nothing goes through!

Protect your samples from unwanted contamination with TipOne[®] filter tips. Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your pipette, but they can contaminate other samples later. TipOne[®] filter tips block over 99% of aerosols*, making them one of the best filters available.

Sterile TipOne® tips are validated according to EN ISO 11137-2:2007.

*See TipOne® brochure in downloads section for details.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.





Starlab (UK) LTD





All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
>	300 µl graduated TipOne® Filter Tip, Natural, Refills (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips) Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1120-9710
- 6	300 µl graduated TipOne® Filter Tip, Natural, Refills (sterile), Case Pack Size: 7,680 Tips (80 Refills × 96 Tips) Create PDF Data Sheet	7,680 Tips (80 Refills × 96 Tips)	S1120-9710-C
	300 µl graduated TipOne® Filter Tip. Natural, Racks (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips) Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1120-9810
	300 µl graduated TipOne® Filter Tip, Natural, Racks (sterile), Case Pack Size: 7,680 Tips (80 Racks × 96 Tips) Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1120-9810-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

