



## **Syringe Filter, PES (Sterile)**



#### **Product Highlights**

- High bubble point membrane provides tighter pore size for greater security
- High burst pressure ensures safe operation
- Large filtration surface area allows you to filter your sample more easily
- Low hold-up volume provides maximum recovery of your solution
- Surfactant-free membrane for cleanliness in downstream applications
- Filter type printed on the filter for ease of identification
- Convenient, larger pack sizes
- Gamma sterilised (SAL 10<sup>-6</sup>)

### **Product Applications**

0.22 µm membranes for sterilization, sterility testing, "final filtration"

 $0.45~\mu m$  membranes for analysis of bacteria, sterility testing, clarification and prefiltration

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD** 5 Tanners Drive

MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







### **General Data**

Art. No.	See variations
Sterile	Yes
Color	Yellow
Membrane	PES Membrane

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







# More informations about Syringe Filter, PES (Sterile)

#### Nearly half the hold up volume of some other syringes on the market!

STARLAB's Syringe Filters are an economical choice due to the maximum recovery of your solution. Manufactured using the latest technology, these 33 mm diameter filters are suitable for volumes 10–200 ml. Use for cell culture media and additives, biological solutions and buffers.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD** 5 Tanners Drive

MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







### **All Variations**

PRODUCT NAME		PACKAGING SIZE	ART. NO.
•	0.22 μm Syringe Filter, PES (Sterile), Yellow,  133 mm  Pore size: 0.22 μm  Create PDF Data Sheet	100 Pcs. (1 Bag × 100 Pcs.)	E4780-1226
0	0.45 μm Syringe Filter, PES (Sterile), Yellow, 33 mm Pore size: 0.45 μm	100 Pcs. (1 Bag × 100 Pcs.)	E4780-1456
	Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

