



## Universal centrifuge adapter, for 5.0 ml tubes, for 15 ml rotor bores



### Product Highlights

- Converts your 15 ml rotor for use with 5.0 ml Centrifuge Tubes
- Suitable for all 15 ml rotor bores or adapters
- Reuseable
- Freezable (-90°C)
- Made of acetal plastic

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E1450-9011
Pack Size	8 Pcs. (1 Pack × 8 Pcs.)
Material	Acetal plastic
Color	Gray
For use with (SL)	5.0 ml tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Universal centrifuge adapter, for 5.0 ml tubes, for 15 ml rotor bores

These adapters convert your centrifuge's 15 ml rotor to fit Starlab's 5.0 ml Centrifuge Tubes. These reusable adapters are compatible with all 15 ml rotor bores.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)