



200 µl Graduated Filter Tip / Wide Orifice, Racked (sterile)



Product Highlights

- Large orifice tip ideal for pipetting macromolecules and viscous samples
- Compatible with a wide range of pipettes
- Certified
- 1.5 mm inside orifice
- Supplied in convenient, press to open rack

Product Applications

Next generation sequencing. Cell culture. Plant sciences.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

| | |
|--------------------|-------------------------------|
| Art. No. | E1011-8420 |
| Pack Size | 960 Tips (10 Racks × 96 Tips) |
| Length | 50 mm |
| Sterile | Yes |
| Autoclavable | No |
| Filter tip | Filter |
| Volume | 200 µl |
| Tip Type | wide orifice |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Free of endotoxins | Yes |
| PCR inhibitor free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Graduations | Yes |
| Color | Natural |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 200 µl Graduated Filter Tip / Wide Orifice, Racked (sterile)

Protect delicate samples from shear force and reduce flow resistance with Starlab's wide orifice tips.

- › Eliminate cell fragmentation caused by mechanical shearing
- › Maintaining the integrity of high molecular weight DNA like gDNA or rDNA
- › Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup, etc.