




## Cell Scrapers



### Product Highlights

- Flexible ABS handle for easy handling in flasks
- Gamma sterilised (SAL 10<sup>-6</sup>)
- Non cytotoxic and non-pyrogenic

## Product Applications

Mechanical harvesting of cells in dishes, plates and flasks

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|               |                              |
|---------------|------------------------------|
| Art. No.      | See variations               |
| Sterile       | Yes                          |
| Pack Size     | 100 Pcs. (1 Pack × 100 Pcs.) |
| DNA free      | Yes                          |
| DNase free    | Yes                          |
| Pyrogen free  | Yes                          |
| RNase free    | Yes                          |
| Non-cytotoxic | Yes                          |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Cell Scrapers







Gently raise/lift cells from the surface with Starlab [cell scrapers](#). Pivoting blades pivot the full 360° around the handle for easy access to flask corners. Flexible ABS handles provide easy handling inside flasks. All scrapers are sterile and supplied individually wrapped.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

| PRODUCT NAME   | PACKAGING SIZE               | ART. NO.    |
|--|------------------------------|-------------|
| <br><a href="#">250 mm Cell Scraper with 10 mm Pivoting Blade</a><br><b>Blade Type:</b> Pivoting<br><b>For use with (SL):</b> small flasks<br><br> <a href="#">Create PDF Data Sheet</a>             | 100 Pcs. (1 Pack × 100 Pcs.) | CC7600-0250 |
| <br><a href="#">220 mm Cell Scraper with 11 mm Fixed Blade</a><br><b>Blade Type:</b> Fixed<br><b>For use with (SL):</b> plates, dishes<br><br> <a href="#">Create PDF Data Sheet</a>                 | 100 Pcs. (1 Pack × 100 Pcs.) | CC7600-0220 |
| <br><a href="#">300 mm Cell Scraper with 20 mm Pivoting Blade (Sterile)</a><br><b>Blade Type:</b> Pivoting<br><b>For use with (SL):</b> large flasks<br><br> <a href="#">Create PDF Data Sheet</a> | 150 Pcs. (1 Pack × 150 Pcs.) | CC7600-0300 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)