




## 5 ml Gilson-Style Graduated Macro Filter Tip (sterile)



### Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified
- DNase free, RNase free, DNA free, Pyrogen free

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                |  |
|----------------|--|
| Art. No.       | I1050-0718   |
| Pack Size      | 500 Tips (10 Racks × 50 Tips)  |
| Sterile        | Yes  |
| Filter tip     | Filter   |
| Volume         | 5 ml   |
| Material       | Polypropylene  |
| DNA free       | Yes  |
| DNase free     | Yes  |
| Human DNA free | Yes  |
| Pyrogen free   | Yes  |
| RNase free     | Yes  |
| Graduations    | No   |
| Design Tips    | Gilson-Style   |
| For use with   | <b>Abimed:</b> Discovery Comfort, LABMATE Optima<br><b>Eppendorf:</b> Reference 2, Research plus, Research pro, Xplorer<br><b>Gilson:</b> Pipetman Concept, Pipetman P, Pipetman Ultra<br><b>HTL:</b> Discovery Comfort, Labmate Pro, Optipette<br><b>Nichiryo:</b> 31N, NPLO, NPX<br><b>Rainin:</b> Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus<br><b>Socorex:</b> Acura electro 935<br><b>STARLAB:</b> ErgoOne, StarPet<br><b>Thermo Fisher Scientific Inc.:</b> Finn Focus, Finnpipette F1/F2 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 5 ml Gilson-Style Graduated Macro Filter Tip (sterile)











Large Volume Filter Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)





## All Variations

| PRODUCT NAME   | PACKAGING SIZE                | ART. NO.   |
|--|-------------------------------|------------|
| <br><a href="#">5 ml Gilson-Style Graduated Macro Filter Tip (sterile)</a><br><b>Design Tips:</b> Gilson-Style<br><b>Sterile:</b> Yes<br><b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)<br> <a href="#">Create PDF Data Sheet</a>         | 500 Tips (10 Racks × 50 Tips) | I1050-0718 |
| <br><a href="#">5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile)</a><br><b>Design Tips:</b> Eppendorf-Style<br><b>Sterile:</b> Yes<br><b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)<br> <a href="#">Create PDF Data Sheet</a>   | 500 Tips (10 Racks × 50 Tips) | I1053-0718 |
| <br><a href="#">5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile)</a><br><b>Design Tips:</b> Eppendorf-Style<br><b>Sterile:</b> Yes<br><b>Pack Size:</b> 100 Tips (2 Racks × 50 Tips)<br> <a href="#">Create PDF Data Sheet</a> | 100 Tips (2 Racks × 50 Tips)  | I1053-0728 |
| <br><a href="#">5 ml Gilson-Style Graduated Macro Filter Tip (sterile)</a><br><b>Design Tips:</b> Gilson-Style<br><b>Sterile:</b> Yes<br><b>Pack Size:</b> 100 Tips (2 Racks × 50 Tips)<br> <a href="#">Create PDF Data Sheet</a>      | 100 Tips (2 Racks × 50 Tips)  | I1050-0728 |
| <br><a href="#">5.0 ml ErgoOne E/Finn-Style Filter Tip (sterile), Racked</a><br><b>Design Tips:</b> Finn-Style<br><b>Sterile:</b> Yes<br><b>Pack Size:</b> 50 Tips (1 Rack × 50 Tips)<br> <a href="#">Create PDF Data Sheet</a>        | 50 Tips (1 Rack × 50 Tips)    | E1009-5718 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



| PRODUCT NAME  | PACKAGING SIZE              | ART. NO.   |
|---|-----------------------------|------------|
| <br><a href="#">10 ml Gilson-Style Graduated Macro Filter Tip (sterile)</a><br><b>Design Tips:</b> Gilson-Style<br><b>Sterile:</b> Yes<br><b>Pack Size:</b> 50 Tips (2 Racks × 25 Tips)<br> <a href="#">Create PDF Data Sheet</a> | 50 Tips (2 Racks × 25 Tips) | I1054-0028 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)