




10 µl Round Gel-Loading Tip, Racked (non-sterile)



Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I1010-3600
Pack Size	1,020 Tips (5 Racks × 204 Tips)
Length	56.62 mm
Sterile	No
Filter tip	Non-filter
Volume	10 µl
Tip Type	round gel-loading
Material	Polypropylene
DNase free	Yes
RNase free	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 10 µl Round Gel-Loading Tip, Racked (non-sterile)

All tips are certified RNase and DNase free. Packaging of these tips differs from the TipOne® range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk