



## Tube Adapter for 15 x 1.5/ 2.0ml Tubes



### Product Highlights

- For use with Microplate Mixer N2400-8040
- Holds 15 x 1.5/2.0 ml Microcentrifuge Tubes
- Tube Adapter can only be used in conjunction Adapter Plate N2400-8043
- Also compatible for use on Vortexer N2400-6110

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	N2400-6012
Pack Size	1 Piece (1 Box × 1 Piece)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Tube Adapter for 15 x 1.5/ 2.0ml Tubes

Convert your Microplate Mixer to mix tubes! This tube adapter holds 15 x 1.5 ml or 2.0 ml [microcentrifuge tubes](#). Please note that to use this Adapter Plate N2400-8043 is also required.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)