




Thermoblock 24 x 1.5 ml Microcentrifuge Tubes



Product Highlights

- Thermoblock for Starlab Mixer HC
- Automatic thermal block recognition simplifies block exchange
- Includes ISO-Rack

Product Applications

Sample Preparation

Purification of plasmids, DNA, RNA

cDNA synthesis

Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)

Denaturation of DNA, RNA, proteins

Labeling of DNA, RNA, proteins

Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates

Transformation of bacteria

Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	S8012-0011
Pack Size	1 Piece (1 Box × 1 Piece)
Operating temperature	4 to 35 °C
Sample capacity	24 × 1.5 ml microcentrifuge tubes
Speed Range	300 – 1,400 rpm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Thermoblock 24 x 1.5 ml Microcentrifuge Tubes

Exchangeable thermoblock for use on the Starlab Mixer HC. Includes ISO-Rack. The Iso-Rack is ideal for time-critical applications and high throughput sample preparation. Prepare your samples in your tubes, using the Iso-Rack as your bench rack. You can then transfer all 24 tubes to the thermoblock at the same time, in seconds!

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk