


8) Remove feed tubing from stand and place into clean water container

9) Draw approximately 500ml of distilled water through the pedestal or cradle ring stand and into vacuum trap to rinse the system



CytoOne® Filtration Docking Units



Pedestal shown with
500 ml Bottle Top
Filtration Unit



Cradle Ring shown with
50 ml Centrifuge Tube
Filtration System
(Stand not included)

Perfect docking!

All Bottle Top Filtration Units and both the Centrifuge Tube Filtration Systems can be used with either the Pedestal or the Cradle Ring for hands-free filtration. The Pedestal is made from powder-coated aluminum and can be left under the hood. The Cradle Ring is made from steel.

CytoOne® Filtration Docking Systems

Cat. No	Description
CC6000-1116	CytoOne® Quick-Connect™ Vacuum Filtration Pedestal
CC6000-1122	CytoOne® Vacuum Filtration Cradle Ring



CC6000-1116



CC6000-1122
(Stand not included)

www.starlab.click/filtration

