

## Breathable Films



Natural, opaque, porous self-adhesive seal. Allows effective gas exchange for cellular and bacterial cultivation, while preventing contamination Available sterile or non-sterile.

### Specifications

- ✦ Suitable for use with plate types:
  - polypropylene
  - polystyrene
- ✦ Suitable for short term storage
- ✦ Standard SBS format to fit all plates
- ✦ Temperature use: 4°C to 40°C
- ✦ Certified RNase, DNase and Pyrogen free

### Storage & Shelf Life

Store at room temperature in a dry environment, out of direct sunlight. Shelf life is three years. However, at optimal storage conditions (room temperature, low humidity) the seals will last well beyond three years.

Date of manufacture can be determined by the Lot No. date on the packaging (yymmdd). This is a 6-digit number of a total of 9 which literally represents Year (yy), Month (mm), Day (dd).

### Ordering Information

Cat. No.	Description	Pack Size
E2796-3005	Breathable Film	100
E2796-3015	Breathable Film (Sterile)	2 x 25

**STARLAB International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Germany

T: +49 (0)40 675 99 39 0  
F: +49 (0)40 675 99 39 20  
info@starlab.de  
www.starlabgroup.com

**STARLAB GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Germany

T: +49 (0)40 675 99 39 0  
F: +49 (0)40 675 99 39 20  
info@starlab.de  
www.starlab.de

**STARLAB (UK), Ltd**  
4 Tanners Drive  
Milton Keynes MK14 5NA  
United Kingdom

T: +44 (0)1908 283800  
F: +44 (0)1908 283802  
info@starlab.co.uk  
www.starlab.co.uk

**STARLAB FRANCE SARL**  
30 rue Jean Rostand  
91400 ORSAY  
France

T: +33 (0)1 60 1371 70  
F: +33 (0)1 69 4148 65  
info@starlab-france.com  
www.starlab-france.com

**STARLAB S.r.l.**  
Via Pinturicchio 1  
20133 Milano  
Italy

T: +39 02 7020 1040  
F: +39 02 7020 1033  
info@starlabitalia.it  
www.starlabitalia.it

[www.starlabgroup.com](http://www.starlabgroup.com)

