

Freezer Racking Guide

– Eppendorf (New Brunswick) Freezers –

The freezer racking listed on this flyer is suitable for use with any standard storage box up to 136 x 136 x 55mm. See below.



UPRIGHT FREEZERS

Eppendorf (New Brunswick) U360VIP

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
1. Stainless Steel, Side Access	I2602-4700	28 (4x7)	9	252
2. Aluminium, Side Access	I2612-4700	28 (4x7)	9	252
3. Stainless Steel, Pull Out Shelves	I2632-4700	28 (4x7)	9	252

Eppendorf (New Brunswick) U410

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
1. Stainless Steel, Side Access	I2602-4400	16 (4x4)	15	240
2. Aluminium, Side Access	I2612-4400	16 (4x4)	15	240
3. Stainless Steel, Pull Out Shelves	I2632-4400	16 (4x4)	15	240

Eppendorf (New Brunswick) U535VIP

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
1. Stainless Steel, Side Access	I2602-4700	28 (4x7)	12	336
2. Aluminium, Side Access	I2612-4700	28 (4x7)	12	336
3. Stainless Steel, Pull Out Shelves	I2632-4700	28 (4x7)	12	336

Eppendorf (New Brunswick) U570

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
1. Stainless Steel, Side Access	I2602-4400	16 (4x4)	25	400
2. Aluminium, Side Access	I2612-4400	16 (4x4)	25	400
3. Stainless Steel, Pull Out Shelves	I2632-4400	16 (4x4)	25	400

Eppendorf (New Brunswick) U725VIP

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
1. Stainless Steel, Side Access	I2602-4700	28 (4x7)	18	504
2. Aluminium, Side Access	I2612-4700	28 (4x7)	18	504
3. Stainless Steel, Pull Out Shelves	I2632-4700	28 (4x7)	18	504

CHEST FREEZERS

Eppendorf (New Brunswick) C340

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
4. Stainless Steel	I2602-1300	13	20	260
5. Stainless Steel, Spring Lock	I2613-1300	13	20	260
6. Aluminium	I2612-1300	13	20	260

Eppendorf (New Brunswick) C585VIP

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
4. Stainless Steel	I2602-1300	13	32	416
5. Stainless Steel, Spring Lock	I2613-1300	13	32	416
6. Aluminium	I2612-1300	13	32	416

Eppendorf (New Brunswick) C660

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
4. Stainless Steel	I2602-1300	13	40	520
5. Stainless Steel, Spring Lock	I2613-1300	13	40	520
6. Aluminium	I2612-1300	13	40	520

Eppendorf (New Brunswick) C760VIP

Freezer Rack Type	Freezer Rack Cat. No.	Boxes per Rack	Racks per Freezer	Boxes per Freezer
4. Stainless Steel	I2602-1300	13	44	572
5. Stainless Steel, Spring Lock	I2613-1300	13	44	572
6. Aluminium	I2612-1300	13	44	572

Freezer racking listed on this flyer is suitable for use with any standard storage box up to 136 x 136 x 55mm, and the following STARLAB storage boxes:

50 mm high cardboard boxes with 133x133 or 136x136 footprint

StarStore 81, StarStore 49, StarStore

81-place and 100-place polycarbonate coloured cryogenic boxes

StarStore 81Z

81-place and 100-place white polycarbonate cryogenic storage boxes

Expanding Work Rack

StarRack 64

Freezer rack images shown to illustrate style of freezer racking only.

www.starlab.co.uk



STARLAB (UK) LTD, 5 Tanners Drive, Blakelands, Milton Keynes, MK14 5BU

T: 01908 283800 F: 01908 283802 infoline@starlab.co.uk www.starlab.co.uk