



## StarChill PCR



### Product Highlights

- 96-well chiller rack with lid for use with 0.2ml single or strip PCR tubes, or one 0.2ml 96-well PCR plate
- Temperature colour change: 7 °C
- Estimated time at 4 °C: + 3 hours\*
- StarChill Rack must be frozen for a minimum of 8 hours at -20°C before use
- Freeze the StarChill Rack upside down in order for the gel to freeze evenly

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	141 mm × 99 mm × 44 mm
Dimensions (W × D)	141 mm × 99 mm
Autoclavable	No
Number of wells	96 wells
Rack Type	Chiller Rack
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about StarChill PCR





### StarChill high performance [chiller racks](#)

Maintain a fridge environment on the benchtop! The centre section of StarChill racks changes colour when frozen. It then gradually turns back to its original colour as it warms up, thus giving an instant, visual-indication of the temperature.

\*Time will vary, depending on ambient temperature. All StarChill racks must be frozen for a minimum of 8 hours at -20 °C before use. Not Autoclavable.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">StarChill PCR, Pink</a> Color: Pink (warm)/purple (cold)  <a href="#">Create PDF Data Sheet</a>	2 Pcs. (1 Pack × 2 Pcs.)	E2396-2004
 <a href="#">StarChill PCR, Yellow</a> Color: Yellow (warm)/green (cold)  <a href="#">Create PDF Data Sheet</a>	2 Pcs. (1 Pack × 2 Pcs.)	E2396-2006

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)