




StarChill SBS



Product Highlights

- 96-well chiller rack without lid for use on automated systems with 0.2 ml single or strip PCR tubes, or one 0.2 ml 96-well PCR plate
- Temperature colour change: 7 °C
- Estimated time at 4 °C: + 2 hours*
- StarChill Rack must be frozen for a minimum of 8 hours at -20°C before use
- Freeze the StarChill Rack upside down in order for the gel to freeze evenly

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Dimensions (W × D × H)	128 mm × 85 mm × 26 mm
Dimensions (W × D)	128 mm × 85 mm
Autoclavable	No
Number of wells	96 wells
Rack Type	Chiller Rack
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about StarChill SBS





StarChill high performance chiller racks

Maintain a fridge environment during automation! The centre section of StarChill racks changes colour when frozen. It then gradually turns back to its original colour as it warms up, thus giving an instant, visual-indication of the temperature.

*Time will vary, depending on ambient temperature. All StarChill racks must be frozen for a minimum of 8 hours at -20 °C before use. Not Autoclavable.



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	StarChill SBS, Pink Color: Pink (warm)/purple (cold)	2 Pcs. (1 Pack × 2 Pcs.)	E2396-5654
	 Create PDF Data Sheet		
	StarChill SBS, Green Color: Yellow (warm)/green (cold)	2 Pcs. (1 Pack × 2 Pcs.)	E2396-5656
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de