



# 10/20 µl XL graduated TipOne® Filter Tip (sterile)



#### **Product Highlights**

- Graduated at 1.0  $\mu$ l, 5.0  $\mu$ l, 10  $\mu$ l and 20  $\mu$ l
- Highly effective HDPE filter that contains no additives to interfere with samples
- Universal tip fit
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- All components, excluding tips and shrink wrap, are recyclable

# **Product Applications**

DNA applications (e.g. PCR)
RNA applications (e.g. Gene expression analysis)
Protein applications (e.g. Antibody Research)
Cell Culture applications (e.g. Media)
Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH** Neuer Höltigbaum 38

22143 Hamburg Email: info@starlab.de







### **General Data**

Art. No.	See variations
Material	Polypropylene
Autoclavable	No
Sterile	Yes
Graduations	Yes
Filter tip	Filter
Volume range	0.1 - 20 μl
Color	Natural
Volume	20 μl

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH** Neuer Höltigbaum 38

22143 Hamburg Email: info@starlab.de







# More informations about $10/20~\mu I$ XL graduated TipOne® Filter Tip (sterile)

"Nothing goes through!

Protect your samples from unwanted contamination with TipOne<sup>\*</sup> filter tips. Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your pipette, but they can contaminate other samples later. TipOne<sup>®</sup> filter tips block over 99% of aerosols\*, making them one of the best filters available.

\*See TipOne® brochure in downloads section for details."

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







# **All Variations**

PRODUCT NAME		PACKAGING SIZE	ART. NO.
<b>5</b>	10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1120-3710
<b>=</b> 6	10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile), Case Pack Size: 7,680 Tips (80 Refills × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Refills × 96 Tips)	S1120-3710-C
	10/20 µl XL graduated TipOne® Filter Tip, Racks (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1120-3810
	10/20 µl XL graduated TipOne® Filter Tip, Racks (sterile), Case Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1120-3810-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH** Neuer Höltigbaum 38

22143 Hamburg Email: info@starlab.de

