




## 20 µl Tips for Beckman Coulter®



### Product Highlights

- Tips suitable for Biomek® i-Series • Biomek® FX • Biomek® NX • Biomek® 3000 • Biomek® 4000 • Multimek®
- Engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit.
- Certified RNase, DNase, DNA and PCR Inhibitor free
- Compatible with AP96, SPAN-8, MP20
- No special set-up required!

## Product Applications

PCR Setup

Assay-Setup

Serial Dilutions

Distributing Reagents

Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Filter tip	Non-filter
Volume	20 µl
Rack Style	Pale Blue SBS Rack
Head	AP96, SPAN-8, MP20
For use with (SL)	Beckman Coulter®

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 20 µl Tips for Beckman Coulter®

Racked, [robotic tips](#) suitable for Biomek® i-Series • Biomek® FX • Biomek® NX • Biomek® 3000 • Biomek® 4000 • Multimek®

Engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit. Don't risk having your time and results compromised with defective materials and poor performance!

Certified free of RNase, DNase, DNA and PCR Inhibitors. Tested Pyrogen free.





Quality control is the essential factor which gives you the freedom to walk away from your automated workstation and get on with other tasks. You need the assurance that the tip will not let you down. STARLAB [robotic tips](#) are put through stringent and thorough testing. Every lot of tips is subject to an array of demanding tests, from the geometry of the tip itself to automation-specific requirements such as clean ejection from the liquid handling head. Every lot is required to meet or exceed the parameters necessary to give you that confidence to walk away.

Tips are tested on the original robotic equipment so there is no doubt they will perform to the same requirements as the disposables from the robot manufacturer themselves.

No special set-up required! STARLAB [robotic tips](#) require no instrument hardware or software adjustment, thereby providing a simple, safe and direct high-performance substitute.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">20 µl Tip for Beckman Coulter®, Racked (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	E1076-4400
 <a href="#">20 µl Tip for Beckman Coulter®, Racked (Sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	E1076-4410

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)