



20 µl Filter Tips for Beckman Coulter®



Product Highlights

Tips suitable for Biomek[®] i-Series • Biomek[®] FX •
Biomek[®] NX • Biomek[®] 3000 • Biomek[®] 4000 •

Multimek®

- Engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit.
- Certified RNase, DNase, DNA and PCR Inhibitor free
- Compatible with AP96, SPAN-8, MP20
- No special set-up required!

Product Applications

PCR Setup Assay-Setup Serial Dilutions Distributing Reagents Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.







General Data

Art. No.	See variations
Sterile	Yes
Filter tip	Filter
Volume	20 µl
Rack Style	Pale Blue SBS Rack
Head	AP96, SPAN-8, MP20
For use with	Beckman Coulter [®]

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.







More informations about 20 µl Filter Tips for Beckman Coulter®

Racked, robotic tips suitable for Biomek[®] i-Series • Biomek[®] FX • Biomek[®] NX • Biomek[®] 3000 • Biomek[®] 4000 • Multimek[®]

Engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit. Don't risk having your time and results compromised with defective materials and poor performance!

Certified free of RNase, DNAse, DNA and PCR Inhibitors. Tested Pyrogen free.

Quality control is the essential factor which gives you the freedom to walk away from your automated workstation and get on with other tasks. You need the assurance that the tip will not let you down. STARLAB robotic tips are put through stringent and thorough testing. Every lot of tips is subject to an array of demanding tests, from the geometry of the tip itself to automation-specific requirements such as clean ejection from the liquid handling head. Every lot is required to meet or exceed the parameters necessary to give you that confidence to walk away.

Tips are tested on the original robotic equipment so there is no doubt they will perform to the same requirements as the disposables from the robot manufacturer themselves.

No special set-up required! STARLAB robotic tips require no instrument hardware of software adjustment, thereby providing a simple, safe and direct high-performance substitute.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.







All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
Y	20 µl Filter Tip for Beckman Coulter®, Racked (Sterile) Pack Size: 960 Tips (10 Racks × 96 Tips) Profile Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	E1076-4418
Y	20 µl Filter Tip for Beckman Coulter®, Racked (Sterile), Case Pack Size: 4,800 Tips (50 Racks × 96 Tips)	4,800 Tips (50 Racks × 96 Tips)	E1076-4418-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

