




Strip Caps Microdilution Tube System



Product Highlights

- Strip Caps Microdilution Tube System
- For use with the 1.2 ml tubes of the Microdilution Tube System
- Temperature use -20 °C to 100 °C
- Certified RNase, DNase, DNA and Pyrogen free

Product Applications

Multichannel and automated pipetting

Preparation of serial dilutions

Mixing

Storage and harvesting of cells

Cell growth assays

DNA screening

Long-term storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Autoclavable	No
Material	Polyethylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Operating temperature	-20°C to 100°C
For use with (SL)	1.2 ml Microtubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Strip Caps Microdilution Tube System

Organised storage!









Microdilution tube system comprising racks containing either individual, 8-strip or 12-strip 1.2 ml tubes. Ideal for serial dilutions, mixing, storage and harvesting of cells, cell growth assays, DNA screening and long-term storage.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 8-Strip Caps for 1.2 ml Reaction Tubes No. of tubes: 8 Sterile: No  Create PDF Data Sheet	125 Pcs. (1 Box × 125 Pcs.)	E1702-8400
 8-Strip Caps for 1.2 ml Reaction Tubes (Sterile) No. of tubes: 8 Sterile: Yes  Create PDF Data Sheet	125 Pcs. (1 Box × 125 Pcs.)	E1702-8410
 12-Strip Caps for 1.2 ml Reaction Tubes No. of tubes: 12 Sterile: No  Create PDF Data Sheet	80 Pcs. (1 Box × 80 Pcs.)	E1703-9400
 12-Strip Caps for 1.2 ml Reaction Tubes (Sterile) No. of tubes: 12 Sterile: Yes  Create PDF Data Sheet	80 Pcs. (1 Box × 80 Pcs.)	E1703-9410

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de