



1000 µl XL Graduated Filter Tip / Wide Orifice, Racked (sterile)



Product Highlights


- Large orifice tip ideal for pipetting macromolecules and viscous samples
- Compatible with a wide range of pipettes
- Certified
- 1.0 mm inside orifice
- Supplied in convenient, press to open rack

Product Applications

Next generation sequencing. Cell culture. Plant sciences.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E1011-9420
Length	86.4 mm
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume	1000 µl
Pack Size	960 Tips (10 Racks × 96 Tips)
Tip Type	XL wide orifice
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 1000 µl XL Graduated Filter Tip / Wide Orifice, Racked (sterile)

- › Protect delicate samples from shear force and reduce flow resistance with Starlab's [wide orifice tips](#).
- › Eliminate cell fragmentation caused by mechanical shearing
- › Maintaining the integrity of high molecular weight DNA like gDNA or rDNA
- › Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup, etc.