



## 200 µl XL Wide Orifice Tips



### Product Highlights

- Large orifice tip ideal for pipetting macromolecules and viscous samples
- Compatible with a wide range of pipettes
- Certified
- 1.0 mm inside orifice
- Supplied in convenient, press to open rack

## Product Applications

Next generation sequencing. Cell culture. Plant sciences.

Eliminate cell fragmentation caused by mechanical shearing

Maintaining the integrity of high molecular weight DNA like gDNA or rDNA

Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup, etc.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Filter tip	Non-filter
Volume	200 µl
Tip Type	XL wide orifice
Material	Polypropylene
DNA free	Yes
DNase free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)







## More informations about 200 µl XL Wide Orifice Tips

- › Protect delicate samples from shear force and reduce flow resistance
- › Eliminate cell fragmentation caused by mechanical shearing
- › Maintaining the integrity of high molecular weight DNA like gDNA or rDNA
- › Pipetting of viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup, etc.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl XL Graduated Tip / Wide Orifice, Racked (non-sterile)</a> <b>Sterile:</b> No  <a href="#">Create PDF Data Sheet</a>	192 Tips (2 Racks × 96 Tips)	E1011-7400
 <a href="#">200 µl XL Wide Orifice Tips</a> <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	192 Tips (2 Racks × 96 Tips)	E1011-7410

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)