



## Thermoblock 24 x 1.5 ml Microcentrifuge Tubes



#### **Product Highlights**

- Thermoblock for Starlab Mixer HC
- Automatic thermal block recognition simplifies block exchange
- Includes ISO-Rack

### **Product Applications**

Sample Preparation
Purification of plasmids, DNA, RNA
cDNA synthesis
Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)
Denaturation of DNA, RNA, proteins
Labeling of DNA, RNA, proteins
Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates
Transformation of bacteria
Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







### **General Data**

Art. No.	S8012-0011
Pack Size	1 Piece (1 Box × 1 Piece)
Operating temperature	4 to 35 °C
Sample capacity	24 × 1.5 ml microcentrifuge tubes
Speed Range	300 – 1,400 rpm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH** Neuer Höltigbaum 38

22143 Hamburg Email: info@starlab.de







# More informations about Thermoblock 24 $\times$ 1.5 ml Microcentrifuge Tubes

Exchangeable thermoblock for use on the Starlab Mixer HC. Includes ISO-Rack. The Iso-Rack is ideal for time-critical applications and high throughput sample preparation. Prepare your samples in your tubes, using the Iso-Rack as your bench rack. You can then transfer all 24 tubes to the thermoblock at the same time, in seconds!

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de

