



## Thermoblock 24 x 2.0 ml Microcentrifuge Tubes Ø 12 mm



### Product Highlights

- Thermoblock for Starlab Mixer HC
- Automatic thermal block recognition simplifies block exchange
- For tubes with V-Shaped bottoms

## Product Applications

Sample Preparation

Purification of plasmids, DNA, RNA

cDNA synthesis

Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)

Denaturation of DNA, RNA, proteins

Labeling of DNA, RNA, proteins

Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates

Transformation of bacteria

Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	S8012-0013
Pack Size	1 Piece (1 Box × 1 Piece)
Operating temperature	4 to 35 °C
Sample capacity	24 × micro test tubes, diameter 12 mm
Speed Range	300 – 1,400 rpm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Thermoblock 24 x 2.0 ml Microcentrifuge Tubes Ø 12 mm

Exchangeable thermoblock for use on the Starlab Mixer HC.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)