



Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue



Product Highlights

- Blue polypropylene storage / transparent box for centrifuge tubes
- Holds 16 x 50ml tubes
- Transparent, lift-off lid
- Base can be removed for use in waterbath
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 145 x 145 x 127 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E9025-1650
Pack Size	2 Pcs. (1 Box × 2 Pcs.)
Dimensions (W × D × H)	145 mm × 145 mm × 127 mm
Dimensions (W × D)	145 mm × 145 mm
Autoclavable	Yes
Rack Type	Tube Rack
Material	Polypropylene
Color	Blue
Max. centrifuge tubes or screwcap microtubes	50ml tubes
Places	16

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue

Polypropylene storage box that holds 16 x 50ml tubes.





Tubes not included.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue Places: 16  Create PDF Data Sheet	2 Pcs. (1 Box × 2 Pcs.)	E9025-1650
 Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green Places: 36  Create PDF Data Sheet	2 Pcs. (1 Box × 2 Pcs.)	E9025-3615

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk