



## PCR Lock-Rack



### Product Highlights

- Dimensions (WxDxH): 133 x 107 x 56 mm (2 sections with lid)

## Product Applications

Sample storage

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E2344-2099
Pack Size	4 Pcs. (1 Pack × 4 Pcs.)
Dimensions (W × D × H)	133 mm × 107 mm × 56 mm
Dimensions (W × D)	133 mm × 107 mm
Autoclavable	Yes
Rack Type	PCR Lock-Rack
Color	Mixed
Max. centrifuge tubes or screwcap microtubes	Each section holds 24 × 1.5 / 2.0 ml tubes, 16 × 0.5 ml tubes and 60 × 0.2 ml tubes.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about PCR Lock-Rack

Polypropylene racks which lock together to make larger racks. Transparent lid fits over two sections locked together. Each section holds 24 x 1.5 / 2.0 ml tubes, 16 x 0.5 ml tubes and 60 x 0.2 ml tubes. Two sections locked together can hold a 96-Well PCR plate.