




9-Place Stainless Steel Rack for Chest Freezer



Product Highlights

- While Stocks last
- The racks for chest freezers have a handle on the top, plus a locking bar to prevent boxes from slipping out when lifting the rack in or out of the freezer
- Suitable for use with the STARLAB boxes listed above, or any other box up to (WxDxH): 135 x 135 x 55 mm.

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I2602-9000
Pack Size	1 Piece (1 Rack × 1 Piece)
Dimensions (W × D × H)	140 mm × 140 mm × 499 mm
Dimensions (W × D)	140 mm × 140 mm
External width	140 mm
Rack Type	Chest Freezer Rack
Material	Stainless Steel
Places	9

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 9-Place Stainless Steel Rack for Chest Freezer





These racks for chest freezers are designed for use standard storage boxes up to 136 x 136 x 54mm.

Suitable for use with STARLAB's:

- › StarStore 130 x 130 mm range of storage boxes
- › 50 mm high cardboard storage boxes
- › Cryogenic storage boxes (excl. 81-Tall)
- › StarRack 64



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	9-Place Stainless Steel Rack for Chest Freezer Places: 9	1 Piece (1 Rack × 1 Piece)	I2602-9000
	 Create PDF Data Sheet		
	13-High Aluminium Rack Places: 13	1 Piece (1 Rack × 1 Piece)	I2612-1300
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk