



## 5-Tray SS Rack for 2200 x 1.5/2.0ml Tubes



### Product Highlights

- 5-drawer stainless steel racking for upright freezers with pull out drawers
- Each drawer contains polypropylene dividers for the direct storage of tubes up to 12.5mm diameter
- Holds a total of 2200 tubes
- Max tube height: 53 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	I2670-2200
Dimensions (W × D × H)	140 mm × 595 mm × 280 mm
Dimensions (W × D)	140 mm × 595 mm
Pack Size	1 Piece (1 Rack × 1 Piece)
Rack Type	Upright Freezer Rack
Material	Stainless Steel
Shelves	5

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 5-Tray SS Rack for 2200 x 1.5/2.0ml Tubes

These racks have pull out trays which contain polypropylene dividers for the direct storage of tubes up to 12.5mm diameter, eg. 1.5/2.0ml [microcentrifuge](#) or [screw cap tubes](#).

Racks with polypropylene dividers for chest freezers are also available. Contact STARLAB for details.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)