



1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile)



Product Highlights

- Graduated at 100 µl, 250 µl, 500 µl and 1000 µl
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- Certified RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids

Dispensing liquids


Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	S1112-1830
Pack Size	960 Tips (10 Racks × 96 Tips)
Color Code Wafer	Blue
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume range	100 - 1,250 µl
Volume	1250 µl
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with	Abimed : BioPette-E, Discovery Comfort, LABMATE Optima Biohit : ePet, mLine, Proline BRAND : Transferpette, Transferpette Electronic, Transferpette S Capp : CappAero Comfort, CappAero Universal, Capptronic Single, ecopipette, Singlechannel Eppendorf : Reference, Reference 2, Research, Research plus, Research pro, Xplorer Gilson : Pipetman Concept, Pipetman Neo, Pipetman P, Pipetman Ultra Hamilton : SoftGrip Pipette, SofTouch HTL : Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette Nichiryo : 31N, 50F, 50V, NLE, NPLO, NPX Rainin : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus Sartorius : ePet, mLine, Proline, Proline Plus, Tacta Socorex : Acura 815 manual, Acura 825, Acura electro 925, Calibra digital STARLAB : ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E Thermo Fisher Scientific Inc. : Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile)

Improved by you - New System Components Overview

The perfect rack

- › Simple! The system comprises only two racks (one small, one large) for less complexity
- › Small rack for use with tip volumes 10 µl to 300 µl; large rack for use with tip volumes 1000 µl and 1,25 ml XL
- › Quick and easy identification of tips by colour-coded wafers with laser printed information and transparent lid
- › Stackable
- › Reusable
- › Autoclavable

The clever refill

- › The new refill is divided into two separate, individually-wrapped stacks of five wafers
- › Completely sealed, the tips are safe from contamination by the tight fitting top shell and unique spacers
- › Economical, space-saving design
- › Autoclavable when racked

The self-standing bag

- › New, self-standing bag for easy and convenient access and storage
- › Tamper-proof, easy-open perforation strip
- › Resealable zip-lock for contamination-free storage
- › Recyclable polyethylene bag
- › Economical – ideal for teaching labs or non-critical work

The compact stack rack

- › Next generation design with wide base for great stability
- › Transparent lid covers the tips when not in use
- › 960 tips for high throughput use, ideal for multichannel users
- › The complete stack rack is autoclavable and recyclable

The unique filter refill













- › A brand new, economical and ecofriendly TipOne® filter tip refill system
- › Fast, easy and safe: The filter tips are refilled with the top and bottom shell intact, thus providing a fast and completely protected refill system that will not compromise sterility
- › Using sterile filter refills to reload racks directly during periods of long use in hoods and PCR-workstations saves valuable space
- › Not only does the refill system generate less waste, the bottom shells are nestable so your waste takes up less storage space too

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk











All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1,250 µl XL graduated TipOne® Tip, Natural, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1112-1820
 1,250 µl XL graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 3,840 Tips (40 Racks × 96 Tips)  Create PDF Data Sheet	3,840 Tips (40 Racks × 96 Tips)	S1112-1820-C
 1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1112-1830
 1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 3,840 Tips (40 Racks × 96 Tips)  Create PDF Data Sheet	3,840 Tips (40 Racks × 96 Tips)	S1112-1830-C
 1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1112-1720
 1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1112-1720-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1,250 µl XL graduated TipOne® Tip, Natural, Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1112-1730
 1,250 µl XL graduated TipOne® Tip, Natural, Refills (sterile), Case Sterile: Yes Pack Size: 3,840 Tips (40 Refills × 96 Tips)  Create PDF Data Sheet	3,840 Tips (40 Refills × 96 Tips)	S1112-1730-C
 1,250 µl XL graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	S1112-1020
 1,250 µl XL graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 16,000 Tips (16 Bags × 1000 Tips)  Create PDF Data Sheet	16,000 Tips (16 Bags × 1000 Tips)	S1112-1020-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk